
Wadsworth Center

New York State Department of Health

TRACE ELEMENTS IN WHOLE BLOOD

Event #1, 2010

March 15, 2010

March 15, 2010

**Trace Elements in Whole Blood
Event #1, 2010**

Dear Laboratory Director:

Results from the first proficiency test (PT) event in 2010 for Trace Elements in Whole Blood have been tabulated and summarized. Target values for Arsenic, Cadmium, Mercury and Lead in whole blood have been established along with acceptable ranges. Results are graded using element-specific criteria as indicated in each narrative section. A laboratory with an unacceptable significant analytical bias relative to the target value will be expected to investigate the source of the error. A confidential three-digit code number assigned by the PT program identifies participant laboratories. The data for blood lead were previously reported in the Blood Lead PT Report issued March 2, 2010, and are reproduced here for completeness.

PT Materials

Test materials for the first event were prepared from caprine (goat) whole blood obtained from animals dosed with lead acetate to create physiologically bound lead (Pb). A total of five blood pools were supplemented with different amounts of arsenic (as As^{3+} , As^{5+} , DMA, MMA and arsenobetaine), cadmium (as Cd^{2+}) and mercury as both inorganic (Hg^{2+}) and as methylmercury (CH_3Hg^+) species. In addition to As, Cd, Pb and Hg, blood pools were supplemented with the trace elements manganese, thallium and tin.


Assignment of Target Values for Trace Elements

Except for blood lead, we will implement robust statistics for assigning target values for all trace element panels. Method specific and additional trace element data will continue to be calculated utilizing traditional statistics. The use of robust statistics for assigning target values for proficiency testing pools is one approach that is acceptable under ISO 13528. In collaboration with other trace element PT scheme organizers, we have conducted an evaluation of robust statistics. As a result of our evaluation, we have elected to introduce this approach in our program.

The next PT event for trace elements in whole blood is scheduled to be mailed April 28th, 2010. Please inform our laboratory staff at (518) 474-4484 if the test materials have not arrived within five days of the scheduled mail out date. **The deadline for reporting results is May 26th, 2010.**

Thank you for your participation in this event.

Yours sincerely,



Patrick J. Parsons, Ph.D.
Section Head, Trace Elements PT Program

New York State Department of Health
Event #1, 2010

Whole Blood Arsenic

Test materials for arsenic were prepared from caprine (goat) whole blood preserved with K₂EDTA anticoagulant. A total of five pools were supplemented with different arsenic species: inorganic As³⁺ and As⁵⁺, monomethylarsonic acid (MMA), dimethylarsinic acid (DMA), and arsenobetaine.

Sample	Arsenic species added
BE10-01	As ³⁺
BE10-02	As ⁵⁺ and MMA
BE10-03	As ³⁺ , MMA, and DMA
BE10-04	As ³⁺ , As ⁵⁺ , MMA, DMA, arsenobetaine
BE10-05	As ⁵⁺ , DMA and arsenobetaine

The Target Value assigned for each PT material is the robust mean of the results reported by all participants in this event. The robust statistics were obtained utilizing algorithms based on those presented in **ISO 13528:2005E** Statistical methods for use in proficiency testing by interlaboratory comparisons. Values for whole blood arsenic range from 20.0 µg/L (0.27 µmol/L) to 53.3 µg/L (0.71 µmol/L).

Acceptable range: The acceptable range for arsenic is set at ±6 µg/L or ±20%, whichever is greater. Thus, it is fixed at ± 6 µg/L for concentrations below 30 µg/L.

Discussion: Based upon the above criteria, 96.8% of test results reported were judged as satisfactory, with one out of 19 participant laboratories (5.3%) reporting 2 or more of the 5 results outside the acceptable ranges.

New York State Department of Health
Blood Arsenic Test Results, 2010 Event #1
ROBUST STATISTICAL SUMMARY

TARGET VALUE ASSIGNMENT AND STATISTICS

Results ($\mu\text{g/L}$ whole blood)

	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
--	---------	---------	---------	---------	---------

Robust Mean	15.2	20.0	31.1	53.3	36.7
Robust Standard Deviation	1.8	2.5	1.3	4.7	2.9
Standard Uncertainty	0.5	0.7	0.4	1.4	0.8
RSD (%)	12.1	12.5	4.1	8.9	7.8
Acceptable Range:					
Upper Limit	21.2	26.0	37.3	64.0	44.0
Lower Limit	9.2	14.0	24.9	42.6	29.4

notes: Results reported as less than the method detection limit are excluded from statistical calculations.

New York State Department of Health
Blood Arsenic Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES

Lab Code	Method	Results (µg/L whole blood)					Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05	
Target Values:		15.2	20.0	31.1	53.3	36.7	
103	DRC/CC-ICP-MS	14.1	18.8	30.3	54.8	36.1	Info
110	DRC/CC-ICP-MS	14.6	19.8	31.1	53.3	38.3	
114	ICP-MS	20	25	36	59	42	
147	ICP-MS	13.4	16.8	28.3	50.4	33.4	Info
156	ICP-MS	11.7	16.6	27.7	43	34.6	
159	ICP-MS	16	22	33	57	41	
164	ICP-MS	23 ↑	30 ↑	38 ↑	62	43	
179	ICP-MS	17	21	31	55	37	
197	DRC/CC-ICP-MS	15	19	31	54	36	
200	ICP-MS	15.1	19.4	31.3	53.1	36.4	Info
206	ICP-MS	16.5	24.4	30.7	56.5	39.3	
208	ICP-MS	17.8	21.3	31.7	51.8	36.1	
293	DRC/CC-ICP-MS	15.7	21.1	32.9	60.1	38.8	Info
305	DRC/CC-ICP-MS	12.7	17.7	29.3	51.2	33	
312	DRC/CC-ICP-MS	17	21	31	54	37	
324	DRC/CC-ICP-MS	14.2	22.2	32.1	47.1	36.6	Info
359	ICP-MS	14.6	19.4	30.4	48.5	34.3	
391	DRC/CC-ICP-MS	13.9	16.3	26.9	46.9	31.5	Info
395	DRC/CC-ICP-MS	13.9	18.7	31.6	53.4	36.1	

Percent satisfactory results for all participants: 96.8 %

notes: ↑ reported outside upper limit
↓ reported outside lower limit
▼ Unacceptable result

notes: Results reported as less than the method detection limit are excluded from statistical calculations.
Info only: results included for informational purposes only.

New York State Department of Health
Blood Arsenic Test Results, 2010 Event #1
STATISTICAL SUMMARY BY METHOD

Results ($\mu\text{g/L}$ whole blood)					
	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
DRC/CC-ICP-MS					
Number of Sample Measurements:	9	9	9	9	9
Mean:	14.6	19.4	30.7	52.8	35.9
Standard Deviation:	1.2	1.8	1.8	4.0	2.3
RSD (%):	8.4	9.5	5.7	7.7	6.5
ICP-MS					
Number of Sample Measurements:	10	10	10	10	10
Mean:	16.5	21.6	31.8	53.6	37.7
Standard Deviation:	3.2	4.1	3.2	5.5	3.4
RSD (%):	19.7	18.8	10.0	10.3	9.0
All Laboratories					
Number of Sample Measurements:	19	19	19	19	19
Mean:	15.6	20.6	31.3	53.2	36.9
Standard Deviation:	2.6	3.3	2.6	4.8	3.0
RSD (%):	16.9	16.1	8.3	8.9	8.2

notes: ? Insufficient data for calculation.

New York State Department of Health
Event #1, 2010

Whole Blood Cadmium

Test materials for cadmium were prepared from caprine (goat) whole blood preserved with K₂EDTA anticoagulant. A total of five blood pools were supplemented with different amounts of cadmium (as Cd²⁺).

The Target Value assigned for each PT material is the robust mean of the results reported by all participants in this event. The robust statistics were obtained utilizing algorithms based on those presented in **ISO 13528:2005E** Statistical methods for use in proficiency testing by interlaboratory comparisons. Values for whole blood cadmium range from 1.9 µg/L (17 nmol/L) to 10.0 µg/L (89 nmol/L).

Acceptable ranges are based on the OSHA criteria of ±15%, or ±1 µg/L around the target value, whichever is greater. So, the range is fixed at ±1 µg/L for concentrations below 6.6 µg/L, where above 6.6 µg/L, it is ±15%.

Discussion: Based upon the above criteria, 92.5%^a of the results reported by all participants were satisfactory, with three out of 32 laboratories (9.4%) reporting 2 or more of the 5 results outside the acceptable ranges.

^a Amended report issued 3/24/10. The number of participating laboratories that reported results was increased to 32 from 30. The percent of results reported as satisfactory increased to 92.5% from 92.0%. There was no change in the number of laboratories reporting 2 or more of the 5 results outside the acceptable ranges. No changes were made to either robust statistical summary or target value assignment.

New York State Department of Health
Blood Cadmium Test Results, 2010 Event #1
ROBUST STATISTICAL SUMMARY

TARGET VALUE ASSIGNMENT AND STATISTICS

Results ($\mu\text{g/L}$ whole blood)

	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
--	---------	---------	---------	---------	---------

Robust Mean	10.0	1.9	2.7	7.2	4.0
Robust Standard Deviation	0.7	0.2	0.2	0.5	0.3
Standard Uncertainty	0.2	0.0	0.1	0.1	0.1
RSD (%)	6.8	8.4	9.1	6.9	8.2
Acceptable Range:					
Upper Limit	11.5	2.9	3.7	8.3	5.0
Lower Limit	8.5	0.9	1.7	6.1	3.0

notes: Results reported as less than the method detection limit are excluded from statistical calculations.

New York State Department of Health
Blood Cadmium Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES

Lab Code	Method	Results (µg/L whole blood)					Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05	
Target Values:		10.0	1.9	2.7	7.2	4.0	
103	DRC/CC-ICP-MS	10.5	2.0	2.9	7.7	4.4	Info
106	ICP-MS	10.4	1.9	2.8	7.5	4.2	Info
107	ICP-MS	9.6	1.7	2.7	6.9	4.0	Info
109	ICP-MS	10.6	2.1	3.0	7.6	4.4	Info
110	ICP-MS	10.3	1.8	2.9	7.5	4.3	
114	ICP-MS	10	2.1	2.8	7.4	4.2	
116	ICP-MS	10.8	1.8	2.8	7.5	4.2	Info
147	ICP-MS	10.2	1.8	2.7	6.9	4.0	Info
156	ICP-MS	9.2	1.7	2.5	11.3 ↑	3.9	
159	ICP-MS	10.2	2	2.9	7.2	4.4	
164	ICP-MS	9.7	1.6	2.6	7.8	3.9	
179	ICP-MS	9.5	2	2.5	6.7	3.8	
197	DRC/CC-ICP-MS	10.3	1.8	2.6	7.2	4	
200	ICP-MS	11.7 ↑	1.9	3.1	8.9 ↑	4.9	Info
206	ICP-MS	10.6	1.9	2.9	7.5	4.4	
208	ICP-MS	8.6	1.5	2.3	6 ↓	3.3	
293	ICP-MS	10.1	1.9	2.7	7.2	4.2	Info
305	ICP-MS	9.9	1.6	2.6	7.1	3.6	
312	ICP-MS	9.8	2	3	7.2	4.1	
324	ICP-MS	10.5	1.9	2.7	7.5	3.8	Info
325	ETAAS-Z	8.1 ↓	0.8 ↓	1.5 ↓	5.8 ↓	3.4	Info
339	HR-ICP-MS	10.1	1.7	2.7	7.1	4.0	Info
359	ICP-MS	8.7	1.8	2.6	6.2	3.6	
366	ETAAS-Z	9.3	1.9	2.4	6.2	2.8 ↓	Info
367	ETAAS-Z	9.7	1.7	2.5	7.1	3.7	Info
383	ETAAS-Z	8.5	2	2.6	6.3	3.7	
385	ICP-MS	9.7	2.0	2.9	7.2	4.1	Info
391	DRC/CC-ICP-MS	11.3	2.8	4.3 ↑	7.8	5.4 ↑	Info
395	ICP-MS	10	1.9	2.5	7.1	4.2	
401	ETAAS-Z	9.0	2.0	2.0	6.0 ↓	4.0	Info
408	ICP-MS	9.6	1.9	2.6	7.0	3.9	Info
410	ICP-MS	10.7	1.8	2.8	7.8	4.3	Info

Percent satisfactory results for all participants: 92.5 %

notes: ↑ reported outside upper limit
↓ reported outside lower limit
▼ Unacceptable result

notes: Results reported as less than the method detection limit are excluded from statistical calculations.
Info only: results included for informational purposes only.

Amended 3/24/10

New York State Department of Health
Blood Cadmium Test Results, 2010 Event #1
STATISTICAL SUMMARY BY METHOD

	Results ($\mu\text{g/L}$ whole blood)				
	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
DRC/CC-ICP-MS					
Number of Sample Measurements:	3	3	3	3	3
Mean:	10.7	2.2	3.3	7.6	4.6
Standard Deviation:	0.5	0.5	0.9	0.3	0.7
RSD (%):	—	—	—	—	—
ETAAS-Z					
Number of Sample Measurements:	5	5	5	5	5
Mean:	8.9	1.7	2.2	6.3	3.5
Standard Deviation:	0.6	0.5	0.5	0.5	0.5
RSD (%):	7.1	30.2	20.6	7.9	12.9
HR-ICP-MS					
Number of Sample Measurements:	1	1	1	1	1
Mean:	10.1	1.7	2.7	7.1	4.0
Standard Deviation:	?	?	?	?	?
RSD (%):	—	—	—	—	—
ICP-MS					
Number of Sample Measurements:	23	23	23	23	23
Mean:	10.0	1.9	2.7	7.4	4.1
Standard Deviation:	0.7	0.2	0.2	1.0	0.3
RSD (%):	6.9	8.4	7.1	13.7	8.2
All Laboratories					
Number of Sample Measurements:	32	32	32	32	32
Mean:	9.9	1.9	2.7	7.3	4.0
Standard Deviation:	0.8	0.3	0.4	1.0	0.5
RSD (%):	8.0	15.8	15.7	13.5	11.5

notes: ? Insufficient data for calculation.

Whole Blood Mercury

Test materials for mercury were prepared from caprine (goat) whole blood preserved with K₂EDTA anticoagulant. A total of five pools were supplemented with different amounts of mercury as both inorganic (Hg²⁺) and organometallic (as methylmercury, CH₃Hg⁺) species.

The Target Value assigned for each PT material is the robust mean of the results reported by all participants in this event. The robust statistics were obtained utilizing algorithms based on those presented in **ISO 13528:2005E** Statistical methods for use in proficiency testing by interlaboratory comparisons. Values for whole blood mercury range from 1.3 µg/L (6 nmol/L) to 45.1 µg/L (225 nmol/L).

Acceptable ranges were fixed at ±30%, or ±3 µg/L around the target value, whichever is greater. That is, the range is fixed at ±3 µg/L for concentrations below 10 µg/L, while above 10 µg/L, it is ±30%.

Discussion: Based on the above criteria, 96.0% of results reported by all participants were satisfactory, with one of the 30 laboratories (3.3%) reporting 2 or more of the 5 results outside the acceptable ranges. Note: Some methods based on cold vapor generation (e.g., CV-AAS) may only detect inorganic Hg thus leading to a low bias compared to methods based on total Hg measurement (e.g., ICP-MS, and CV-AAS combined with on-line microwave digestion). See Barbosa et al. (2004) JAAS (1) for more details on total Hg in blood using CV-AAS.

(1) Barbosa F, Palmer CD, Krug FJ, Parsons PJ. Determination of total mercury in whole blood by flow injection cold vapor atomic absorption spectrometry with room temperature digestion using tetramethylammonium hydroxide. Journal of Analytical Atomic Spectrometry 2004;19(8):1000-1005.

New York State Department of Health
Blood Mercury Test Results, 2010 Event #1
ROBUST STATISTICAL SUMMARY

TARGET VALUE ASSIGNMENT AND STATISTICS

Results ($\mu\text{g/L}$ whole blood)

	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
--	---------	---------	---------	---------	---------

Robust Mean	3.1	1.3	9.4	45.1	11.8
Robust Standard Deviation	0.4	0.3	0.8	4.0	1.3
Standard Uncertainty	0.1	0.1	0.2	0.9	0.3
RSD (%)	12.5	22.7	8.2	8.8	10.9
Acceptable Range:					
Upper Limit	6.1	4.3	12.4	58.6	15.3
Lower Limit	0.1	0.0	6.4	31.6	8.3

notes: Results reported as less than the method detection limit are excluded from statistical calculations.

New York State Department of Health
Blood Mercury Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES

Lab Code	Method	Results (µg/L whole blood)					Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05	
Target Values:		3.1	1.3	9.4	45.1	11.8	
103	DRC/CC-ICP-MS	3.2	1.4	9.8	47.0	12.5	Info
106	ICP-MS	3.2	1.3	9.8	47.3	12.5	Info
107	ICP-MS	3.3	1.5	10.3	48.0	12.9	Info
109	ICP-MS	3.1	1.3	9.6	45.4	12.1	Info
110	ICP-MS	3.1	1.3	9.8	47.9	11.9	
114	ICP-MS	3.1	1.4	9.4	44.9	11.6	
116	ICP-MS	3.1	1.2	10.0	48.6	12.8	Info
147	ICP-MS	2.8	1.3	9.1	44.1	11.9	Info
156	ICP-MS	<3.0	<3.0	8.9	44	11.6	
159	ICP-MS	3	<2.0	10	44	13	
164	ICP-MS	2.6	1	8.4	42	11	
179	ICP-MS	3	1	9	42	11	
197	DRC/CC-ICP-MS	<5.0	<5.0	11	53	18 ↑	
200	ICP-MS	3.5	1.5	10.5	47.6	12.4	Info
206	ICP-MS	3	<3.0	9	34.5	10	
208	ICP-MS	3.3	1.1	9.4	46.6	10.3	
293	ICP-MS	2.8	11.4 ↑	8.9	0.6 ↓	11.5	Info
305	ICP-MS	6.4 ↑	<2.0	8.8	41.1	10.6	
312	ICP-MS	2.5	0.9	8.7	42	10	
324	CV-AAS	2.8	1.2	8.7	42.3	11.1	Info
339	HR-ICP-MS	2.8	1.1	9.2	44.4	11.5	Info
359	ICP-MS	2.5	1	8	38.2	10.3	
366	ICP-MS	5.3	2.5	11.6	59.0 ↑	14.6	Info
367	CV-AAS	2.8	0.8	9.4	48.6	13.4	Info
385	ICP-MS	3.5	<2.0	9.9	46.1	12.6	Info
391	CV-AAS	6.7 ↑	1.7	10.0	55.7	9.9	Info
395	ICP-MS	2.9	1.3	8.9	44.3	12	
401	CV-AAS	3.0	1.7	7.4	40.8	10.0	Info
408	ICP-MS	3.3	1.5	9.5	44.2	11.7	Info
410	ICP-MS	3.3	1.4	10.6	50.1	13.1	Info

Percent satisfactory results for all participants: 96.0 %

notes: ↑ reported outside upper limit
↓ reported outside lower limit
▼ Unacceptable result

notes: Results reported as less than the method detection limit are excluded from statistical calculations.
Info only: results included for informational purposes only.

New York State Department of Health
Blood Mercury Test Results, 2010 Event #1
STATISTICAL SUMMARY BY METHOD

Results ($\mu\text{g/L}$ whole blood)					
	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
CV-AAS					
Number of Sample Measurements:	3	4	4	4	4
Mean:	2.9	1.4	8.9	46.9	11.1
Standard Deviation:	0.1	0.4	1.1	6.8	1.6
RSD (%):	—	32.3	12.6	14.5	14.7
DRC/CC-ICP-MS					
Number of Sample Measurements:	1	1	2	2	2
Mean:	3.2	1.4	10.4	50.0	15.3
Standard Deviation:	?	?	0.8	4.2	3.9
RSD (%):	—	—	—	—	—
HR-ICP-MS					
Number of Sample Measurements:	1	1	1	1	1
Mean:	2.8	1.1	9.2	44.4	11.5
Standard Deviation:	?	?	?	?	?
RSD (%):	—	—	—	—	—
ICP-MS					
Number of Sample Measurements:	21	16	23	22	23
Mean:	3.2	1.3	9.5	45.1	11.8
Standard Deviation:	0.6	0.2	0.8	4.7	1.1
RSD (%):	18.1	15.7	8.5	10.5	9.7
All Laboratories					
Number of Sample Measurements:	26	22	30	29	30
Mean:	3.1	1.3	9.5	45.6	11.9
Standard Deviation:	0.5	0.2	0.9	4.9	1.6
RSD (%):	16.8	19.1	9.2	10.8	13.7

notes: ? Insufficient data for calculation.

New York State Department of Health
Event #1, 2010

Whole Blood Lead

Test materials for lead were prepared from caprine (goat) whole blood obtained from animals dosed with lead acetate to create physiologically-bound Pb. Whole blood was collected into collection bags containing K₂EDTA anticoagulant.

Target values were established as the mean of 20 measurements performed by 19 reference laboratories using ICP-MS, ETAAS and ASV methods. Values range from 3 µg/dL to 43 µg/dL. Among the reference group, imprecision (SD) varied from 0.6 - 1.9 µg/dL, increasing with Pb concentration.

Acceptable ranges are based on the CLIA '88 criteria (Federal Register Volume 57, Number 40, §§ 493.2 and 493.937, February 28, 1992). The criteria are set at ±10% or ±4 µg/dL, whichever is greater.

Discussion Based on the CLIA '88 criteria, 92.0% of results reported by all participants were judged as satisfactory, with 13 out of 100 participant laboratories (13.0%) reporting 2 or more of the 5 results outside the acceptable ranges.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES**

Lab Code	Method	Results (µg/dL whole blood)					Normalized Mean	Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05		
Target Values:		17	21	43	38	3		
103	ASV-LeadCare	18	21	46	40	3	1.05	Info
103	DRC/CC-ICP-MS	17	20	41	36	3	0.96	
104	ETAAS-Z	18	22	44	39	3	1.04	
106	ICP-MS	17	21	43	38	3	1.00	Info
107	ICP-MS	18	22	44	39	3	1.04	
107	ASV-LeadCare	14	17	37 ↓	30 ↓	<3	0.82	Info
108	ETAAS-Z	16	21	39	36	5	0.95	
109	ETAAS-Z	19	23	46	41	4	1.09	Info
109	ASV-LeadCare	18	20	42	36	3	0.98	
109	ASV-LeadCare	16	19	39	34	<3	0.91	
110	ETAAS-Z	18	22	44	38	4	1.03	Info
110	ICP-MS	18	21	43	38	3	1.01	
110	ASV-LeadCare	18	21	44	38	3	1.02	
110	ASV-LeadCare	18	20	46	36	<3	1.01	Info
112	ASV-3010	15	19	39	34	<2	0.90	
114	ETAAS-Z	19	24	47	40	6	1.10	Info
115	ETAAS-Z	17	20	43	30 ↓	2	0.94	
116	ICP-MS	18	22	45	39	3	1.04	
121	ETAAS-Z	20	24	52 ↑	60 ↑	3	1.28	Info
123	ETAAS-Z	16	20	40	35	3	0.94	
126	ETAAS-Z	17	22	43	39	3	1.02	Info
131	ETAAS-Z	19	22	43	32 ↓	7	1.00	
132	ETAAS-Z	17	20	41	36	3	0.96	
143	ETAAS-Z	16	20	41	36	3	0.95	Info
144	ETAAS-Z	17	21	43	37	3	0.99	
146	ETAAS-Z	16	20	41	35	2	0.94	
147	ICP-MS	17	20	41	36	3	0.96	Info
150	ASV-LeadCare	17	21	43	38	3	1.00	
156	ICP-MS	18	25	56 ↑	44 ↑	3	1.18	
158	ETAAS-Z	18	21	44	38	4	1.02	Info

notes: ↑ reported value outside upper limit
↓ reported value outside lower limit

Normalized mean: The average of each reported result divided by the corresponding target value. It measures bias.
Info only: results included for informational purposes only.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES**

Lab Code	Method	Results (µg/dL whole blood)					Normalized Mean	Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05		
Target Values:		17	21	43	38	3		
159	ICP-MS	17	20	42	36	3	0.97	
160	ETAAS-Z	19	25	56 ↑	51 ↑	4	1.24	
164	ICP-MS	18	22	45	39	3	1.04	
166	ASV-3010	17	20	41	36	2	0.96	
168	ETAAS-Z	18	23	45	39	4	1.06	
170	ETAAS Other	17	21	43	36	3	0.99	
179	ICP-MS	17	22	44	37	3	1.01	
197	ICP-MS	15	19	38 ↓	33 ↓	3	0.88	
198	ETAAS-Z	17	22	43	38	3	1.01	
199	ETAAS-Z	18	22	43	39	3	1.03	
200	ETAAS-Z	17	21	45	38	2	1.01	
204	ASV-3010	17	18	39	37	2	0.93	
206	ICP-MS	17	21	43	36	3	0.99	
208	ETAAS-Z	21	25	52 ↑	43 ↑	<3	1.19	
215	ETAAS-Z	20	23	46	41	6	1.11	
221	ETAAS-Z	19	21	41	36	4	1.00	
232	ASV-3010	18	22	46	40	3	1.06	
237	ETAAS-Z	18	22	45	40	3	1.05	
243	ASV-3010	17	20	42	37	2	0.98	
249	ASV-3010	17	20	40	33 ↓	3	0.94	
254	ETAAS-Z	18	21	42	38	4	1.01	
255	ETAAS-Z	17	22	42	37	3	1.00	
261	ETAAS-Z	16	20	39	34	3	0.92	
269	ETAAS-Z	18	20	42	36	3	0.98	
271	ASV-3010	15	18	40	35	1	0.90	
272	ETAAS-Z	17	21	44	37	3	1.00	
279	ETAAS-Z	15	18	33 ↓	28 ↓	3	0.81	
282	ASV-3010	13	15 ↓	34 ↓	30 ↓	1	0.76	
286	ASV-LeadCare	16	18	39	31 ↓	<3	0.88	
290	ICP-MS	18	21	43	38	3	1.01	

notes: ↑ reported value outside upper limit
↓ reported value outside lower limit

Normalized mean: The average of each reported result divided by the corresponding target value. It measures bias.
Info only: results included for informational purposes only.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES**

Lab Code	Method	Results (µg/dL whole blood)					Normalized Mean	Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05		
Target Values:		17	21	43	38	3		
291	ASV-3010	14	16 ↓	28 ↓	26 ↓	17 ↑	1.72	
293	ICP-MS	17	22	44	38	3	1.02	
295	ASV-3010	15	19	40	35	3	0.91	
300	ASV-3010	18	21	44	38	<3	1.02	
301	ETAAS-Z	17	21	42	36	3	0.98	
305	ETAAS-Z	17	21	41	36	3	0.98	
312	ICP-MS	17	21	44	38	3	1.01	
317	ETAAS-Z	19	23	47	40	4	1.09	
324	ICP-MS	18	22	42	39	3	1.03	
325	ETAAS-Z	15	18	40	35	3	0.90	
333	ETAAS-Z	18	22	43	38	4	1.03	
339	HR-ICP-MS	16	19	40	35	3	0.92	Info
340	ETAAS-Z	17	20	40	35	3	0.95	
343	ASV-LeadCare	17	18	42	35	3	0.94	Info
348	ETAAS-Z	17	21	41	36	3	0.98	
349	ETAAS-Z	17	22	43	39	3	1.02	
350	ASV-3010	18	23	46	41	4	1.08	
352	ASV-3010	18	23	46	41	4	1.08	
353	ETAAS-Z	16	20	43	37	<2	0.97	
359	ICP-MS	15	18	38 ↓	31 ↓	3	0.86	
365	ETAAS-Z	19	24	49 ↑	42	3	1.13	
366	ETAAS-Z	17	23	46	41	3	1.06	Info
367	ETAAS-Z	17	21	42	39	3	1.00	Info
368	ASV-3010	15	16 ↓	40	31 ↓	2	0.85	
369	ASV-3010	17	19	40	33 ↓	3	0.93	
374	ASV-3010	15	18	40	34	2	0.89	
376	ASV-LeadCare	14	20	40	33 ↓	3	0.89	
383	ETAAS-Z	17	20	42	35	3	0.96	
384	ASV-3010	17	20	42	38	4	0.98	
385	ICP-MS	17	21	41	35	4	0.97	Info

notes: ↑ reported value outside upper limit
↓ reported value outside lower limit

Normalized mean: The average of each reported result divided by the corresponding target value. It measures bias.
Info only: results included for informational purposes only.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
PERFORMANCE OF PARTICIPATING LABORATORIES**

Lab Code	Method	Results (µg/dL whole blood)					Normalized Mean	Info Only
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05		
Target Values:		17	21	43	38	3		
388	ASV-LeadCare	16	22	46	34	2	0.99	
389	ETAAS-Z	18	21	43	38	3	1.01	
391	ETAAS Other	16	20	40	32 ↓	6	0.92	Info
395	ICP-MS	17	21	43	38	3	1.00	
401	ETAAS-Z	19	23	46	39	4	1.08	Info
408	ICP-MS	16	20	40	35	3	0.94	Info
410	ICP-MS	18	22	45	39	3	1.04	Info
449	ASV-LeadCare	17	16 ↓	47	37	3	0.96	
455	ASV-LeadCare	14	15 ↓	35 ↓	33 ↓	<3	0.81	
456	ASV-LeadCare	13	16 ↓	37 ↓	32 ↓	3	0.81	
Percent satisfactory results for all participants:							92.0 %	

notes: ↑ reported value outside upper limit
↓ reported value outside lower limit

Normalized mean: The average of each reported result divided by the corresponding target value. It measures bias.
Info only: results included for informational purposes only.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
STATISTICAL SUMMARY**

TARGET VALUE ASSIGNMENT AND STATISTICS						
Lab Code	Method	Results (µg/dL whole blood)				
		BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
103	DRC/CC-ICP-MS	17	20	41	36	3
104	ETAAS-Z	18	22	44	39	3
107	ICP-MS	18	22	44	39	3
109	ETAAS-Z	19	23	46	41	4
110	ETAAS-Z	18	22	44	38	4
110	ICP-MS	18	21	43	38	3
112	ASV-3010	15	19	39	34	<2
147	ICP-MS	17	20	41	36	3
159	ICP-MS	17	20	42	36	3
164	ICP-MS	18	22	45	39	3
166	ASV-3010	17	20	41	36	2
179	ICP-MS	17	22	44	37	3
198	ETAAS-Z	17	22	43	38	3
199	ETAAS-Z	18	22	43	39	3
200	ETAAS-Z	17	21	45	38	2
243	ASV-3010	17	20	42	37	2
293	ICP-MS	17	22	44	38	3
324	ICP-MS	18	22	42	39	3
325	ETAAS-Z	15	18	40	35	3
350	ASV-3010	18	23	46	41	4
Number of Sample Measurements:		20	20	20	20	19
Mean (target value):		17	21	43	38	3
Standard Deviation:		1.0	1.3	1.9	1.8	0.6
RSD (%):		5.7	6.4	4.5	4.9	19.2
Acceptable Range:						
Upper Limit:		21	25	47	42	7
Lower Limit:		13	17	39	34	0

notes: Results reported as less than the detection limits are treated as zero for statistical and grading purposes.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
STATISTICAL SUMMARY BY METHOD**

	Results (µg/dL whole blood)				
	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
ASV-3010					
Number of Sample Measurements:	17	17	17	17	14
Mean:	16.2	19.2	40.4	35.2	2.6
Standard Deviation:	1.6	2.3	4.4	4.0	1.0
RSD (%):	9.6	12.0	10.9	11.3	39.5
ASV-LeadCare					
Number of Sample Measurements:	14	14	14	14	9
Mean:	16.1	18.9	41.6	34.8	2.9
Standard Deviation:	1.7	2.2	3.9	2.9	0.3
RSD (%):	10.8	11.7	9.4	8.3	11.5
DRC/CC-ICP-MS					
Number of Sample Measurements:	1	1	1	1	1
Mean:	17.0	20.0	41.0	36.0	3.0
Standard Deviation:	?	?	?	?	?
RSD (%):	—	—	—	—	—
ETAAS Other					
Number of Sample Measurements:	2	2	2	2	1
Mean:	16.5	20.5	41.5	34.0	3.0
Standard Deviation:	0.7	0.7	2.1	2.8	?
RSD (%):	—	—	—	—	—
ETAAS-Z					
Number of Sample Measurements:	45	45	45	45	40
Mean:	17.6	21.5	43.5	38.1	3.2
Standard Deviation:	1.3	1.6	3.8	4.9	0.6
RSD (%):	7.4	7.3	8.7	12.8	19.2
HR-ICP-MS					
Number of Sample Measurements:	1	1	1	1	1
Mean:	16.0	19.0	40.0	35.0	3.0
Standard Deviation:	?	?	?	?	?
RSD (%):	—	—	—	—	—
ICP-MS					
Number of Sample Measurements:	20	20	20	20	20
Mean:	17.2	21.2	43.2	37.3	3.1
Standard Deviation:	0.9	1.4	3.7	2.7	0.2
RSD (%):	5.4	6.7	8.5	7.2	7.3
All Laboratories					
Number of Sample Measurements:	100	100	100	100	86
Mean:	17.0	20.6	42.6	36.8	3.0
Standard Deviation:	1.4	2.1	4.0	4.2	0.6
RSD (%):	8.5	10.0	9.3	11.4	21.1

notes: ? Insufficient data for calculation.

**New York State Department of Health
Blood Lead Test Results, 2010 Event #1
STATISTICAL SUMMARY BY CLASS**

	Results ($\mu\text{g/dL}$ whole blood)				
	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
Evaluated					
Number of Sample Measurements:	62	62	62	62	53
Mean:	16.9	20.4	42.3	36.3	3.0
Standard Deviation:	1.6	2.3	4.5	4.0	0.7
RSD (%):	9.4	11.2	10.7	10.9	23.9
Info					
Number of Sample Measurements:	18	18	18	18	14
Mean:	17.2	20.7	43.1	37.8	3.1
Standard Deviation:	1.4	1.8	3.5	6.2	0.4
RSD (%):	7.9	8.6	8.2	16.5	11.6
Reference					
Number of Sample Measurements:	20	20	20	20	19
Mean:	17.3	21.2	43.0	37.7	3.0
Standard Deviation:	1.0	1.3	1.9	1.8	0.6
RSD (%):	5.7	6.4	4.5	4.9	19.2
All Laboratories					
Number of Sample Measurements:	100	100	100	100	86
Mean:	17.0	20.6	42.6	36.8	3.0
Standard Deviation:	1.4	2.1	4.0	4.2	0.6
RSD (%):	8.5	10.0	9.3	11.4	21.1

notes: ? Insufficient data for calculation.

Additional Trace Elements Reported in Whole Blood

For this interlaboratory study, we requested that participant laboratories report their analytical results for any additional trace elements (other than As, Cd, Pb and Hg) that are routinely reported so that a more complete characterization can be recorded for these materials. Results for additional trace elements are reported here, but no target value is implied nor are any acceptable ranges provided. These data are provided solely for educational and informational purposes.

In addition to As, Cd, Pb and Hg, blood pools were supplemented with additional trace elements as indicated below.

Additional Elements

Mn, Sn, Tl

New York State Department of Health
Whole Blood Additional Elements, 2010 Event #1
Page 1

Blood Barium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	9.2	16.2	11.4	14.9	13.6
147	ICP-MS	9.54	15.65	11.51	14.83	12.65
197	ICP-MS	9.7	17.5	12.6	15.5	13.4
312	ICP-MS	9.0	15.4	11.5	15.4	13.0
Mean, n=4		9.4	16.2	11.8	15.2	13.2
SD		0.3	0.9	0.6	0.3	0.4

Blood Beryllium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	< 0.14	0.96	0.65	< 0.14	0.44
197	ICP-MS	< 0.2	< 0.2	< 0.2	< 0.2	< 0.2

Blood Cobalt Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	0.46	0.47	0.47	< 0.36	< 0.36
147	ICP-MS	0.5	0.43	0.53	0.39	0.47
159	ICP-MS	< 1.0	< 1.0	< 1.0	< 1.0	< 1.0
197	ICP-MS	< 1.0	< 1.0	< 1.0	< 1.0	< 1.0
312	ICP-MS	0.5	0.4	0.5	0.4	0.4
391	DRC/CC-ICP-MS	0.4	0.3	0.4	0.3	0.3
Mean, n=6		0.47	0.43	0.50	0.40	0.39
SD		0.05	0.09	0.09	0.10	0.09

Blood Chromium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
197	DRC/CC-ICP-MS	< 1.0	< 1.0	< 1.0	< 1.0	< 1.0
312	DRC/CC-ICP-MS	0.5	0.5	0.6	0.5	0.4

Blood Cesium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	0.56	0.54	0.71	0.52	0.46

Blood Copper Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	1060	1072	1097	1059	1013
147	ICP-MS	1054.64	1080.05	1124.52	1099.11	1054.64
197	ICP-MS	104	102	114	99	94
312	ICP-MS	1030	1060	1130	1070	1000
Mean, n=4		1048	1071	1117	1076	1023
SD		16	10	18	21	29

New York State Department of Health
Whole Blood Additional Elements, 2010 Event #1
Page 2

Blood Manganese Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
107	DRC/CC-ICP-MS	51.6	53.6	26.1	24.3	38.5
110	ETAAS-Z	46.4	48.8	23.4	21.8	34.0
147	ICP-MS	46.1	49.18	24.84	24.51	35.44
159	ICP-MS	55	59	32	29	45
179	ETAAS-Z	49.4	51.1	25.1	22.7	37.2
197	DRC/CC-ICP-MS	53.0	52.6	27.2	25.0	39.7
293	ICP-MS	44.1	46	25.3	17.6	31.7
305	ICP-MS	48.5	47.8	24.9	21.4	33.3
312	DRC/CC-ICP-MS	51	53	27	25	38
391	DRC/CC-ICP-MS	41.6	41.7	20.5	20.2	30.1
Mean, n=10		48.7	50.3	25.6	23.2	36.3
SD		4.2	4.7	2.9	3.1	4.4

Blood Molybdenum Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
147	ICP-MS	21.79	23.13	26.97	27.26	17.47

Blood Nickel Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
147	ICP-MS	0.86	0.79	0.91	0.93	0.7

Blood Antimony Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	< 0.10	0.12	< 0.10	< 0.10	< 0.10

Blood Selenium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
147	ICP-MS	401.26	409.16	328.59	308.85	331.75
197	ICP-MS	428	394	326	295	341
305	ICP-MS	465	441	380	313	406
312	ICP-MS	452	435	374	345	396
391	DRC/CC-ICP-MS	450.5	423.3	354.8	327.9	373.3
Mean, n=5		438	420	353	318	370
SD		25	19	25	19	33

Blood Tin Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	9.9	1.6	2.5	7.5	3.9
147	ICP-MS	9.64	1.22	2.1	6.76	3.59
197	ICP-MS	9.5	< 5.0	< 5.0	7.3	< 5.0
Mean, n=3		9.7	1.4	2.3	7.2	3.7
SD		0.2	0.3	0.3	0.4	0.2

New York State Department of Health
Whole Blood Additional Elements, 2010 Event #1
Page 3

Blood Thallium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	0.79	1.46	2.34	4.28	3.22
147	ICP-MS	0.7	1.35	1.97	3.9	2.92
159	ICP-MS	< 1.0	1.4	2.1	4.4	3.2
197	ICP-MS	< 1.0	1.2	1.8	3.8	2.7
312	ICP-MS	0.7	1.4	2.2	4.3	3.1
Mean, n=5		0.7	1.4	2.1	4.1	3.0
SD		0.1	0.1	0.2	0.3	0.2

Blood Uranium Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	< 0.02	< 0.02	< 0.02	< 0.02	< 0.02
312	ICP-MS	< 0.1	< 0.1	< 0.1	< 0.1	< 0.1

Blood Zinc Results (µg/L)						
Lab Code	Method	BE10-01	BE10-02	BE10-03	BE10-04	BE10-05
110	ICP-MS	1717	2082	1967	1688	3379
147	ICP-MS	1712.42	2294.12	2169.93	2039.22	3457.52
312	ICP-MS	1850	2250	2290	1990	3730
Mean, n=3		1760	2209	2142	1906	3522
SD		78	112	163	190	184

New York State Department of Health
Trace Elements in Whole Blood
METHOD NOTES

ATOMIC SPECTROMETRY METHODS

- A-1 ETAAS-Z (Electrothermal atomic absorption spectrometry with Zeeman background correction)
- A-2 ETAAS other (i.e., D₂, S-H background correction)
- A-3 FAAS (Flame atomic absorption spectrometry)
- A-4 CV-AAS (Cold vapor atomic absorption spectrometry)
- A-5 HG-AAS (Hydride generation atomic absorption spectrometry)
- A-6 AFS (Atomic fluorescence spectrometry)
- A-7 Other

INDUCTIVELY COUPLED PLASMA

- P-1 ICP-MS (Inductively coupled plasma - mass spectrometry)
- P-2 DRC/CC-ICP-MS (ICP-MS used in the Dynamic Reaction Cell or Collision Cell mode)
- P-3 ICP-AES/OES (ICP atomic/optical emission spectrometry)
- P-4 HR-ICP-MS (High resolution ICP-MS)
- P-5 ETV-ICP-MS (Electrothermal vaporization ICP-MS)
- P-6 ID-ICP-MS (Isotope dilution ICP-MS)
- P-7 Other

ELECTROCHEMICAL METHODS

- E-1 ASV (Anodic stripping voltammetry without digestion)
- E-2 ASV-LeadCare® (Anodic stripping voltammetry using the ESA LeadCare® system)
- E-3 Fluoride specific electrode
- E-4 Other

MOLECULAR FLUORIMETRY

- F-1 EtOAc (Ethyl acetate-acetic acid extraction method for determination of erythrocyte protoporphyrin)
- F-2 Aviv hematofluorometry (for determination of EP at hematocrit 35)
- F-3 Helena ZPP (for determination of zinc protoporphyrin in $\mu\text{mol ZPP/mol heme}$)
- F-4 Other

OTHER METHODS

If your method is not listed in the above list, please describe it briefly.
