
Howard A. Zucker, M.D., J.D.
Acting Commissioner of Health

NEW YORK
state department of
HEALTH

Sally Dreslin, M.S., R.N.
Executive Deputy Commissioner

Statistical Report
Urinalysis/Urine Pregnancy Testing
Mail out date: October 27, 2014

This report summarizes data from the urinalysis/urine pregnancy proficiency testing of October 27, 2014.

If you have any questions regarding this report, please contact Regina Saulsbery at 518-474-6156.
You may also contact us by e-mail: clinchem@wadsworth.org

Results are shown for each analyte arranged by instrument and reagent system. Results for specific gravity and pH are listed as the mean $\pm$ 1SD for each of the instruments and reagents. All other results are shown as the number of participant results for that instrument or reagent.

The use of brand and/or trade names in this report does not constitute an endorsement of the products on the part of the Wadsworth Center or the New York State Department of Health.

PLEASE NOTE: The evaluation report and summary of the educational quantitative urine clinical chemistry proficiency test will be mailed separately if your laboratory participated in that proficiency test.

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Specific Gravity

Specimen: U26	Specimen: U27	Specimen: U28	Specimen: U29	Specimen: U30	Number	Instrument or Reagent System
-----	-----	-----	-----	-----	-----	-----
1.020 ± 0.004	1.017 ± 0.007	1.017 ± 0.005	1.017 ± 0.004	1.014 ± 0.005	n = 359	All Methods & Instruments
1.020	1.015	1.016	1.017	1.012	n = 359	Median value, All Methods & Instruments <Instruments>
1.018 ± 0.003	1.018 ± 0.003	1.012 ± 0.003	1.013 ± 0.003	1.009 ± 0.003	n = 18	Visual examination
1.019 ± 0.002	1.013 ± 0.003	1.016 ± 0.002	1.018 ± 0.002	1.012 ± 0.002	n = 96	Siemens Atlas
1.025 ± 0.002	1.023 ± 0.003	1.024 ± 0.002	1.020 ± 0.001	1.021 ± 0.002	n = 42	Siemens Clinitek Status
1.021 ± 0.002	1.023 ± 0.004	1.020 ± 0.000	1.014 ± 0.002	1.015 ± 0.000	n = 5	Siemens Clinitek 50
1.024 ± 0.003	1.024 ± 0.003	1.024 ± 0.003	1.023 ± 0.003	1.023 ± 0.003	n = 4	Siemens Clinitek 100
1.023 ± 0.003	1.024 ± 0.002	1.021 ± 0.001	1.019 ± 0.002	1.018 ± 0.003	n = 29	Siemens Clinitek 500
1.023 ± 0.003	1.023 ± 0.003	1.020 ± 0.003	1.018 ± 0.003	1.018 ± 0.003	n = 61	Siemens Clinitek Advantus
1.014 ± 0.002	1.016 ± 0.002	1.011 ± 0.002	1.013 ± 0.003	1.008 ± 0.003	n = 6	Roche Chemstrip Criterion
1.016 ± 0.002	1.016 ± 0.004	1.012 ± 0.003	1.015 ± 0.002	1.011 ± 0.002	n = 28	Roche Urisys
1.015 ± 0.000	1.018 ± 0.003	1.010 ± 0.000	1.013 ± 0.003	1.014 ± 0.002	n = 5	Roche Cobas U411
1.015 ± 0.002	1.007 ± 0.002	1.012 ± 0.002	1.013 ± 0.002	1.009 ± 0.002	n = 25	IRIS AUTION MAX
1.016 ± 0.001	1.007 ± 0.000	1.013 ± 0.001	1.012 ± 0.001	1.010 ± 0.001	n = 25	IRIS iChem Velocity
1.017 ± 0.002	1.010 ± 0.002	1.014 ± 0.002	1.016 ± 0.002	1.011 ± 0.001	n = 15	Refractometer
						<Reagents>
1.023 ± 0.004	1.022 ± 0.004	1.020 ± 0.004	1.018 ± 0.004	1.018 ± 0.004	n = 165	Siemens Reagent Strips
1.019 ± 0.002	1.013 ± 0.003	1.016 ± 0.002	1.018 ± 0.002	1.012 ± 0.002	n = 87	Siemens Atlas Reagent
1.015 ± 0.002	1.018 ± 0.003	1.011 ± 0.002	1.013 ± 0.002	1.011 ± 0.003	n = 31	Roche Chemstrips
1.017 ± 0.002	1.011 ± 0.001	1.014 ± 0.002	1.016 ± 0.002	1.011 ± 0.002	n = 11	Roche Urisys Cassette
1.015 ± 0.002	1.007 ± 0.002	1.012 ± 0.002	1.013 ± 0.002	1.009 ± 0.002	n = 16	ARKRAY Aution Sticks
1.016 ± 0.001	1.007 ± 0.001	1.013 ± 0.000	1.012 ± 0.001	1.010 ± 0.001	n = 30	IRIS Reagent Strips
1.017 ± 0.002	1.010 ± 0.003	1.014 ± 0.003	1.015 ± 0.002	1.011 ± 0.002	n = 18	None, no reagent or strip used

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

pH

Specimen: U26	Specimen: U27	Specimen: U28	Specimen: U29	Specimen: U30	Number	Instrument or Reagent System
6.95 ± 0.17	5.98 ± 0.18	7.50 ± 0.23	6.92 ± 0.22	7.95 ± 0.33	n = 364	All Methods & Instruments
7.00	6.00	7.50	7.00	8.00	n = 364	Median value, All Methods & Instruments <Instruments>
6.84 ± 0.35	5.96 ± 0.43	7.58 ± 0.19	6.90 ± 0.25	8.04 ± 0.29	n = 25	Visual examination
7.00 ± 0.00	6.01 ± 0.05	7.49 ± 0.05	6.99 ± 0.05	7.86 ± 0.23	n = 96	Siemens Atlas
6.99 ± 0.07	6.01 ± 0.19	7.50 ± 0.00	7.00 ± 0.00	8.21 ± 0.41	n = 47	Siemens Clinitek Status
6.90 ± 0.22	6.00 ± 0.00	7.40 ± 0.22	6.90 ± 0.22	7.70 ± 0.45	n = 5	Siemens Clinitek 50
7.00 ± 0.00	6.00 ± 0.00	7.60 ± 0.22	7.00 ± 0.00	8.20 ± 0.27	n = 5	Siemens Clinitek 100
6.92 ± 0.19	6.00 ± 0.00	7.48 ± 0.09	6.83 ± 0.27	7.85 ± 0.30	n = 30	Siemens Clinitek 500
6.92 ± 0.18	5.98 ± 0.14	7.49 ± 0.11	6.87 ± 0.30	7.89 ± 0.30	n = 63	Siemens Clinitek Advantus
7.00 ± 0.00	6.00 ± 0.00	7.88 ± 0.35	6.94 ± 0.18	8.63 ± 0.52	n = 8	Roche Chemstrip Criterion
6.98 ± 0.09	5.86 ± 0.36	7.54 ± 0.51	7.02 ± 0.21	7.93 ± 0.26	n = 28	Roche Urisys
6.90 ± 0.22	6.00 ± 0.00	7.20 ± 0.45	6.80 ± 0.27	8.00 ± 0.00	n = 5	Roche Cobas U411
6.78 ± 0.25	6.00 ± 0.00	7.48 ± 0.10	6.52 ± 0.10	7.68 ± 0.24	n = 25	IRIS AUTION MAX
7.00 ± 0.00	6.00 ± 0.00	7.38 ± 0.50	7.00 ± 0.00	8.00 ± 0.00	n = 26	IRIS iChem Velocity <Reagents>
6.94 ± 0.18	6.01 ± 0.15	7.50 ± 0.10	6.91 ± 0.24	7.98 ± 0.35	n = 180	Siemens Reagent Strips
7.00 ± 0.00	6.01 ± 0.05	7.49 ± 0.05	6.99 ± 0.05	7.85 ± 0.23	n = 85	Siemens Atlas Reagent
6.93 ± 0.21	5.81 ± 0.39	7.69 ± 0.46	6.94 ± 0.26	8.13 ± 0.48	n = 35	Roche Chemstrips
7.00 ± 0.00	5.91 ± 0.30	7.27 ± 0.47	7.00 ± 0.00	8.00 ± 0.00	n = 11	Roche Urisys Cassette
6.75 ± 0.26	6.00 ± 0.00	7.53 ± 0.11	6.53 ± 0.11	7.70 ± 0.25	n = 20	ARKRAY Aution Sticks
6.98 ± 0.09	6.00 ± 0.00	7.37 ± 0.45	6.92 ± 0.19	7.94 ± 0.17	n = 31	IRIS Reagent Strips

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Bilirubin	Specimen: U26	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>						
Visual examination		2	2	9	13	55
Siemens Atlas		1	2	87	0	6
Siemens Clinitek Status		0	0	35	10	0
Siemens Clinitek 50		0	0	4	1	0
Siemens Clinitek 100		0	0	1	4	0
Siemens Clinitek 500		0	0	30	0	0
Siemens Clinitek Advantus		0	0	61	1	1
Roche Chemstrip Criterion		0	7	1	0	0
Roche Urisys		0	25	2	0	0
Roche Cobas U411		0	5	0	0	0
IRIS AUTION MAX		1	1	23	0	0
IRIS iChem Velocity		26	0	0	0	0
<Reagents>						
Siemens Reagent Strips		3	2	142	26	3
Siemens Atlas Reagent		0	1	76	0	6
Roche Chemstrips		0	28	4	2	0
Roche Urisys Cassette		0	10	0	0	0
ARKRAY Aution Sticks		1	1	17	0	0
IRIS Reagent Strips		26	0	6	0	5
Siemens Ictotest		0	0	4	1	44
Siemens Hemastix		0	0	1	0	0
Siemens Uristix		0	0	2	0	0
Biorex Ictero-Check		0	0	1	0	4
All Instruments and Reagents		30	42	253	29	62

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Bilirubin	Specimen: U27	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>						
Visual examination	31	0	0	0	0	0
Siemens Atlas	96	0	0	0	0	0
Siemens Clinitek Status	44	0	1	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	29	0	0	0	0	0
Siemens Clinitek Advantus	63	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	27	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	174	0	1	0	0	0
Siemens Atlas Reagent	83	0	0	0	0	0
Roche Chemstrips	33	0	0	0	0	0
Roche Urisys Cassette	10	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0
IRIS Reagent Strips	33	0	0	0	0	0
Siemens Ictotest	9	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0	0
All Instruments and Reagents	364	0	1	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Bilirubin	Specimen: U28	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>						
Visual examination	31	0	0	0	0	0
Siemens Atlas	96	0	0	0	0	0
Siemens Clinitek Status	45	0	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	29	0	0	0	0	0
Siemens Clinitek Advantus	63	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	27	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	175	0	0	0	0	0
Siemens Atlas Reagent	83	0	0	0	0	0
Roche Chemstrips	33	0	0	0	0	0
Roche Urisys Cassette	10	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0
IRIS Reagent Strips	33	0	0	0	0	0
Siemens Ictotest	9	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0	0
All Instruments and Reagents	365	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Bilirubin	Specimen: U29	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>						
Visual examination	31	0	0	0	0	0
Siemens Atlas	96	0	0	0	0	0
Siemens Clinitek Status	45	0	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	29	0	0	0	0	0
Siemens Clinitek Advantus	63	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	27	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	175	0	0	0	0	0
Siemens Atlas Reagent	83	0	0	0	0	0
Roche Chemstrips	33	0	0	0	0	0
Roche Urisys Cassette	10	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0
IRIS Reagent Strips	33	0	0	0	0	0
Siemens Ictotest	9	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0	0
All Instruments and Reagents	365	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Bilirubin	Specimen: U30	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>						
Visual examination	31	0	0	0	0	0
Siemens Atlas	96	0	0	0	0	0
Siemens Clinitek Status	45	0	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	29	0	0	0	0	0
Siemens Clinitek Advantus	63	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	27	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	175	0	0	0	0	0
Siemens Atlas Reagent	83	0	0	0	0	0
Roche Chemstrips	33	0	0	0	0	0
Roche Urisys Cassette	10	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0
IRIS Reagent Strips	33	0	0	0	0	0
Siemens Ictotest	9	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0	0
All Instruments and Reagents	365	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Glucose	Specimen: U26	Negative	Trace 50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination	72		0	0	0	0	0	0
Siemens Atlas	97		0	0	0	0	0	0
Siemens Clinitek Status	47		0	0	0	0	0	0
Siemens Clinitek 50	5		0	0	0	0	0	0
Siemens Clinitek 100	5		0	0	0	0	0	0
Siemens Clinitek 500	30		0	0	0	0	0	0
Siemens Clinitek Advantus	63		0	0	0	0	0	0
Roche Chemstrip Criterion	8		0	0	0	0	0	0
Roche Urisys	28		0	0	0	0	0	0
Roche Cobas U411	5		0	0	0	0	0	0
IRIS AUTION MAX	25		0	0	0	0	0	0
IRIS iChem Velocity	26		0	0	0	0	0	0
Other Chemistry Analyzer	0		0	0	0	0	0	0
Other	1		0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips	182		0	0	0	0	0	0
Siemens Atlas Reagent	85		0	0	0	0	0	0
Roche Chemstrips	34		0	0	0	0	0	0
Roche Urisys Cassette	11		0	0	0	0	0	0
ARKRAY Aution Sticks	20		0	0	0	0	0	0
IRIS Reagent Strips	31		0	0	0	0	0	0
Bayer Clinitest	46		0	0	0	0	0	0
Bayer Diastix	0		0	0	0	0	0	0
Siemens Uristix	3		0	0	0	0	0	0
All Instruments and Reagents	412		0	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Glucose	Specimen: U27	Negative	Trace 50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination		0	0	1	32	29	32	6
Siemens Atlas		0	0	0	77	2	18	0
Siemens Clinitek Status		0	0	2	41	0	4	0
Siemens Clinitek 50		0	0	0	1	0	4	0
Siemens Clinitek 100		0	0	0	1	0	4	0
Siemens Clinitek 500		1	0	0	15	0	16	0
Siemens Clinitek Advantus		0	2	0	27	0	34	0
Roche Chemstrip Criterion		0	0	0	1	0	7	0
Roche Urisys		0	0	0	1	0	27	0
Roche Cobas U411		0	0	0	0	1	4	0
IRIS AUTION MAX		0	0	1	23	1	0	0
IRIS iChem Velocity		1	0	0	24	1	0	0
Other Chemistry Analyzer		0	0	0	0	1	0	0
Other		0	0	0	1	0	0	0
<Reagents>								
Siemens Reagent Strips		1	2	3	104	1	70	3
Siemens Atlas Reagent		0	0	0	66	1	18	0
Roche Chemstrips		0	0	0	4	1	32	0
Roche Urisys Cassette		0	0	0	0	0	11	0
ARKRAY Aution Sticks		0	0	0	19	1	0	0
IRIS Reagent Strips		1	0	1	28	1	0	0
Bayer Clinitest		0	0	0	19	30	18	3
Bayer Diastix		0	0	0	1	0	0	0
Siemens Uristix		0	0	0	3	0	1	0
All Instruments and Reagents		2	2	4	244	35	150	6

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Glucose	Specimen: U28	Negative	Trace 50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination	65	9	0	0	0	0	0	0
Siemens Atlas	89	8	0	0	0	0	0	0
Siemens Clinitek Status	38	9	0	0	0	0	0	0
Siemens Clinitek 50	1	4	0	0	0	0	0	0
Siemens Clinitek 100	3	2	0	0	0	0	0	0
Siemens Clinitek 500	26	5	0	0	0	0	0	0
Siemens Clinitek Advantus	53	10	0	0	0	0	0	0
Roche Chemstrip Criterion	5	3	0	0	0	0	0	0
Roche Urisys	21	7	0	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0	0	0
IRIS AUTION MAX	12	13	0	0	0	0	0	0
IRIS iChem Velocity	0	26	0	0	0	0	0	0
Other Chemistry Analyzer	0	0	0	0	0	0	0	0
Other	1	0	0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips	151	31	0	0	0	0	0	0
Siemens Atlas Reagent	78	7	0	0	0	0	0	0
Roche Chemstrips	25	9	0	0	0	0	0	0
Roche Urisys Cassette	6	5	0	0	0	0	0	0
ARKRAY Aution Sticks	10	10	0	0	0	0	0	0
IRIS Reagent Strips	2	29	0	0	0	0	0	0
Bayer Clinitest	44	4	0	0	0	0	0	0
Bayer Diastix	0	0	0	0	0	0	0	0
Siemens Uristix	3	1	0	0	0	0	0	0
All Instruments and Reagents	319	96	0	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Glucose	Specimen: U29	Negative	Trace 50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination		6	19	42	25	3	3	0
Siemens Atlas		0	2	80	15	0	0	0
Siemens Clinitek Status		0	46	0	1	0	0	0
Siemens Clinitek 50		0	5	0	0	0	0	0
Siemens Clinitek 100		3	2	0	0	0	0	0
Siemens Clinitek 500		11	18	3	0	0	0	0
Siemens Clinitek Advantus		18	45	0	0	0	0	0
Roche Chemstrip Criterion		0	0	6	1	1	0	0
Roche Urisys		1	0	23	0	0	4	0
Roche Cobas U411		0	0	5	0	0	0	0
IRIS AUTION MAX		0	0	22	3	0	0	0
IRIS iChem Velocity		0	0	1	25	0	0	0
Other Chemistry Analyzer		0	0	1	0	0	0	0
Other		0	0	1	0	0	0	0
<Reagents>								
Siemens Reagent Strips		38	126	16	3	0	0	0
Siemens Atlas Reagent		0	2	69	14	0	0	0
Roche Chemstrips		1	1	31	3	1	1	0
Roche Urisys Cassette		0	0	7	0	0	4	0
ARKRAY Aution Sticks		0	0	17	3	0	0	0
IRIS Reagent Strips		0	0	6	25	0	0	0
Bayer Clinitest		0	3	38	22	3	2	0
Bayer Diastix		0	1	0	0	0	0	0
Siemens Uristix		0	4	0	0	0	0	0
All Instruments and Reagents		39	137	184	70	4	7	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Glucose	Specimen: U30	Negative	Trace 50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination	72		1	0	0	0	0	0
Siemens Atlas	97		0	0	0	0	0	0
Siemens Clinitek Status	47		0	0	0	0	0	0
Siemens Clinitek 50	5		0	0	0	0	0	0
Siemens Clinitek 100	5		0	0	0	0	0	0
Siemens Clinitek 500	30		0	0	0	0	0	0
Siemens Clinitek Advantus	63		0	0	0	0	0	0
Roche Chemstrip Criterion	8		0	0	0	0	0	0
Roche Urisys	27		0	1	0	0	0	0
Roche Cobas U411	5		0	0	0	0	0	0
IRIS AUTION MAX	24		1	0	0	0	0	0
IRIS iChem Velocity	3		23	0	0	0	0	0
Other Chemistry Analyzer	0		0	0	0	0	0	0
Other	1		0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips	182		0	0	0	0	0	0
Siemens Atlas Reagent	85		0	0	0	0	0	0
Roche Chemstrips	32		1	1	0	0	0	0
Roche Urisys Cassette	11		0	0	0	0	0	0
ARKRAY Aution Sticks	19		1	0	0	0	0	0
IRIS Reagent Strips	8		23	0	0	0	0	0
Bayer Clinitest	47		0	0	0	0	0	0
Bayer Diastix	0		0	0	0	0	0	0
Siemens Uristix	3		0	0	0	0	0	0
All Instruments and Reagents	387		25	1	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Hemoglobin Specimen: U26	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	25	0	0	0	0	0
Siemens Atlas	94	1	1	0	0	0
Siemens Clinitek Status	46	2	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	30	0	0	0	0	0
Siemens Clinitek Advantus	61	1	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	28	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	180	2	0	0	0	0
Siemens Atlas Reagent	81	1	1	0	0	0
Roche Chemstrips	35	0	0	0	0	0
Roche Urisys Cassette	11	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	31	0	0	0	0	0
Siemens Uristix	0	1	0	0	0	0
None, no reagent or strip use	0	0	0	0	0	0
All Instruments and Reagents	358	4	1	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Hemoglobin Specimen: U27	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	24	0	0	0	0	1
Siemens Atlas	96	0	0	0	0	0
Siemens Clinitek Status	45	2	0	1	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	30	0	0	0	0	0
Siemens Clinitek Advantus	62	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	28	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	179	2	0	1	0	0
Siemens Atlas Reagent	83	0	0	0	0	0
Roche Chemstrips	34	0	0	0	0	1
Roche Urisys Cassette	11	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	31	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
None, no reagent or strip use	0	0	0	0	0	0
All Instruments and Reagents	359	2	0	1	0	1

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Hemoglobin Specimen: U28	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	0	0	0	1	0	27
Siemens Atlas	0	0	0	3	0	93
Siemens Clinitek Status	0	0	0	14	0	34
Siemens Clinitek 50	0	0	0	0	0	5
Siemens Clinitek 100	0	0	0	0	0	5
Siemens Clinitek 500	0	0	0	3	0	28
Siemens Clinitek Advantus	0	0	0	2	0	60
Roche Chemstrip Criterion	0	0	0	0	0	8
Roche Urisys	0	0	0	2	1	25
Roche Cobas U411	0	0	0	0	1	4
IRIS AUTION MAX	0	0	0	21	0	4
IRIS iChem Velocity	0	0	0	26	0	0
<Reagents>						
Siemens Reagent Strips	0	0	0	18	0	167
Siemens Atlas Reagent	0	0	0	3	0	80
Roche Chemstrips	0	0	0	3	2	31
Roche Urisys Cassette	0	0	0	0	0	11
ARKRAY Aution Sticks	0	0	0	18	0	2
IRIS Reagent Strips	0	0	0	29	0	2
Siemens Uristix	0	0	0	1	0	0
None, no reagent or strip use	0	0	0	0	0	0
All Instruments and Reagents	0	0	0	72	2	293

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Hemoglobin Specimen: U29	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	25	0	0	0	0	0
Siemens Atlas	96	0	0	0	0	0
Siemens Clinitek Status	45	2	1	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0
Siemens Clinitek 500	29	1	0	0	0	0
Siemens Clinitek Advantus	61	1	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0
Roche Urisys	27	0	0	0	0	1
Roche Cobas U411	5	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	177	4	1	0	0	0
Siemens Atlas Reagent	83	0	0	0	0	0
Roche Chemstrips	34	0	0	0	0	1
Roche Urisys Cassette	11	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	31	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
None, no reagent or strip use	0	0	0	0	0	0
All Instruments and Reagents	357	4	1	0	0	1

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Hemoglobin Specimen: U30	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	0	0	0	1	0	27
Siemens Atlas	0	0	0	9	0	87
Siemens Clinitek Status	0	0	0	30	0	18
Siemens Clinitek 50	0	0	0	2	0	3
Siemens Clinitek 100	0	0	0	1	0	4
Siemens Clinitek 500	0	0	0	8	0	23
Siemens Clinitek Advantus	0	0	0	7	0	55
Roche Chemstrip Criterion	0	0	0	0	0	8
Roche Urisys	1	0	0	3	5	19
Roche Cobas U411	0	0	0	0	1	4
IRIS AUTION MAX	0	0	0	24	0	1
IRIS iChem Velocity	1	0	0	25	0	0
<Reagents>						
Siemens Reagent Strips	0	0	0	49	0	136
Siemens Atlas Reagent	0	0	0	7	0	76
Roche Chemstrips	1	0	0	2	4	29
Roche Urisys Cassette	0	0	0	2	2	7
ARKRAY Aution Sticks	0	0	0	20	0	0
IRIS Reagent Strips	1	0	0	29	0	1
Siemens Uristix	0	0	0	1	0	0
None, no reagent or strip use	0	0	0	0	0	0
All Instruments and Reagents	2	0	0	110	6	249

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Ketones	Specimen: U26	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>							
Visual examination	30	0	0	0	0	0	0
Siemens Atlas	97	0	0	0	0	0	0
Siemens Clinitek Status	48	0	0	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0	0
Siemens Clinitek 500	30	0	0	0	0	0	0
Siemens Clinitek Advantus	62	0	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0	0
Roche Urisys	28	0	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	182	0	0	0	0	0	0
Siemens Atlas Reagent	84	0	0	0	0	0	0
Roche Chemstrips	36	0	0	0	0	0	0
Roche Urisys Cassette	11	0	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0	0
IRIS Reagent Strips	32	0	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0	0
Siemens Ketostix	0	0	0	0	0	0	0
All Instruments and Reagents	369	0	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Ketones	Specimen: U27	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>							
Visual examination	30	0	0	0	0	0	0
Siemens Atlas	97	0	0	0	0	0	0
Siemens Clinitek Status	47	0	0	0	0	1	0
Siemens Clinitek 50	5	0	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0	0
Siemens Clinitek 500	30	0	0	0	0	0	0
Siemens Clinitek Advantus	62	0	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0	0
Roche Urisys	28	0	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	181	0	0	0	0	1	0
Siemens Atlas Reagent	84	0	0	0	0	0	0
Roche Chemstrips	36	0	0	0	0	0	0
Roche Urisys Cassette	11	0	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0	0
IRIS Reagent Strips	32	0	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0	0
Siemens Ketostix	0	0	0	0	0	0	0
All Instruments and Reagents	368	0	0	0	0	1	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Ketones	Specimen: U28	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>							
Visual examination	30	0	0	0	0	0	0
Siemens Atlas	97	0	0	0	0	0	0
Siemens Clinitek Status	48	0	0	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0	0
Siemens Clinitek 500	30	0	0	0	0	0	0
Siemens Clinitek Advantus	62	0	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0	0
Roche Urisys	28	0	0	0	0	0	0
Roche Cobas U411	5	0	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	182	0	0	0	0	0	0
Siemens Atlas Reagent	84	0	0	0	0	0	0
Roche Chemstrips	36	0	0	0	0	0	0
Roche Urisys Cassette	11	0	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0	0
IRIS Reagent Strips	32	0	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0	0
Siemens Ketostix	0	0	0	0	0	0	0
All Instruments and Reagents	369	0	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Ketones	Specimen: U29	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>							
Visual examination		0	0	0	6	40	3
Siemens Atlas		0	0	0	0	97	0
Siemens Clinitek Status		1	0	0	1	46	0
Siemens Clinitek 50		0	0	0	0	5	0
Siemens Clinitek 100		0	0	0	0	5	0
Siemens Clinitek 500		0	0	0	4	27	0
Siemens Clinitek Advantus		0	0	0	5	57	0
Roche Chemstrip Criterion		0	0	0	1	7	0
Roche Urisys		2	0	0	2	24	0
Roche Cobas U411		0	0	0	0	5	0
IRIS AUTION MAX		0	0	0	0	25	0
IRIS iChem Velocity		0	0	0	13	13	0
<Reagents>							
Siemens Reagent Strips		1	0	0	10	173	0
Siemens Atlas Reagent		0	0	0	0	84	0
Roche Chemstrips		2	0	0	5	30	0
Roche Urisys Cassette		0	0	0	0	11	0
ARKRAY Aution Sticks		0	0	0	0	19	0
IRIS Reagent Strips		0	0	0	13	19	0
Siemens Acetest		0	0	0	4	12	2
Bayer Diastix		0	0	0	0	0	1
Siemens Uristix		0	0	0	0	1	0
Siemens Ketostix		0	0	0	0	2	0
All Instruments and Reagents		3	0	0	32	351	3

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Ketones	Specimen: U30	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>							
Visual examination	30	0	0	0	0	0	0
Siemens Atlas	97	0	0	0	0	0	0
Siemens Clinitek Status	48	0	0	0	0	0	0
Siemens Clinitek 50	5	0	0	0	0	0	0
Siemens Clinitek 100	5	0	0	0	0	0	0
Siemens Clinitek 500	30	0	0	0	0	0	0
Siemens Clinitek Advantus	62	0	0	0	0	0	0
Roche Chemstrip Criterion	8	0	0	0	0	0	0
Roche Urisys	27	0	0	0	0	1	0
Roche Cobas U411	5	0	0	0	0	0	0
IRIS AUTION MAX	25	0	0	0	0	0	0
IRIS iChem Velocity	26	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	182	0	0	0	0	0	0
Siemens Atlas Reagent	84	0	0	0	0	0	0
Roche Chemstrips	35	0	0	0	0	1	0
Roche Urisys Cassette	11	0	0	0	0	0	0
ARKRAY Aution Sticks	19	0	0	0	0	0	0
IRIS Reagent Strips	32	0	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0	0
Siemens Ketostix	0	0	0	0	0	0	0
All Instruments and Reagents	368	0	0	0	0	1	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Protein	Specimen: U26	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination		1	3	44	9	0	0	7
Siemens Atlas		0	1	94	2	0	0	0
Siemens Clinitek Status		1	0	12	34	0	0	0
Siemens Clinitek 50		0	0	2	3	0	0	0
Siemens Clinitek 100		0	0	3	2	0	0	0
Siemens Clinitek 500		1	6	24	0	0	0	0
Siemens Clinitek Advantus		1	6	56	0	0	0	0
Roche Chemstrip Criterion		0	1	7	0	0	0	0
Roche Urisys		0	5	15	8	0	0	0
Roche Cobas U411		0	0	4	1	0	0	0
IRIS AUTION MAX		0	1	24	0	0	0	0
IRIS iChem Velocity		0	0	26	0	0	0	0
Hitachi/Roche Analyzer		0	0	0	0	0	0	1
<Reagents>								
Siemens Reagent Strips		3	11	123	43	0	0	0
Siemens Atlas Reagent		0	1	84	1	0	0	0
Sulfosalicylic Acid		1	3	21	4	0	0	7
Roche Chemstrips		0	4	22	10	0	0	0
Roche Urisys Cassette		0	2	9	0	0	0	0
ARKRAY Aution Sticks		0	1	19	0	0	0	0
IRIS Reagent Strips		0	0	31	0	0	0	0
Siemens Uristix		0	1	2	1	0	0	0
Hitachi/Roche Reagents		0	0	0	0	0	0	1
All Instruments and Reagents		4	23	311	59	0	0	8

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Protein	Specimen: U27	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination		0	0	2	40	16	0	9
Siemens Atlas		1	0	0	91	5	0	0
Siemens Clinitek Status		1	0	0	19	27	0	0
Siemens Clinitek 50		0	0	0	0	5	0	0
Siemens Clinitek 100		0	0	0	2	3	0	0
Siemens Clinitek 500		0	0	1	30	0	0	0
Siemens Clinitek Advantus		0	0	0	62	1	0	0
Roche Chemstrip Criterion		0	0	0	7	1	0	0
Roche Urisys		0	0	0	28	0	0	0
Roche Cobas U411		0	0	0	5	0	0	0
IRIS AUTION MAX		0	0	0	25	0	0	0
IRIS iChem Velocity		0	0	0	26	0	0	0
Hitachi/Roche Analyzer		0	0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips		1	0	1	133	45	0	0
Siemens Atlas Reagent		1	0	0	80	5	0	0
Sulfosalicylic Acid		0	0	2	23	5	0	9
Roche Chemstrips		0	0	0	34	2	0	0
Roche Urisys Cassette		0	0	0	11	0	0	0
ARKRAY Aution Sticks		0	0	0	20	0	0	0
IRIS Reagent Strips		0	0	0	31	0	0	0
Siemens Uristix		0	0	0	3	1	0	0
Hitachi/Roche Reagents		0	0	0	0	0	0	0
All Instruments and Reagents		2	0	3	335	58	0	9

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Protein	Specimen: U28	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination		0	0	1	2	44	12	9
Siemens Atlas		0	0	0	1	73	23	0
Siemens Clinitek Status		1	0	0	0	45	1	0
Siemens Clinitek 50		0	0	0	0	5	0	0
Siemens Clinitek 100		0	0	0	1	4	0	0
Siemens Clinitek 500		0	0	0	3	28	0	0
Siemens Clinitek Advantus		0	0	0	3	60	0	0
Roche Chemstrip Criterion		0	0	0	0	8	0	0
Roche Urisys		0	0	0	0	28	0	0
Roche Cobas U411		0	0	0	0	5	0	0
IRIS AUTION MAX		0	0	0	0	25	0	0
IRIS iChem Velocity		0	0	0	0	26	0	0
Hitachi/Roche Analyzer		0	0	0	0	0	0	1
<Reagents>								
Siemens Reagent Strips		1	0	0	6	169	5	0
Siemens Atlas Reagent		0	0	0	1	63	22	0
Sulfosalicylic Acid		0	0	1	2	18	9	9
Roche Chemstrips		0	0	0	0	36	0	0
Roche Urisys Cassette		0	0	0	0	11	0	0
ARKRAY Aution Sticks		0	0	0	0	20	0	0
IRIS Reagent Strips		0	0	0	0	31	0	0
Siemens Uristix		0	0	0	1	3	0	0
Hitachi/Roche Reagents		0	0	0	0	0	0	1
All Instruments and Reagents		1	0	1	10	351	36	10

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Protein	Specimen: U29	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination		0	1	32	25	0	0	8
Siemens Atlas		0	0	23	74	0	0	0
Siemens Clinitek Status		1	0	1	45	0	0	0
Siemens Clinitek 50		0	0	0	5	0	0	0
Siemens Clinitek 100		0	0	1	4	0	0	0
Siemens Clinitek 500		0	0	29	2	0	0	0
Siemens Clinitek Advantus		0	0	56	7	0	0	0
Roche Chemstrip Criterion		0	0	6	2	0	0	0
Roche Urisys		0	0	7	20	1	0	0
Roche Cobas U411		0	0	0	5	0	0	0
IRIS AUTION MAX		0	1	22	2	0	0	0
IRIS iChem Velocity		0	0	4	22	0	0	0
Hitachi/Roche Analyzer		0	0	0	0	0	0	1
<Reagents>								
Siemens Reagent Strips		1	0	98	81	0	0	0
Siemens Atlas Reagent		0	0	18	68	0	0	0
Sulfosalicylic Acid		0	1	19	10	0	0	8
Roche Chemstrips		0	0	16	19	1	0	0
Roche Urisys Cassette		0	0	1	10	0	0	0
ARKRAY Aution Sticks		0	1	18	1	0	0	0
IRIS Reagent Strips		0	0	8	23	0	0	0
Siemens Uristix		0	0	3	1	0	0	0
Hitachi/Roche Reagents		0	0	0	0	0	0	1
All Instruments and Reagents		1	2	181	213	1	0	9

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Protein	Specimen: U30	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination		0	0	1	2	50	6	9
Siemens Atlas		0	0	0	1	89	7	0
Siemens Clinitek Status		1	0	0	0	46	0	0
Siemens Clinitek 50		0	0	0	0	5	0	0
Siemens Clinitek 100		0	0	0	1	4	0	0
Siemens Clinitek 500		0	0	0	12	19	0	0
Siemens Clinitek Advantus		0	0	0	21	42	0	0
Roche Chemstrip Criterion		0	0	0	0	8	0	0
Roche Urisys		0	0	0	2	26	0	0
Roche Cobas U411		0	0	0	0	5	0	0
IRIS AUTION MAX		0	0	0	4	21	0	0
IRIS iChem Velocity		0	0	0	1	25	0	0
Hitachi/Roche Analyzer		0	0	0	0	0	0	1
<Reagents>								
Siemens Reagent Strips		1	0	0	33	146	1	0
Siemens Atlas Reagent		0	0	0	1	78	7	0
Sulfosalicylic Acid		0	0	1	2	22	5	9
Roche Chemstrips		0	0	0	2	34	0	0
Roche Urisys Cassette		0	0	0	0	11	0	0
ARKRAY Aution Sticks		0	0	0	3	17	0	0
IRIS Reagent Strips		0	0	0	2	29	0	0
Siemens Uristix		0	0	0	1	3	0	0
Hitachi/Roche Reagents		0	0	0	0	0	0	1
All Instruments and Reagents		1	0	1	44	340	13	10

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Qualitative hCG report:

Alere hCG Cassette	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	13	13	0	13	0
Indeterminant	0	0	0	0	0
Positive	0	0	13	0	13
Siemens Clinitek Status	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	9	9	0	9	0
Indeterminant	0	0	0	0	0
Positive	0	0	9	0	9
Beckman Coulter ICON II	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	2	2	0	2	0
Indeterminant	0	0	0	0	0
Positive	0	0	2	0	2
Beckman Coulter ICON 25	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	21	21	0	21	0
Indeterminant	0	0	0	0	0
Positive	0	0	21	0	21
BTNX Rapid Response	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	2	2	0	2	0
Indeterminant	0	0	0	0	0
Positive	0	0	2	0	2
Beckman Coulter ICON 20	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	23	23	0	23	0
Indeterminant	0	0	0	0	0
Positive	0	0	23	0	23
Cardinal Health hCG Combo	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	28	28	0	28	0
Indeterminant	0	0	0	0	0
Positive	0	0	28	0	28

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Wampole Clearview hCG Combo	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Cardinal Health hCG Cassette	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	4	4	0	4	0
Indeterminant	0	0	0	0	0
Positive	0	0	4	0	4
Clarity hCG Cassette	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Cen Med Elite	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	16	16	0	16	0
Indeterminant	0	0	0	0	0
Positive	0	0	16	0	16
Sekisui OSOM Combo (Genzyme)	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	16	16	0	16	0
Indeterminant	0	0	0	0	0
Positive	0	0	16	0	16
Henry Schein One Step	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	3	3	0	3	0
Indeterminant	0	0	0	0	0
Positive	0	0	3	0	3
Immunostics hCG Detector	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Inverness Medical BioStar Acceava	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Jant Accutest Value±	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Lifesign Status (serum/urine)	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
McKesson hCG Cassette	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Meridian Diagnostics	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Polymedco Poly stat	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	5	5	0	5	0
Indeterminant	0	0	0	0	0
Positive	0	0	5	0	5
Quidel QuickVue	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	68	68	1	67	1
Indeterminant	0	0	0	0	0
Positive	0	0	67	1	67

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

R. Medical One Step Strip	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	2	2	0	2	0
Indeterminant	0	0	0	0	0
Positive	0	0	2	0	2
SA Scientific	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	1	1	0	1	0
Indeterminant	0	0	0	0	0
Positive	0	0	1	0	1
Stanbio QuPID	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	13	13	0	13	0
Indeterminant	0	0	0	0	0
Positive	0	0	13	0	13
Fisher Sure-Vue serum/urine	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	41	41	0	41	0
Indeterminant	0	0	0	0	0
Positive	0	0	41	0	41
Fisher Sure-Vue STAT serum/urine	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	26	26	0	26	0
Indeterminant	0	0	0	0	0
Positive	0	0	26	0	26
Signify hCG Urine/Serum	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	3	3	0	3	0
Indeterminant	0	0	0	0	0
Positive	0	0	3	0	3
Siemens Immulite	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	5	5	0	5	0
Indeterminant	0	0	0	0	0
Positive	0	0	5	0	5

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

Fisher Sure-View (urine only)	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	8	8	0	8	0
Indeterminant	0	0	0	0	0
Positive	0	0	8	0	8

True Medix	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	5	5	0	5	0
Indeterminant	0	0	0	0	0
Positive	0	0	5	0	5

All Methods & Instruments	Specimen:U26	Specimen:U27	Specimen:U28	Specimen:U29	Specimen:U30
Negative	322	322	1	321	1
Indeterminant	0	0	0	0	0
Positive	0	0	321	1	321

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 27 October 2014

hCG Summary:

[Code], hCG Kit name	U26			U27			U28			U29			U30		
	N	I	P	N	I	P	N	I	P	N	I	P	N	I	P

[AC] Alere hCG Cassette	13	0	0	13	0	0	0	0	13	13	0	0	0	0	13
[BA] Siemens Clinitek Status	9	0	0	9	0	0	0	0	9	9	0	0	0	0	9
[BC] Beckman Coulter ICON II	2	0	0	2	0	0	0	0	2	2	0	0	0	0	2
[BN] Beckman Coulter ICON 25	21	0	0	21	0	0	0	0	21	21	0	0	0	0	21
[BR] BTNX Rapid Response	2	0	0	2	0	0	0	0	2	2	0	0	0	0	2
[BT] Beckman Coulter ICON 20	23	0	0	23	0	0	0	0	23	23	0	0	0	0	23
[CA] Cardinal Health hCG Combo	28	0	0	28	0	0	0	0	28	28	0	0	0	0	28
[CC] Wampole Clearview hCG Combo	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[CH] Cardinal Health hCG Cassette	4	0	0	4	0	0	0	0	4	4	0	0	0	0	4
[CL] Clarity hCG Cassette	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[CM] Cen Med Elite	16	0	0	16	0	0	0	0	16	16	0	0	0	0	16
[GO] Sekisui OSOM Combo (Genzyme)	16	0	0	16	0	0	0	0	16	16	0	0	0	0	16
[HS] Henry Schein One Step	3	0	0	3	0	0	0	0	3	3	0	0	0	0	3
[ID] Immunostics hCG Detector	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[IM] Inverness Medical BioStar Acceava	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[JA] Jant Accutest Value+	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[LS] Lifesign Status (serum/urine)	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[MC] McKesson hCG Cassette	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[MD] Meridian Diagnostics	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[PM] Polymedco Poly stat	5	0	0	5	0	0	0	0	5	5	0	0	0	0	5
[QQ] Quidel QuickVue	68	0	0	68	0	0	1	0	67	67	0	1	1	0	67
[RM] R. Medical One Step Strip	2	0	0	2	0	0	0	0	2	2	0	0	0	0	2
[SA] SA Scientific	1	0	0	1	0	0	0	0	1	1	0	0	0	0	1
[SB] Stanbio QuPID	13	0	0	13	0	0	0	0	13	13	0	0	0	0	13
[SE] Fisher Sure-Vue serum/urine	41	0	0	41	0	0	0	0	41	41	0	0	0	0	41
[SF] Fisher Sure-Vue STAT serum/urine	26	0	0	26	0	0	0	0	26	26	0	0	0	0	26
[SG] Signify hCG Urine/Serum	3	0	0	3	0	0	0	0	3	3	0	0	0	0	3
[SI] Siemens Immulite	5	0	0	5	0	0	0	0	5	5	0	0	0	0	5
[SU] Fisher Sure-Vue (urine only)	8	0	0	8	0	0	0	0	8	8	0	0	0	0	8
[TX] True Medix	5	0	0	5	0	0	0	0	5	5	0	0	0	0	5

[99] All hCG Methods	322	0	0	322	0	0	1	0	321	321	0	1	1	0	321