

New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004

Specific Gravity	N	U21	U22	U23	U24	U25
[01] Bayer Multistix (Manual)	51	1.017 ± 0.004	1.019 ± 0.003	1.016 ± 0.004	1.019 ± 0.004	1.023 ± 0.003
[02] Bayer Clinitek	172	1.018 ± 0.003	1.018 ± 0.003	1.022 ± 0.003	1.019 ± 0.003	1.021 ± 0.004
[03] Bayer Atlas	71	1.019 ± 0.002	1.015 ± 0.002	1.016 ± 0.002	1.014 ± 0.002	1.016 ± 0.002
[10] Roche Manual (Chemstrip)	17	1.010 ± 0.001	1.015 ± 0.001	1.010 ± 0.002	1.015 ± 0.001	1.017 ± 0.002
[11] Roche UA Reader/Super UA Reade	27	1.011 ± 0.002	1.015 ± 0.002	1.013 ± 0.003	1.016 ± 0.003	1.017 ± 0.003
[12] Roche Mini UA Reader	5	1.010 ± 0.000	1.015 ± 0.000	1.012 ± 0.003	1.015 ± 0.000	1.017 ± 0.003
[13] Roche Chemstrip Criterion	23	1.010 ± 0.003	1.015 ± 0.001	1.012 ± 0.003	1.015 ± 0.003	1.017 ± 0.003
[16] IRIS Mass Gravity Meter	6	1.020 ± 0.003	1.016 ± 0.001	1.018 ± 0.001	1.015 ± 0.001	1.017 ± 0.001
[17] IRIS BMC UA Reader	2	1.016 ± 0.000	1.011 ± 0.000	1.013 ± 0.003	1.011 ± 0.006	1.014 ± 0.001
[18] IRIS Clinitek Reader	1	1.021 ± 0.000	1.016 ± 0.000	1.018 ± 0.000	1.015 ± 0.000	1.017 ± 0.000
[20] Bayer Clinitest	3	1.020 ± 0.000	1.017 ± 0.003	1.020 ± 0.005	1.017 ± 0.003	1.017 ± 0.003
[23] Refractometer	33	1.019 ± 0.001	1.014 ± 0.002	1.016 ± 0.001	1.013 ± 0.001	1.015 ± 0.001
[25] Hitachi/Roche Analyzer	2	1.019 ± 0.000	1.013 ± 0.000	1.016 ± 0.000	1.013 ± 0.000	1.015 ± 0.000
[30] Other Chemistry Analyzer	11	1.013 ± 0.004	1.015 ± 0.005	1.013 ± 0.003	1.015 ± 0.006	1.018 ± 0.004
[90] Method not listed	7	1.014 ± 0.005	1.014 ± 0.002	1.013 ± 0.003	1.014 ± 0.003	1.017 ± 0.003
[99] All methods	431	1.017 ± 0.004	1.017 ± 0.003	1.017 ± 0.005	1.017 ± 0.003	1.019 ± 0.004

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – October 25, 2004**

pH	N	U21	U22	U23	U24	U25
[01] Bayer Multistix (Manual)	64	6.91 ± 0.36	6.28 ± 0.38	7.59 ± 0.29	6.36 ± 0.31	5.20 ± 0.39
[02] Bayer Clinitek	179	7.02 ± 0.29	6.39 ± 0.30	7.64 ± 0.26	6.44 ± 0.21	5.17 ± 0.26
[03] Bayer Atlas	70	6.99 ± 0.13	6.50 ± 0.00	7.52 ± 0.10	6.51 ± 0.06	5.46 ± 0.19
[10] Roche Manual (Chemstrip)	27	7.00 ± 0.00	6.00 ± 0.24	7.87 ± 0.30	6.00 ± 0.24	5.00 ± 0.00
[11] Roche UA Reader/Super UA Reade	29	7.00 ± 0.00	6.52 ± 0.21	7.97 ± 0.19	6.47 ± 0.19	5.07 ± 0.26
[12] Roche Mini UA Reader	4	7.00 ± 0.00	6.50 ± 0.00	8.00 ± 0.00	6.50 ± 0.00	5.00 ± 0.00
[13] Roche Chemstrip Criterion	25	7.00 ± 0.00	6.44 ± 0.22	7.88 ± 0.33	6.34 ± 0.28	5.08 ± 0.28
[17] IRIS BMC UA Reader	7	7.00 ± 0.00	6.14 ± 0.38	7.86 ± 0.24	6.00 ± 0.50	5.07 ± 0.19
[18] IRIS Clinitek Reader	3	7.00 ± 0.00	6.00 ± 0.00	7.83 ± 0.29	6.00 ± 0.00	5.00 ± 0.00
[20] Bayer Clinitest	3	7.33 ± 0.29	6.50 ± 0.00	7.83 ± 0.29	6.67 ± 0.29	5.17 ± 0.29
[24] pH Meter	2	6.90 ± 0.14	6.35 ± 0.07	7.65 ± 0.07	6.40 ± 0.14	5.35 ± 0.07
[25] Hitachi/Roche Analyzer	3	7.00 ± 0.00	6.50 ± 0.00	7.67 ± 0.58	6.33 ± 0.29	5.00 ± 0.00
[29] Chimera UA Perfect	1	7.00 ± 0.00	6.30 ± 0.00	7.20 ± 0.00	6.30 ± 0.00	5.80 ± 0.00
[30] Other Chemistry Analyzer	12	7.00 ± 0.00	6.29 ± 0.26	7.79 ± 0.33	6.33 ± 0.25	5.00 ± 0.00
[90] Method not listed	7	7.00 ± 0.00	6.43 ± 0.19	7.93 ± 0.19	6.29 ± 0.27	5.14 ± 0.38
[95] Method not specified	2	7.00 ± 0.00	6.25 ± 0.35	7.50 ± 0.71	6.25 ± 0.35	5.00 ± 0.00
[99] All methods	438	7.00 ± 0.24	6.37 ± 0.30	7.68 ± 0.29	6.39 ± 0.26	5.19 ± 0.29

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Bilirubin - Specimen: U21		Number of Laboratories Reporting Results				
		Negative	Small	Moderate	Large	Positive
[01]	Bayer Multistix (Manual)	59	0	0	0	0
[02]	Bayer Clinitek	176	0	0	0	0
[03]	Bayer Atlas	71	0	0	0	0
[10]	Roche Manual (Chemstrip)	21	0	0	0	0
[11]	Roche UA Reader/Super UA	29	0	0	0	0
[12]	Roche Mini UA Reader	5	0	0	0	0
[13]	Roche Chemstrip Criterion	26	0	0	0	0
[17]	IRIS BMC UA Reader	6	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0
[20]	Bayer Clinitest	3	0	0	0	0
[21]	Bayer Ictotest	18	0	0	0	0
[25]	Hitachi/Roche Analyzer	3	0	0	0	0
[30]	Other Chemistry Analyzer	12	0	0	0	0
[90]	Method not listed	7	0	0	0	0
[95]	Method not specified	2	0	0	0	0
[99]	All methods	441	0	0	0	0

Bilirubin - Specimen: U22		Number of Laboratories Reporting Results				
		Negative	Small	Moderate	Large	Positive
[01]	Bayer Multistix (Manual)	1	3	34	19	1
[02]	Bayer Clinitek	2	11	133	31	0
[03]	Bayer Atlas	0	0	61	10	0
[10]	Roche Manual (Chemstrip)	4	5	12	1	0
[11]	Roche UA Reader/Super UA	2	26	0	0	1
[12]	Roche Mini UA Reader	0	5	0	0	0
[13]	Roche Chemstrip Criterion	21	5	0	0	0
[17]	IRIS BMC UA Reader	1	5	0	0	0
[18]	IRIS Clinitek Reader	0	3	0	0	0
[20]	Bayer Clinitest	0	1	2	0	0
[21]	Bayer Ictotest	2	1	19	3	149
[25]	Hitachi/Roche Analyzer	0	2	1	0	0
[30]	Other Chemistry Analyzer	0	10	2	0	0
[90]	Method not listed	1	6	0	0	0
[95]	Method not specified	0	0	1	1	2
[99]	All methods	34	83	265	65	153

New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004
Number of Laboratories Reporting Results

Bilirubin - Specimen: U23		Negative	Small	Moderate	Large	Positive
[01] Bayer Multistix (Manual)	59	0	0	0	0	0
[02] Bayer Clinitek	176	0	0	0	0	0
[03] Bayer Atlas	71	0	0	0	0	0
[10] Roche Manual (Chemstrip)	21	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0
[17] IRIS BMC UA Reader	6	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0
[20] Bayer Clinitest	3	0	0	0	0	0
[21] Bayer Ictotest	18	0	0	0	0	0
[25] Hitachi/Roche Analyzer	3	0	0	0	0	0
[30] Other Chemistry Analyzer	12	0	0	0	0	0
[90] Method not listed	7	0	0	0	0	0
[95] Method not specified	1	0	0	0	0	0
[99] All methods	440	0	0	0	0	0

Bilirubin - Specimen: U24		Negative	Small	Moderate	Large	Positive
[01] Bayer Multistix (Manual)	59	0	0	0	0	0
[02] Bayer Clinitek	176	0	0	0	0	0
[03] Bayer Atlas	71	0	0	0	0	0
[10] Roche Manual (Chemstrip)	21	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0
[17] IRIS BMC UA Reader	6	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0
[20] Bayer Clinitest	3	0	0	0	0	0
[21] Bayer Ictotest	18	0	0	0	0	0
[25] Hitachi/Roche Analyzer	3	0	0	0	0	0
[30] Other Chemistry Analyzer	12	0	0	0	0	0
[90] Method not listed	7	0	0	0	0	0
[95] Method not specified	1	0	0	0	0	0
[99] All methods	440	0	0	0	0	0

Bilirubin - Specimen: U25		Negative	Small	Moderate	Large	Positive
[01] Bayer Multistix (Manual)	59	0	0	0	0	0
[02] Bayer Clinitek	175	1	0	0	0	0
[03] Bayer Atlas	71	0	0	0	0	0
[10] Roche Manual (Chemstrip)	21	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0
[17] IRIS BMC UA Reader	6	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0
[20] Bayer Clinitest	3	0	0	0	0	0
[21] Bayer Ictotest	18	0	0	0	0	0
[25] Hitachi/Roche Analyzer	3	0	0	0	0	0
[30] Other Chemistry Analyzer	12	0	0	0	0	0
[90] Method not listed	7	0	0	0	0	0
[95] Method not specified	1	0	0	0	0	0
[99] All methods	439	1	0	0	0	0

New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004

Glucose - Specimen: U21	Negative	50 mg/dL	100 mg/dL	250 mg/dL	500 mg/dL	1000 mg/dL	2000 mg/dL
[01] Bayer Multistix (Manual)	66	0	0	0	0	0	0
[02] Bayer Clinitek	179	0	1	0	0	0	0
[03] Bayer Atlas	71	0	0	0	0	0	0
[10] Roche Manual (Chemstrip)	28	0	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0	0
[17] IRIS BMC UA Reader	7	0	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0	0
[20] Bayer Clinitest	23	0	0	0	0	0	0
[25] Hitachi/Roche Analyzer	4	0	0	0	0	0	0
[26] Bayer Diastix	1	0	0	0	0	0	0
[30] Other Chemistry Analyzer	11	0	1	0	0	0	0
[90] Method not listed	12	0	0	0	0	0	0
[99] All methods	465	0	2	0	0	0	0

Glucose - Specimen: U22	Negative	50 mg/dL	100 mg/dL	250 mg/dL	500 mg/dL	1000 mg/dL	2000 mg/dL
[01] Bayer Multistix (Manual)	0	0	1	30	26	9	0
[02] Bayer Clinitek	0	0	5	102	66	7	0
[03] Bayer Atlas	0	0	1	69	2	0	0
[10] Roche Manual (Chemstrip)	0	0	0	15	14	6	0
[11] Roche UA Reader/Super UA	0	0	0	26	2	1	0
[12] Roche Mini UA Reader	0	0	0	5	0	0	0
[13] Roche Chemstrip Criterion	0	0	0	17	9	0	0
[17] IRIS BMC UA Reader	0	0	0	1	6	0	0
[18] IRIS Clinitek Reader	0	0	0	0	3	0	0
[20] Bayer Clinitest	0	0	3	49	46	2	1
[25] Hitachi/Roche Analyzer	0	0	1	4	0	0	0
[26] Bayer Diastix	0	0	0	4	0	0	0
[30] Other Chemistry Analyzer	0	0	0	7	4	2	0
[90] Method not listed	0	0	1	8	4	0	0
[99] All methods	0	0	12	337	182	27	1

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Glucose - Specimen: U23	Negative	50 mg/dL	100 mg/dL	250 mg/dL	500 mg/dL	1000 mg/dL	2000 mg/dL
[01] Bayer Multistix (Manual)	0	0	1	41	21	3	0
[02] Bayer Clinitek	1	0	3	134	40	2	0
[03] Bayer Atlas	0	0	0	16	56	0	0
[10] Roche Manual (Chemstrip)	0	0	0	1	10	24	0
[11] Roche UA Reader/Super UA	0	0	0	1	1	26	1
[12] Roche Mini UA Reader	0	0	0	1	0	4	0
[13] Roche Chemstrip Criterion	0	0	0	0	1	25	0
[17] IRIS BMC UA Reader	0	0	0	0	1	6	0
[18] IRIS Clinitek Reader	0	0	0	0	0	3	0
[20] Bayer Clinitest	0	0	0	32	48	21	1
[25] Hitachi/Roche Analyzer	0	0	0	1	1	3	0
[26] Bayer Diastix	0	0	0	3	1	0	0
[30] Other Chemistry Analyzer	0	0	0	0	3	10	0
[90] Method not listed	0	0	0	3	4	6	0
[99] All methods	1	0	4	233	187	133	2

Glucose - Specimen: U24	Negative	50 mg/dL	100 mg/dL	250 mg/dL	500 mg/dL	1000 mg/dL	2000 mg/dL
[01] Bayer Multistix (Manual)	0	0	1	0	4	38	23
[02] Bayer Clinitek	0	0	0	2	17	161	0
[03] Bayer Atlas	1	0	0	0	3	67	1
[10] Roche Manual (Chemstrip)	0	0	0	0	0	29	5
[11] Roche UA Reader/Super UA	0	0	0	0	0	28	1
[12] Roche Mini UA Reader	0	0	0	1	0	4	0
[13] Roche Chemstrip Criterion	0	0	0	0	1	25	0
[17] IRIS BMC UA Reader	0	0	0	0	0	7	0
[18] IRIS Clinitek Reader	0	0	0	0	0	3	0
[20] Bayer Clinitest	0	0	0	0	6	76	24
[25] Hitachi/Roche Analyzer	0	0	0	0	1	4	0
[26] Bayer Diastix	0	0	0	0	1	2	1
[30] Other Chemistry Analyzer	0	0	0	0	0	11	2
[90] Method not listed	0	0	0	0	0	12	1
[99] All methods	1	0	1	3	33	467	58

Glucose - Specimen: U25	Negative	50 mg/dL	100 mg/dL	250 mg/dL	500 mg/dL	1000 mg/dL	2000 mg/dL
[01] Bayer Multistix (Manual)	66	0	0	0	0	0	0
[02] Bayer Clinitek	178	0	2	0	0	0	0
[03] Bayer Atlas	71	0	0	0	0	0	0
[10] Roche Manual (Chemstrip)	30	0	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0	0
[17] IRIS BMC UA Reader	7	0	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0	0
[20] Bayer Clinitest	23	0	0	0	0	0	0
[25] Hitachi/Roche Analyzer	4	0	0	0	0	0	0
[26] Bayer Diastix	1	0	0	0	0	0	0
[30] Other Chemistry Analyzer	12	0	0	0	0	0	0
[90] Method not listed	13	0	0	0	0	0	0
[99] All methods	468	0	2	0	0	0	0

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Hemoglobin	Specimen: U21	Negative	Trace	Small	Moderate	Large
[01]	Bayer Multistix (Manual)	0	0	25	37	2
[02]	Bayer Clinitek	0	5	72	102	1
[03]	Bayer Atlas	1	3	43	24	0
[10]	Roche Manual (Chemstrip)	0	3	3	18	6
[11]	Roche UA Reader/Super UA	2	1	6	17	3
[12]	Roche Mini UA Reader	0	0	0	4	1
[13]	Roche Chemstrip Criterion	0	1	3	12	10
[17]	IRIS BMC UA Reader	0	0	3	4	0
[18]	IRIS Clinitek Reader	0	0	0	3	0
[20]	Bayer Clinitest	0	0	2	1	0
[25]	Hitachi/Roche Analyzer	0	0	1	2	0
[27]	Bayer Hemastix	0	0	3	3	0
[29]	Chimera UA Perfect	0	0	1	0	0
[30]	Other Chemistry Analyzer	0	0	4	5	3
[90]	Method not listed	0	2	1	3	1
[99]	All methods	3	15	167	235	27

Hemoglobin	Specimen: U22	Negative	Trace	Small	Moderate	Large
[01]	Bayer Multistix (Manual)	63	1	0	0	0
[02]	Bayer Clinitek	178	1	0	1	0
[03]	Bayer Atlas	71	0	0	0	0
[10]	Roche Manual (Chemstrip)	26	1	1	0	0
[11]	Roche UA Reader/Super UA	29	0	0	0	0
[12]	Roche Mini UA Reader	5	0	0	0	0
[13]	Roche Chemstrip Criterion	26	0	0	0	0
[17]	IRIS BMC UA Reader	7	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0
[20]	Bayer Clinitest	3	0	0	0	0
[25]	Hitachi/Roche Analyzer	3	0	0	0	0
[27]	Bayer Hemastix	2	0	0	0	0
[29]	Chimera UA Perfect	1	0	0	0	0
[30]	Other Chemistry Analyzer	12	0	0	0	0
[90]	Method not listed	7	0	0	0	0
[99]	All methods	436	3	1	1	0

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Hemoglobin	Specimen: U23	Negative	Trace	Small	Moderate	Large
[01]	Bayer Multistix (Manual)	16	41	8	0	0
[02]	Bayer Clinitek	10	157	13	0	0
[03]	Bayer Atlas	1	68	2	0	0
[10]	Roche Manual (Chemstrip)	14	14	1	0	0
[11]	Roche UA Reader/Super UA	16	13	0	0	0
[12]	Roche Mini UA Reader	1	4	0	0	0
[13]	Roche Chemstrip Criterion	20	5	1	0	0
[17]	IRIS BMC UA Reader	7	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0
[20]	Bayer Clinitest	0	2	1	0	0
[25]	Hitachi/Roche Analyzer	1	2	0	0	0
[27]	Bayer Hemastix	1	2	1	0	0
[29]	Chimera UA Perfect	0	0	0	1	0
[30]	Other Chemistry Analyzer	5	7	0	0	0
[90]	Method not listed	6	1	0	0	0
[99]	All methods	101	316	27	1	0

Hemoglobin	Specimen: U24	Negative	Trace	Small	Moderate	Large
[01]	Bayer Multistix (Manual)	0	0	1	43	21
[02]	Bayer Clinitek	0	1	4	175	0
[03]	Bayer Atlas	0	0	0	71	0
[10]	Roche Manual (Chemstrip)	0	2	0	15	13
[11]	Roche UA Reader/Super UA	0	0	0	9	20
[12]	Roche Mini UA Reader	0	0	0	1	4
[13]	Roche Chemstrip Criterion	0	0	1	4	21
[17]	IRIS BMC UA Reader	0	0	1	6	0
[18]	IRIS Clinitek Reader	0	0	0	2	1
[20]	Bayer Clinitest	0	0	0	2	1
[25]	Hitachi/Roche Analyzer	0	0	0	1	2
[27]	Bayer Hemastix	0	0	0	5	1
[29]	Chimera UA Perfect	0	0	0	1	0
[30]	Other Chemistry Analyzer	0	1	1	3	7
[90]	Method not listed	0	1	1	5	0
[99]	All methods	0	5	9	343	91

Hemoglobin	Specimen: U25	Negative	Trace	Small	Moderate	Large
[01]	Bayer Multistix (Manual)	61	3	0	0	0
[02]	Bayer Clinitek	165	14	1	0	0
[03]	Bayer Atlas	71	0	0	0	0
[10]	Roche Manual (Chemstrip)	28	0	0	0	0
[11]	Roche UA Reader/Super UA	28	1	0	0	0
[12]	Roche Mini UA Reader	5	0	0	0	0
[13]	Roche Chemstrip Criterion	26	0	0	0	0
[17]	IRIS BMC UA Reader	7	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0
[20]	Bayer Clinitest	3	0	0	0	0
[25]	Hitachi/Roche Analyzer	3	0	0	0	0
[27]	Bayer Hemastix	2	0	0	0	0
[29]	Chimera UA Perfect	1	0	0	0	0
[30]	Other Chemistry Analyzer	12	0	0	0	0
[90]	Method not listed	7	0	0	0	0
[99]	All methods	422	18	1	0	0

New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004

Ketones - Specimen: U21	Negative	Trace	Small	Moderate	Large	Positive
[01] Bayer Multistix (Manual)	65	0	0	0	0	0
[02] Bayer Clinitek	179	0	0	0	0	0
[03] Bayer Atlas	72	0	0	0	0	0
[10] Roche Manual (Chemstrip)	28	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0
[17] IRIS BMC UA Reader	7	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0
[19] Bayer Acetest	5	0	0	0	0	0
[20] Bayer Clinitest	3	0	0	0	0	0
[25] Hitachi/Roche Analyzer	3	0	0	0	0	0
[30] Other Chemistry Analyzer	12	0	0	0	0	0
[90] Method not listed	7	0	0	0	0	0
[99] All methods	444	0	0	0	0	0

Ketones - Specimen: U22	Negative	Trace	Small	Moderate	Large	Positive
[01] Bayer Multistix (Manual)	65	0	0	0	0	0
[02] Bayer Clinitek	178	0	1	0	0	0
[03] Bayer Atlas	72	0	0	0	0	0
[10] Roche Manual (Chemstrip)	28	0	0	0	0	0
[11] Roche UA Reader/Super UA	29	0	0	0	0	0
[12] Roche Mini UA Reader	5	0	0	0	0	0
[13] Roche Chemstrip Criterion	26	0	0	0	0	0
[17] IRIS BMC UA Reader	7	0	0	0	0	0
[18] IRIS Clinitek Reader	3	0	0	0	0	0
[19] Bayer Acetest	5	0	0	0	0	0
[20] Bayer Clinitest	3	0	0	0	0	0
[25] Hitachi/Roche Analyzer	3	0	0	0	0	0
[30] Other Chemistry Analyzer	12	0	0	0	0	0
[90] Method not listed	7	0	0	0	0	0
[99] All methods	443	0	1	0	0	0

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Ketones	- Specimen: U23	Negative	Trace	Small	Moderate	Large	Positive
[01]	Bayer Multistix (Manual)	0	0	2	53	10	0
[02]	Bayer Clinitek	0	0	67	111	1	0
[03]	Bayer Atlas	0	0	12	61	0	0
[10]	Roche Manual (Chemstrip)	0	0	3	28	0	0
[11]	Roche UA Reader/Super UA	0	0	0	27	2	0
[12]	Roche Mini UA Reader	0	0	1	4	0	0
[13]	Roche Chemstrip Criterion	0	0	4	22	0	0
[17]	IRIS BMC UA Reader	0	0	1	6	0	0
[18]	IRIS Clinitek Reader	0	0	0	3	0	0
[19]	Bayer Acetest	0	1	4	28	0	28
[20]	Bayer Clinitest	0	0	1	2	0	0
[25]	Hitachi/Roche Analyzer	0	0	1	2	0	0
[26]	Bayer Diastix	0	0	0	3	0	0
[30]	Other Chemistry Analyzer	0	1	1	10	0	0
[90]	Method not listed	0	1	1	6	0	0
[99]	All methods	0	3	98	366	13	28

Ketones	- Specimen: U24	Negative	Trace	Small	Moderate	Large	Positive
[01]	Bayer Multistix (Manual)	65	0	0	0	0	0
[02]	Bayer Clinitek	178	0	0	0	1	0
[03]	Bayer Atlas	72	0	0	0	0	0
[10]	Roche Manual (Chemstrip)	27	0	1	0	0	0
[11]	Roche UA Reader/Super UA	29	0	0	0	0	0
[12]	Roche Mini UA Reader	5	0	0	0	0	0
[13]	Roche Chemstrip Criterion	26	0	0	0	0	0
[17]	IRIS BMC UA Reader	7	0	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0	0
[19]	Bayer Acetest	5	0	0	0	0	0
[20]	Bayer Clinitest	3	0	0	0	0	0
[25]	Hitachi/Roche Analyzer	3	0	0	0	0	0
[30]	Other Chemistry Analyzer	12	0	0	0	0	0
[90]	Method not listed	6	0	0	0	0	0
[99]	All methods	441	0	1	0	1	0

Ketones	- Specimen: U25	Negative	Trace	Small	Moderate	Large	Positive
[01]	Bayer Multistix (Manual)	1	0	50	14	0	0
[02]	Bayer Clinitek	1	18	159	1	0	0
[03]	Bayer Atlas	1	1	71	0	0	0
[10]	Roche Manual (Chemstrip)	0	1	27	3	0	0
[11]	Roche UA Reader/Super UA	0	1	26	2	0	0
[12]	Roche Mini UA Reader	0	0	4	1	0	0
[13]	Roche Chemstrip Criterion	0	1	25	0	0	0
[17]	IRIS BMC UA Reader	0	3	4	0	0	0
[18]	IRIS Clinitek Reader	0	0	3	0	0	0
[19]	Bayer Acetest	0	2	29	1	0	28
[20]	Bayer Clinitest	0	1	2	0	0	0
[25]	Hitachi/Roche Analyzer	0	0	3	0	0	0
[26]	Bayer Diastix	0	1	2	0	0	0
[30]	Other Chemistry Analyzer	0	0	11	1	0	0
[90]	Method not listed	0	0	8	0	0	0
[95]	Method not specified	0	0	1	0	0	0
[99]	All methods	3	29	425	23	0	28

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Protein	- Specimen: U21	Negative	Trace	30 mg/dL	100 mg/dL	300 mg/dL	500 mg/dL	Positive
[01]	Bayer Multistix (Manual)	66	0	0	0	0	0	0
[02]	Bayer Clinitek	179	0	0	1	0	0	0
[03]	Bayer Atlas	71	0	0	0	0	0	0
[10]	Roche Manual (Chemstrip)	27	0	0	0	0	0	0
[11]	Roche UA Reader/Super UA	29	0	0	0	0	0	0
[12]	Roche Mini UA Reader	5	0	0	0	0	0	0
[13]	Roche Chemstrip Criterion	26	0	0	0	0	0	0
[17]	IRIS BMC UA Reader	7	0	0	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0	0	0
[20]	Bayer Clinitest	3	0	0	0	0	0	0
[22]	Sulfosalicylic Acid	15	0	0	0	0	0	0
[25]	Hitachi/Roche Analyzer	3	0	0	0	0	0	0
[30]	Other Chemistry Analyzer	11	1	0	0	0	0	0
[90]	Method not listed	11	0	0	0	0	0	0
[99]	All methods	456	1	0	1	0	0	0

Protein	- Specimen: U22	Negative	Trace	30 mg/dL	100 mg/dL	300 mg/dL	500 mg/dL	Positive
[01]	Bayer Multistix (Manual)	0	22	44	0	0	0	0
[02]	Bayer Clinitek	3	125	52	0	0	0	0
[03]	Bayer Atlas	0	44	28	0	0	0	0
[10]	Roche Manual (Chemstrip)	1	19	9	0	0	0	0
[11]	Roche UA Reader/Super UA	5	8	16	0	0	0	0
[12]	Roche Mini UA Reader	0	0	5	0	0	0	0
[13]	Roche Chemstrip Criterion	3	22	1	0	0	0	0
[17]	IRIS BMC UA Reader	2	5	0	0	0	0	0
[18]	IRIS Clinitek Reader	1	2	0	0	0	0	0
[20]	Bayer Clinitest	0	0	3	0	0	0	0
[22]	Sulfosalicylic Acid	3	24	14	1	0	0	56
[25]	Hitachi/Roche Analyzer	0	1	2	0	0	0	0
[30]	Other Chemistry Analyzer	1	5	5	2	0	0	0
[90]	Method not listed	2	8	2	0	0	0	0
[95]	Method not specified	0	0	0	0	0	0	1
[99]	All methods	21	285	181	3	0	0	57

**New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004**

Protein	- Specimen: U23	Negative	Trace	30 mg/dL	100 mg/dL	300 mg/dL	500 mg/dL	Positive
[01]	Bayer Multistix (Manual)	0	0	0	2	50	14	0
[02]	Bayer Clinitek	0	0	0	14	157	9	0
[03]	Bayer Atlas	0	0	0	0	68	4	0
[10]	Roche Manual (Chemstrip)	0	0	0	0	17	12	0
[11]	Roche UA Reader/Super UA	0	0	0	17	0	12	0
[12]	Roche Mini UA Reader	0	0	0	2	0	3	0
[13]	Roche Chemstrip Criterion	0	0	0	3	7	16	0
[17]	IRIS BMC UA Reader	0	0	0	2	4	1	0
[18]	IRIS Clinitek Reader	0	0	0	0	3	0	0
[20]	Bayer Clinitest	0	0	0	0	3	0	0
[22]	Sulfosalicylic Acid	0	0	0	5	32	12	68
[25]	Hitachi/Roche Analyzer	0	0	0	2	0	2	0
[30]	Other Chemistry Analyzer	0	0	0	3	5	5	0
[90]	Method not listed	0	0	0	3	6	3	0
[95]	Method not specified	0	0	0	0	0	0	1
[99]	All methods	0	0	0	53	352	93	69

Protein	- Specimen: U24	Negative	Trace	30 mg/dL	100 mg/dL	300 mg/dL	500 mg/dL	Positive
[01]	Bayer Multistix (Manual)	0	0	1	55	10	0	0
[02]	Bayer Clinitek	0	0	0	150	28	2	0
[03]	Bayer Atlas	0	0	0	72	0	0	0
[10]	Roche Manual (Chemstrip)	0	0	2	26	1	0	0
[11]	Roche UA Reader/Super UA	0	0	6	23	0	0	0
[12]	Roche Mini UA Reader	0	0	0	5	0	0	0
[13]	Roche Chemstrip Criterion	0	0	2	24	0	0	0
[17]	IRIS BMC UA Reader	0	0	2	5	0	0	0
[18]	IRIS Clinitek Reader	0	0	0	3	0	0	0
[20]	Bayer Clinitest	0	0	0	3	0	0	0
[22]	Sulfosalicylic Acid	0	0	4	38	5	0	64
[25]	Hitachi/Roche Analyzer	0	0	1	3	0	0	0
[30]	Other Chemistry Analyzer	0	0	2	11	0	0	0
[90]	Method not listed	0	0	1	11	0	0	0
[95]	Method not specified	0	0	0	0	0	0	1
[99]	All methods	0	0	21	429	44	2	65

Protein	- Specimen: U25	Negative	Trace	30 mg/dL	100 mg/dL	300 mg/dL	500 mg/dL	Positive
[01]	Bayer Multistix (Manual)	66	0	0	0	0	0	0
[02]	Bayer Clinitek	180	0	0	0	0	0	0
[03]	Bayer Atlas	71	0	0	0	0	0	0
[10]	Roche Manual (Chemstrip)	27	0	0	0	0	0	0
[11]	Roche UA Reader/Super UA	29	0	0	0	0	0	0
[12]	Roche Mini UA Reader	5	0	0	0	0	0	0
[13]	Roche Chemstrip Criterion	26	0	0	0	0	0	0
[17]	IRIS BMC UA Reader	7	0	0	0	0	0	0
[18]	IRIS Clinitek Reader	3	0	0	0	0	0	0
[20]	Bayer Clinitest	3	0	0	0	0	0	0
[22]	Sulfosalicylic Acid	15	0	0	0	0	0	0
[25]	Hitachi/Roche Analyzer	3	0	0	0	0	0	0
[30]	Other Chemistry Analyzer	13	0	0	0	0	0	0
[90]	Method not listed	11	0	0	0	0	0	0
[99]	All methods	459	0	0	0	0	0	0

New York State Department of Health - Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing - October 25, 2004

Qualitative hCG - Specimen: U21

hCG Content in Sample: 0 mIU/mL

Number of Laboratories Reporting Results

0 = Positive

0 = Indeterminate

382 = Negative

Qualitative hCG - Specimen: U22

hCG Content in Sample: 0 mIU/mL

Number of Laboratories Reporting Results

0 = Positive

0 = Indeterminate

382 = Negative

Qualitative hCG - Specimen: U23

hCG Content in Sample: 65 mIU/mL

Number of Laboratories Reporting Results

374 = Positive

3 = Indeterminate

4 = Negative

Qualitative hCG - Specimen: U24

hCG Content in Sample: 65 mIU/mL

Number of Laboratories Reporting Results

375 = Positive

1 = Indeterminate

5 = Negative

Qualitative hCG - Specimen: U25

hCG Content in Sample: 0 mIU/mL

Number of Laboratories Reporting Results

0 = Positive

0 = Indeterminate

382 = Negative