



Department of Health

ANDREW M. CUOMO
Governor

HOWARD A. ZUCKER, M.D., J.D.
Commissioner

SALLY DRESLIN, M.S., R.N.
Executive Deputy Commissioner

Statistical Report Urinalysis/Urine Pregnancy Testing Mail out date: October 26, 2015

This report summarizes data from the urinalysis/urine pregnancy proficiency testing of October 26, 2015.

If you have any questions regarding this report, please contact Regina Saulsbery at 518-474-6156.
You may also contact us by e-mail: clinchem@wadsworth.org

Results are shown for each analyte arranged by instrument and reagent system. Results for specific gravity and pH are listed as the mean $\pm$ 1SD for each of the instruments and reagents. All other results are shown as the number of participant results for that instrument or reagent.

If your laboratory participated in the educational quantitative urine clinical chemistry proficiency test, the evaluation report and summary are separate and included in this mailing.

The use of brand and/or trade names in this report does not constitute an endorsement of the products on the part of the Wadsworth Center or the New York State Department of Health.

NOTE: An unusually wide dispersion of results for glucose on sample U40 were observed. Since the participant consensus was less than 80%, the sample was not evaluated (not gradable).

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Specific Gravity:

Specimen: U36	Specimen: U37	Specimen: U38	Specimen: U39	Specimen: U40	Number	Instrument or Reagent System
1.014 ± 0.003	1.017 ± 0.003	1.019 ± 0.005	1.018 ± 0.003	1.021 ± 0.004	n = 350	All Methods & Instruments
1.015	1.016	1.019	1.018	1.021	n = 350	Median value, All Methods & Instruments
<Instruments>						
1.010 ± 0.003	1.014 ± 0.003	1.014 ± 0.004	1.015 ± 0.004	1.018 ± 0.003	n = 14	Visual examination
1.015 ± 0.002	1.017 ± 0.002	1.019 ± 0.002	1.018 ± 0.002	1.022 ± 0.002	n = 88	Siemens Atlas
1.018 ± 0.003	1.020 ± 0.003	1.024 ± 0.003	1.021 ± 0.002	1.023 ± 0.003	n = 46	Siemens Clinitek Status
1.015 ± 0.004	1.018 ± 0.006	1.026 ± 0.005	1.021 ± 0.006	1.023 ± 0.005	n = 4	Siemens Clinitek 50
1.015 ± 0.004	1.018 ± 0.004	1.022 ± 0.003	1.020 ± 0.003	1.022 ± 0.003	n = 23	Siemens Clinitek 500
1.015 ± 0.003	1.018 ± 0.003	1.022 ± 0.003	1.019 ± 0.003	1.023 ± 0.003	n = 67	Siemens Clinitek Advantus
1.009 ± 0.002	1.015 ± 0.000	1.010 ± 0.000	1.013 ± 0.003	1.015 ± 0.000	n = 4	Roche Chemstrip Criterion
1.012 ± 0.002	1.015 ± 0.001	1.015 ± 0.003	1.016 ± 0.001	1.019 ± 0.003	n = 19	Roche Urisys
1.010 ± 0.003	1.015 ± 0.001	1.012 ± 0.003	1.015 ± 0.001	1.015 ± 0.001	n = 11	Roche Cobas U411
1.010 ± 0.001	1.014 ± 0.001	1.014 ± 0.002	1.016 ± 0.002	1.019 ± 0.002	n = 21	IRIS AUTION MAX
1.010 ± 0.001	1.014 ± 0.001	1.014 ± 0.001	1.016 ± 0.001	1.019 ± 0.001	n = 33	IRIS iChem Velocity
1.014 ± 0.001	1.015 ± 0.001	1.017 ± 0.001	1.017 ± 0.001	1.021 ± 0.001	n = 13	Refractometer
1.013 ± 0.005	1.016 ± 0.001	1.018 ± 0.002	1.017 ± 0.004	1.021 ± 0.003	n = 5	Other
<Reagents>						
1.015 ± 0.004	1.018 ± 0.004	1.022 ± 0.004	1.019 ± 0.004	1.022 ± 0.003	n = 160	Siemens Reagent Strips
1.015 ± 0.002	1.016 ± 0.002	1.019 ± 0.002	1.018 ± 0.002	1.022 ± 0.002	n = 80	Siemens Atlas Reagent
1.010 ± 0.003	1.015 ± 0.001	1.012 ± 0.003	1.015 ± 0.002	1.015 ± 0.002	n = 27	Roche Chemstrips
1.014 ± 0.001	1.015 ± 0.000	1.017 ± 0.001	1.017 ± 0.001	1.022 ± 0.001	n = 9	Roche Urisys Cassette
1.010 ± 0.001	1.014 ± 0.001	1.014 ± 0.002	1.016 ± 0.002	1.019 ± 0.002	n = 15	ARKRAY Aution Sticks
1.010 ± 0.000	1.014 ± 0.001	1.014 ± 0.001	1.016 ± 0.000	1.018 ± 0.001	n = 34	IRIS Reagent Strips
1.013 ± 0.002	1.015 ± 0.001	1.017 ± 0.003	1.017 ± 0.001	1.020 ± 0.005	n = 18	None, no reagent or strip used
1.014 ± 0.005	1.017 ± 0.002	1.019 ± 0.003	1.017 ± 0.004	1.022 ± 0.003	n = 6	Other

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

pH:

Specimen: U36	Specimen: U37	Specimen: U38	Specimen: U39	Specimen: U40	Number	Instrument or Reagent System
7.13 ± 0.23	6.82 ± 0.29	7.51 ± 0.29	6.93 ± 0.18	6.88 ± 0.25	n = 354	All Methods & Instruments
7.00	7.00	7.50	7.00	7.00	n = 354	Median value, All Methods & Instruments
<Instruments>						
7.35 ± 0.29	6.85 ± 0.37	7.68 ± 0.24	6.90 ± 0.26	6.93 ± 0.18	n = 20	Visual examination
7.05 ± 0.14	6.80 ± 0.25	7.49 ± 0.05	6.97 ± 0.12	6.89 ± 0.21	n = 88	Siemens Atlas
7.07 ± 0.17	6.99 ± 0.07	7.50 ± 0.24	7.00 ± 0.00	7.00 ± 0.00	n = 51	Siemens Clinitek Status
7.00 ± 0.00	6.75 ± 0.29	7.50 ± 0.00	6.88 ± 0.25	7.00 ± 0.00	n = 4	Siemens Clinitek 50
7.29 ± 0.25	6.88 ± 0.22	7.52 ± 0.38	6.92 ± 0.19	6.98 ± 0.10	n = 24	Siemens Clinitek 500
7.33 ± 0.24	6.83 ± 0.31	7.54 ± 0.20	6.94 ± 0.16	6.92 ± 0.20	n = 69	Siemens Clinitek Advantus
7.00 ± 0.00	6.30 ± 0.76	7.80 ± 0.45	6.80 ± 0.27	6.70 ± 0.27	n = 5	Roche Chemstrip Criterion
7.00 ± 0.00	6.76 ± 0.26	7.61 ± 0.49	6.84 ± 0.24	6.82 ± 0.25	n = 19	Roche Urisys
7.00 ± 0.00	6.59 ± 0.20	7.73 ± 0.47	6.77 ± 0.26	6.68 ± 0.25	n = 11	Roche Cobas U411
7.00 ± 0.00	6.52 ± 0.19	7.40 ± 0.20	6.64 ± 0.23	6.60 ± 0.20	n = 21	IRIS AUTION MAX
7.00 ± 0.00	6.91 ± 0.29	7.27 ± 0.45	7.00 ± 0.00	6.85 ± 0.36	n = 33	IRIS iChem Velocity
7.40 ± 0.22	6.90 ± 0.22	7.70 ± 0.27	6.90 ± 0.22	6.90 ± 0.22	n = 5	Other
<Reagents>						
7.23 ± 0.26	6.89 ± 0.25	7.53 ± 0.25	6.95 ± 0.14	6.96 ± 0.15	n = 175	Siemens Reagent Strips
7.04 ± 0.13	6.79 ± 0.25	7.50 ± 0.00	6.97 ± 0.11	6.89 ± 0.21	n = 78	Siemens Atlas Reagent
7.00 ± 0.00	6.52 ± 0.37	7.79 ± 0.40	6.75 ± 0.25	6.70 ± 0.25	n = 28	Roche Chemstrips
7.00 ± 0.00	6.94 ± 0.17	7.33 ± 0.50	7.00 ± 0.00	6.94 ± 0.17	n = 9	Roche Urisys Cassette
7.00 ± 0.00	6.53 ± 0.21	7.39 ± 0.21	6.64 ± 0.23	6.58 ± 0.19	n = 18	ARKRAY Aution Sticks
7.00 ± 0.00	6.88 ± 0.30	7.29 ± 0.44	6.97 ± 0.12	6.83 ± 0.36	n = 36	IRIS Reagent Strips
7.36 ± 0.24	6.86 ± 0.24	7.64 ± 0.24	6.79 ± 0.39	6.86 ± 0.24	n = 7	Other

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Bilirubin Specimen: U36

	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>					
Visual examination	27	0	0	0	0
Siemens Atlas	88	0	0	0	0
Siemens Clinitek Status	49	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0
Siemens Clinitek 500	22	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0
Roche Urisys	18	0	0	0	0
Roche Cobas U411	11	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0
Other	5	0	0	0	0
<Reagents>					
Siemens Reagent Strips	169	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0
Roche Chemstrips	28	0	0	0	0
Roche Urisys Cassette	8	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0
Siemens Ictotest	9	0	0	0	0
Siemens Hemastix	1	0	0	0	0
Siemens Uristix	1	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0
Other	7	0	0	0	0
All Instruments and Reagents	354	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Bilirubin Specimen: U37

	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>					
Visual examination	2	1	0	22	58
Siemens Atlas	0	0	0	86	2
Siemens Clinitek Status	0	0	0	49	0
Siemens Clinitek 50	0	0	0	4	0
Siemens Clinitek 100	0	0	0	2	0
Siemens Clinitek 500	0	0	3	18	1
Siemens Clinitek Advantus	0	0	15	53	1
Roche Chemstrip Criterion	0	0	5	0	0
Roche Urisys	0	1	15	2	0
Roche Cobas U411	0	0	10	1	0
IRIS AUTION MAX	2	1	9	9	0
IRIS iChem Velocity	33	0	0	0	0
Other	5	0	1	0	0
<Reagents>					
Siemens Reagent Strips	0	0	17	149	4
Siemens Atlas Reagent	0	0	0	74	2
Roche Chemstrips	0	1	23	5	0
Roche Urisys Cassette	0	0	7	1	0
ARKRAY Aution Sticks	3	1	7	6	0
IRIS Reagent Strips	32	0	2	3	0
Siemens Ictotest	2	0	0	4	54
Siemens Hemastix	0	0	0	1	0
Siemens Uristix	0	0	0	1	0
Biorex Ictero-Check	0	1	0	1	2
Other	5	0	2	1	0
All Instruments and Reagents	42	3	58	246	62

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Bilirubin Specimen: U38

	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>					
Visual examination	27	0	0	0	0
Siemens Atlas	88	0	0	0	0
Siemens Clinitek Status	49	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0
Siemens Clinitek 500	22	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0
Roche Urisys	18	0	0	0	0
Roche Cobas U411	11	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0
Other	5	0	0	0	0
<Reagents>					
Siemens Reagent Strips	169	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0
Roche Chemstrips	28	0	0	0	0
Roche Urisys Cassette	8	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0
Siemens Ictotest	9	0	0	0	0
Siemens Hemastix	1	0	0	0	0
Siemens Uristix	1	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0
Other	7	0	0	0	0
All Instruments and Reagents	354	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Bilirubin Specimen: U39

	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>					
Visual examination	27	0	0	0	0
Siemens Atlas	88	0	0	0	0
Siemens Clinitek Status	49	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0
Siemens Clinitek 500	22	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0
Roche Urisys	18	0	0	0	0
Roche Cobas U411	11	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0
Other	5	0	0	0	0
<Reagents>					
Siemens Reagent Strips	169	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0
Roche Chemstrips	28	0	0	0	0
Roche Urisys Cassette	8	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0
Siemens Ictotest	9	0	0	0	0
Siemens Hemastix	1	0	0	0	0
Siemens Uristix	1	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0
Other	7	0	0	0	0
All Instruments and Reagents	354	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Bilirubin **Specimen: U40**

	Negative	Small 1+ 1mg/dL	Moderate 2+ 3mg/dL	Large 3+ 6mg/dL	Positive
<Instruments>					
Visual examination	27	0	0	0	0
Siemens Atlas	88	0	0	0	0
Siemens Clinitek Status	49	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0
Siemens Clinitek 500	22	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0
Roche Urisys	18	0	0	0	0
Roche Cobas U411	11	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0
Other	5	0	0	0	0
<Reagents>					
Siemens Reagent Strips	169	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0
Roche Chemstrips	28	0	0	0	0
Roche Urisys Cassette	8	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0
Siemens Ictotest	9	0	0	0	0
Siemens Hemastix	1	0	0	0	0
Siemens Uristix	1	0	0	0	0
Biorex Ictero-Check	1	0	0	0	0
Other	7	0	0	0	0
All Instruments and Reagents	354	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Glucose	Specimen: U36							
	Negative	Trace	50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination	1	2		10	62	9	5	1
Siemens Atlas	0	0		2	86	0	1	0
Siemens Clinitek Status	2	1		5	38	0	4	1
Siemens Clinitek 50	0	1		0	2	0	1	0
Siemens Clinitek 100	0	0		0	1	0	1	0
Siemens Clinitek 500	0	0		2	23	0	0	0
Siemens Clinitek Advantus	1	2		6	58	0	2	0
Roche Chemstrip Criterion	0	0		2	1	0	2	0
Roche Urisys	0	0		4	1	0	14	0
Roche Cobas U411	0	0		10	0	0	1	0
IRIS AUTION MAX	0	0		0	19	2	0	0
IRIS iChem Velocity	0	0		0	33	0	0	0
Other Chemistry Analyzer	0	0		0	0	1	0	0
Other	0	0		0	4	0	0	0
<Reagents>								
Siemens Reagent Strips	3	4		16	142	0	10	2
Siemens Atlas Reagent	0	0		2	74	0	1	0
Roche Chemstrips	0	0		16	6	0	10	0
Roche Urisys Cassette	0	0		0	0	0	9	0
ARKRAY Aution Sticks	0	0		0	16	2	0	0
IRIS Reagent Strips	0	0		0	36	0	0	0
Bayer Clinitest	1	2		6	43	10	1	0
Bayer Diastix	0	0		0	1	0	0	0
Siemens Uristix	0	0		0	4	0	0	0
Other	0	0		1	6	0	0	0
All Instruments and Reagents	4	6		41	328	12	31	2

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Glucose		Specimen: U37							
		Negative	Trace	50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>									
Visual examination	64	0		0		0	0	0	0
Siemens Atlas	89	0		0		0	0	0	0
Siemens Clinitek Status	51	0		0		0	0	0	0
Siemens Clinitek 50	4	0		0		0	0	0	0
Siemens Clinitek 100	2	0		0		0	0	0	0
Siemens Clinitek 500	23	1		0		0	0	0	0
Siemens Clinitek Advantus	68	1		0		0	0	0	0
Roche Chemstrip Criterion	5	0		0		0	0	0	0
Roche Urisys	19	0		0		0	0	0	0
Roche Cobas U411	11	0		0		0	0	0	0
IRIS AUTION MAX	21	0		0		0	0	0	0
IRIS iChem Velocity	33	0		0		0	0	0	0
Other Chemistry Analyzer	0	0		0		0	0	0	0
Other	4	0		0		0	0	0	0
<Reagents>									
Siemens Reagent Strips	175	2		0		0	0	0	0
Siemens Atlas Reagent	77	0		0		0	0	0	0
Roche Chemstrips	28	0		0		0	0	0	0
Roche Urisys Cassette	9	0		0		0	0	0	0
ARKRAY Aution Sticks	18	0		0		0	0	0	0
IRIS Reagent Strips	36	0		0		0	0	0	0
Bayer Clinitest	40	0		0		0	0	0	0
Bayer Diastix	0	0		0		0	0	0	0
Siemens Uristix	4	0		0		0	0	0	0
Other	7	0		0		0	0	0	0
All Instruments and Reagents	394	2		0		0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Glucose		Specimen: U38									
		Negative	Trace	50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL		
<Instruments>											
Visual examination		1	2		53	33	0	1		0	
Siemens Atlas		0	0		68	21	0	0		0	
Siemens Clinitek Status		0	0		28	22	0	0		1	
Siemens Clinitek 50		0	0		2	2	0	0		0	
Siemens Clinitek 100		0	0		0	2	0	0		0	
Siemens Clinitek 500		0	0		16	9	0	0		0	
Siemens Clinitek Advantus		0	2		38	27	0	2		0	
Roche Chemstrip Criterion		0	0		5	0	0	0		0	
Roche Urisys		0	1		16	0	0	2		0	
Roche Cobas U411		0	0		11	0	0	0		0	
IRIS AUTION MAX		0	0		16	3	2	0		0	
IRIS iChem Velocity		0	0		3	30	0	0		0	
Other Chemistry Analyzer		0	0		1	0	0	0		0	
Other		0	0		0	4	0	0		0	
<Reagents>											
Siemens Reagent Strips		0	3		100	70	0	3		1	
Siemens Atlas Reagent		0	0		58	19	0	0		0	
Roche Chemstrips		0	1		29	2	0	0		0	
Roche Urisys Cassette		0	0		7	0	0	2		0	
ARKRAY Aution Sticks		0	0		13	3	2	0		0	
IRIS Reagent Strips		0	0		6	30	0	0		0	
Bayer Clinitest		1	1		38	23	0	0		0	
Bayer Diastix		0	0		1	0	0	0		0	
Siemens Uristix		0	0		4	0	0	0		0	
Other		0	0		1	6	0	0		0	
All Instruments and Reagents		1	5		257	153	2	5		1	

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Glucose Specimen: U39		Negative	Trace	50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>									
Visual examination	64	0		0		0	0	0	0
Siemens Atlas	89	0		0		0	0	0	0
Siemens Clinitek Status	51	0		0		0	0	0	0
Siemens Clinitek 50	4	0		0		0	0	0	0
Siemens Clinitek 100	2	0		0		0	0	0	0
Siemens Clinitek 500	24	0		0		0	0	0	0
Siemens Clinitek Advantus	69	0		0		0	0	0	0
Roche Chemstrip Criterion	5	0		0		0	0	0	0
Roche Urisys	19	0		0		0	0	0	0
Roche Cobas U411	11	0		0		0	0	0	0
IRIS AUTION MAX	21	0		0		0	0	0	0
IRIS iChem Velocity	33	0		0		0	0	0	0
Other Chemistry Analyzer	0	0		0		0	0	0	0
Other	4	0		0		0	0	0	0
<Reagents>									
Siemens Reagent Strips	177	0		0		0	0	0	0
Siemens Atlas Reagent	77	0		0		0	0	0	0
Roche Chemstrips	28	0		0		0	0	0	0
Roche Urisys Cassette	9	0		0		0	0	0	0
ARKRAY Aution Sticks	18	0		0		0	0	0	0
IRIS Reagent Strips	36	0		0		0	0	0	0
Bayer Clinitest	40	0		0		0	0	0	0
Bayer Diastix	0	0		0		0	0	0	0
Siemens Uristix	4	0		0		0	0	0	0
Other	7	0		0		0	0	0	0
All Instruments and Reagents	396	0		0		0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Glucose	Specimen: U40							
	Negative	Trace	50-100mg/dL	1/4% 150-300mg/dL	1/2% 500mg/dL	3/4% 750mg/dL	1% >=1000mg/dL	2% >=2000mg/dL
<Instruments>								
Visual examination	13	16		44	11	0	0	0
Siemens Atlas	0	1		87	1	0	0	0
Siemens Clinitek Status	11	40		0	0	0	0	0
Siemens Clinitek 50	1	3		0	0	0	0	0
Siemens Clinitek 100	0	2		0	0	0	0	0
Siemens Clinitek 500	16	7		1	0	0	0	0
Siemens Clinitek Advantus	43	26		0	0	0	0	0
Roche Chemstrip Criterion	0	0		5	0	0	0	0
Roche Urisys	0	3		15	1	0	0	0
Roche Cobas U411	0	3		8	0	0	0	0
IRIS AUTION MAX	0	0		17	2	2	0	0
IRIS iChem Velocity	0	0		16	17	0	0	0
Other Chemistry Analyzer	0	0		1	0	0	0	0
Other	0	0		3	1	0	0	0
<Reagents>								
Siemens Reagent Strips	80	85		12	0	0	0	0
Siemens Atlas Reagent	0	1		75	1	0	0	0
Roche Chemstrips	1	6		21	1	0	0	0
Roche Urisys Cassette	0	0		9	0	0	0	0
ARKRAY Aution Sticks	0	0		15	1	2	0	0
IRIS Reagent Strips	0	0		18	18	0	0	0
Bayer Clinitest	2	4		43	11	0	0	0
Bayer Diastix	0	0		0	0	0	0	0
Siemens Uristix	0	4		0	0	0	0	0
Other	1	1		4	1	0	0	0
All Instruments and Reagents	84	101		197	33	2	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Hemoglobin Specimen: U36

	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	0	0	0	0	0	22
Siemens Atlas	0	0	0	0	0	88
Siemens Clinitek Status	0	0	0	6	0	46
Siemens Clinitek 50	0	1	0	0	0	3
Siemens Clinitek 100	0	0	0	0	0	2
Siemens Clinitek 500	0	0	0	1	0	23
Siemens Clinitek Advantus	1	0	0	1	0	66
Roche Chemstrip Criterion	0	0	0	0	0	5
Roche Urisys	0	0	0	1	0	18
Roche Cobas U411	0	0	0	0	0	11
IRIS AUTION MAX	0	0	0	5	0	16
IRIS iChem Velocity	0	0	0	8	0	25
Other	0	0	0	0	0	5
<Reagents>						
Siemens Reagent Strips	1	1	0	8	0	168
Siemens Atlas Reagent	0	0	0	0	0	76
Roche Chemstrips	0	0	0	1	0	28
Roche Urisys Cassette	0	0	0	0	0	9
ARKRAY Aution Sticks	0	0	0	5	0	13
IRIS Reagent Strips	0	0	0	8	0	28
Siemens Uristix	0	0	0	0	0	1
Other	0	0	0	0	0	7
All Instruments and Reagents	1	1	0	22	0	330

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Hemoglobin Specimen: U37

	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	19	1	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	51	1	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	23	0	0	1	0	0
Siemens Clinitek Advantus	66	0	1	0	0	1
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	18	1	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	32	0	1	0	0	0
Other	5	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	172	2	1	1	0	1
Siemens Atlas Reagent	76	0	0	0	0	0
Roche Chemstrips	27	1	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	18	0	0	0	0	0
IRIS Reagent Strips	35	0	1	0	0	0
Siemens Uristix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	345	3	2	1	0	1

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Hemoglobin Specimen: U38

<Instruments>

	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
Visual examination	0	0	0	0	0	22
Siemens Atlas	0	0	0	1	0	87
Siemens Clinitek Status	0	0	0	1	0	51
Siemens Clinitek 50	0	1	0	0	0	3
Siemens Clinitek 100	0	0	0	0	0	2
Siemens Clinitek 500	0	0	0	0	0	24
Siemens Clinitek Advantus	0	0	0	1	0	67
Roche Chemstrip Criterion	0	0	0	0	1	4
Roche Urisys	0	0	0	1	0	18
Roche Cobas U411	0	0	0	0	0	11
IRIS AUTION MAX	0	0	0	3	0	18
IRIS iChem Velocity	0	0	0	7	0	26
Other	0	0	0	0	0	5

<Reagents>

Siemens Reagent Strips	0	1	0	3	0	174
Siemens Atlas Reagent	0	0	0	0	0	76
Roche Chemstrips	0	0	0	1	1	27
Roche Urisys Cassette	0	0	0	0	0	9
ARKRAY Aution Sticks	0	0	0	3	0	15
IRIS Reagent Strips	0	0	0	7	0	29
Siemens Uristix	0	0	0	0	0	1
Other	0	0	0	0	0	7
All Instruments and Reagents	0	1	0	14	1	338

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Hemoglobin Specimen: U39

	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
<Instruments>						
Visual examination	20	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	48	4	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	23	1	0	0	0	0
Siemens Clinitek Advantus	68	0	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	19	0	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	32	0	1	0	0	0
Other	5	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	172	5	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0
Roche Chemstrips	28	0	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	18	0	0	0	0	0
IRIS Reagent Strips	35	0	1	0	0	0
Siemens Uristix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	346	5	1	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Hemoglobin Specimen: U40

<Instruments>

	Negative	Trace 5-10ery/uL	Small 1+ 25ery/uL	Moderate 2+ 50ery/uL	150ery/uL	Large 3+ 250ery/uL
Visual examination	20	0	0	0	0	0
Siemens Atlas	87	0	1	0	0	0
Siemens Clinitek Status	49	3	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	24	0	0	0	0	0
Siemens Clinitek Advantus	66	1	1	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	19	0	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0	0
Other	5	0	0	0	0	0

<Reagents>

Siemens Reagent Strips	172	4	1	0	0	0
Siemens Atlas Reagent	75	0	1	0	0	0
Roche Chemstrips	28	0	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	18	0	0	0	0	0
IRIS Reagent Strips	36	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	346	4	2	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Ketones **Specimen: U36**

	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>						
Visual examination	23	1	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	51	0	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	24	0	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	19	0	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0	0
Other	5	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	176	1	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0
Roche Chemstrips	28	0	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Siemens Ketostix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	355	1	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Ketones **Specimen: U37**

	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>						
Visual examination	24	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	51	0	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	24	0	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	19	0	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0	0
Other	5	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	177	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0
Roche Chemstrips	28	0	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Siemens Ketostix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	356	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Ketones Specimen: U38

<Instruments>

	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
Visual examination	24	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	51	0	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	24	0	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	19	0	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0	0
Other	5	0	0	0	0	0

<Reagents>

Siemens Reagent Strips	177	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0
Roche Chemstrips	28	0	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Siemens Ketostix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	356	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Ketones **Specimen: U39**

	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>						
Visual examination	24	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	51	0	0	0	0	0
Siemens Clinitek 50	4	0	0	0	0	0
Siemens Clinitek 100	2	0	0	0	0	0
Siemens Clinitek 500	24	0	0	0	0	0
Siemens Clinitek Advantus	69	0	0	0	0	0
Roche Chemstrip Criterion	5	0	0	0	0	0
Roche Urisys	19	0	0	0	0	0
Roche Cobas U411	11	0	0	0	0	0
IRIS AUTION MAX	21	0	0	0	0	0
IRIS iChem Velocity	33	0	0	0	0	0
Other	5	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	177	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0
Roche Chemstrips	28	0	0	0	0	0
Roche Urisys Cassette	9	0	0	0	0	0
ARKRAY Aution Sticks	17	0	0	0	0	0
IRIS Reagent Strips	37	0	0	0	0	0
Siemens Acetest	3	0	0	0	0	0
Siemens Uristix	1	0	0	0	0	0
Siemens Ketostix	1	0	0	0	0	0
Other	7	0	0	0	0	0
All Instruments and Reagents	356	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Ketones	Specimen: U40					
	Negative	Trace 5mg/dL	Small 1+ 15mg/dL	Moderate 2+ 40-50mg/dL	Large 3+ >=80mg/dL	Positive
<Instruments>						
Visual examination	0	0	2	7	29	2
Siemens Atlas	0	0	0	75	13	0
Siemens Clinitek Status	0	0	0	11	40	0
Siemens Clinitek 50	0	0	0	0	4	0
Siemens Clinitek 100	0	0	0	1	1	0
Siemens Clinitek 500	0	0	0	20	4	0
Siemens Clinitek Advantus	0	0	3	46	20	0
Roche Chemstrip Criterion	0	0	0	1	4	0
Roche Urisys	0	0	0	4	15	0
Roche Cobas U411	0	0	0	4	7	0
IRIS AUTION MAX	0	0	1	11	9	0
IRIS iChem Velocity	0	0	7	12	14	0
Other	0	0	0	1	4	0
<Reagents>						
Siemens Reagent Strips	0	0	2	89	86	0
Siemens Atlas Reagent	0	0	0	65	11	0
Roche Chemstrips	0	0	0	8	21	0
Roche Urisys Cassette	0	0	0	1	8	0
ARKRAY Aution Sticks	0	0	1	8	8	0
IRIS Reagent Strips	0	0	7	15	15	0
Siemens Acetest	0	0	2	5	8	2
Siemens Uristix	0	0	0	0	1	0
Siemens Ketostix	0	0	0	0	2	0
Other	0	0	1	2	4	0
All Instruments and Reagents	0	0	13	193	164	2

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Protein	Specimen: U36						
	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>							
Visual examination	0	0	0	4	40	4	7
Siemens Atlas	0	0	0	0	86	2	0
Siemens Clinitek Status	0	0	0	0	51	0	0
Siemens Clinitek 50	0	0	0	0	4	0	0
Siemens Clinitek 100	0	0	0	0	2	0	0
Siemens Clinitek 500	0	0	0	11	13	0	0
Siemens Clinitek Advantus	0	0	0	23	46	0	0
Roche Chemstrip Criterion	0	0	0	0	5	0	0
Roche Urisys	0	0	0	1	18	0	0
Roche Cobas U411	0	0	0	0	11	0	0
IRIS AUTION MAX	0	0	0	1	20	0	0
IRIS iChem Velocity	0	0	0	5	28	0	0
Hitachi/Roche Analyzer	0	0	0	0	0	0	0
Other	0	0	0	0	5	0	0
<Reagents>							
Siemens Reagent Strips	0	0	0	33	140	3	0
Siemens Atlas Reagent	0	0	0	0	77	1	0
Sulfosalicylic Acid	0	0	0	4	17	2	7
Roche Chemstrips	0	0	0	1	28	0	0
Roche Urisys Cassette	0	0	0	0	9	0	0
ARKRAY Aution Sticks	0	0	0	1	17	0	0
IRIS Reagent Strips	0	0	0	5	31	0	0
Siemens Uristix	0	0	0	0	4	0	0
Hitachi/Roche Reagents	0	0	0	0	0	0	0
Other	0	0	0	1	6	0	0
All Instruments and Reagents	0	0	0	45	329	6	7

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Protein	Specimen: U37							
	Negative	Trace	10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination	0	5		36	6	0	0	7
Siemens Atlas	0	0		81	7	0	0	0
Siemens Clinitek Status	0	0		5	46	0	0	0
Siemens Clinitek 50	0	0		0	4	0	0	0
Siemens Clinitek 100	0	0		0	2	0	0	0
Siemens Clinitek 500	0	0		24	0	0	0	0
Siemens Clinitek Advantus	0	1		68	0	0	0	0
Roche Chemstrip Criterion	0	1		4	0	0	0	0
Roche Urisys	0	2		6	11	0	0	0
Roche Cobas U411	0	0		6	5	0	0	0
IRIS AUTION MAX	0	0		20	1	0	0	0
IRIS iChem Velocity	0	0		31	2	0	0	0
Hitachi/Roche Analyzer	0	0		0	0	0	0	0
Other	0	0		4	1	0	0	0
<Reagents>								
Siemens Reagent Strips	0	1		118	56	0	0	0
Siemens Atlas Reagent	0	0		73	5	0	0	0
Sulfosalicylic Acid	0	5		15	3	0	0	7
Roche Chemstrips	0	2		18	9	0	0	0
Roche Urisys Cassette	0	1		1	7	0	0	0
ARKRAY Aution Sticks	0	0		18	0	0	0	0
IRIS Reagent Strips	0	0		33	3	0	0	0
Siemens Uristix	0	0		3	1	0	0	0
Hitachi/Roche Reagents	0	0		0	0	0	0	0
Other	0	0		6	1	0	0	0
All Instruments and Reagents	0	9		285	85	0	0	7

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Protein	Specimen: U38						
	Negative	Trace 10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>							
Visual examination	0	0	0	7	40	2	7
Siemens Atlas	0	0	0	0	87	1	0
Siemens Clinitek Status	0	0	0	0	51	0	0
Siemens Clinitek 50	0	0	0	0	4	0	0
Siemens Clinitek 100	0	0	0	0	2	0	0
Siemens Clinitek 500	0	0	0	11	13	0	0
Siemens Clinitek Advantus	0	0	0	21	48	0	0
Roche Chemstrip Criterion	0	0	0	0	5	0	0
Roche Urisys	0	0	0	1	18	0	0
Roche Cobas U411	0	0	0	1	10	0	0
IRIS AUTION MAX	0	0	0	3	18	0	0
IRIS iChem Velocity	0	0	0	6	27	0	0
Hitachi/Roche Analyzer	0	0	0	0	0	0	0
Other	0	0	0	0	5	0	0
<Reagents>							
Siemens Reagent Strips	0	0	0	32	143	1	0
Siemens Atlas Reagent	0	0	0	0	77	1	0
Sulfosalicylic Acid	0	0	0	4	19	1	7
Roche Chemstrips	0	0	0	4	25	0	0
Roche Urisys Cassette	0	0	0	0	9	0	0
ARKRAY Aution Sticks	0	0	0	3	15	0	0
IRIS Reagent Strips	0	0	0	6	30	0	0
Siemens Uristix	0	0	0	0	4	0	0
Hitachi/Roche Reagents	0	0	0	0	0	0	0
Other	0	0	0	1	6	0	0
All Instruments and Reagents	0	0	0	50	328	3	7

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Protein	Specimen: U39							
	Negative	Trace	10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination	0	0		7	32	7	1	7
Siemens Atlas	0	0		0	87	1	0	0
Siemens Clinitek Status	0	0		0	41	10	0	0
Siemens Clinitek 50	0	0		0	1	3	0	0
Siemens Clinitek 100	0	0		0	1	1	0	0
Siemens Clinitek 500	0	0		0	24	0	0	0
Siemens Clinitek Advantus	0	0		0	69	0	0	0
Roche Chemstrip Criterion	0	0		0	5	0	0	0
Roche Urisys	0	0		1	18	0	0	0
Roche Cobas U411	0	0		0	11	0	0	0
IRIS AUTION MAX	0	0		1	20	0	0	0
IRIS iChem Velocity	0	0		0	33	0	0	0
Hitachi/Roche Analyzer	0	0		0	0	0	0	0
Other	0	0		0	5	0	0	0
<Reagents>								
Siemens Reagent Strips	0	0		1	155	18	1	0
Siemens Atlas Reagent	0	0		0	77	1	0	0
Sulfosalicylic Acid	0	0		5	16	2	0	7
Roche Chemstrips	0	0		2	27	0	0	0
Roche Urisys Cassette	0	0		0	9	0	0	0
ARKRAY Aution Sticks	0	0		1	17	0	0	0
IRIS Reagent Strips	0	0		0	36	0	0	0
Siemens Uristix	0	0		0	4	0	0	0
Hitachi/Roche Reagents	0	0		0	0	0	0	0
Other	0	0		0	6	1	0	0
All Instruments and Reagents	0	0		9	347	22	1	7

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Protein	Specimen: U40							
	Negative	Trace	10-15mg/dL	1+ 25-30mg/dL	2+ 75-150mg/dL	3+ 300-600mg/dL	4+ >=1000mg/dL	Positive
<Instruments>								
Visual examination	0	0		4	29	14	1	7
Siemens Atlas	0	0		0	65	23	0	0
Siemens Clinitek Status	1	0		0	7	43	0	0
Siemens Clinitek 50	0	0		0	0	4	0	0
Siemens Clinitek 100	0	0		0	0	2	0	0
Siemens Clinitek 500	0	0		0	24	0	0	0
Siemens Clinitek Advantus	0	0		0	69	0	0	0
Roche Chemstrip Criterion	0	0		0	5	0	0	0
Roche Urisys	0	0		1	18	0	0	0
Roche Cobas U411	0	0		0	10	1	0	0
IRIS AUTION MAX	0	0		0	21	0	0	0
IRIS iChem Velocity	0	0		0	32	1	0	0
Hitachi/Roche Analyzer	0	0		0	0	0	0	0
Other	0	0		0	0	5	0	0
<Reagents>								
Siemens Reagent Strips	1	0		1	113	60	1	0
Siemens Atlas Reagent	0	0		0	57	21	0	0
Sulfosalicylic Acid	0	0		2	17	4	0	7
Roche Chemstrips	0	0		2	26	1	0	0
Roche Urisys Cassette	0	0		0	9	0	0	0
ARKRAY Aution Sticks	0	0		0	18	0	0	0
IRIS Reagent Strips	0	0		0	35	1	0	0
Siemens Uristix	0	0		0	3	1	0	0
Hitachi/Roche Reagents	0	0		0	0	0	0	0
Other	0	0		0	2	5	0	0
All Instruments and Reagents	1	0		5	280	93	1	7

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Qualitative hCG report:

Other	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	5	0	5	5
Indeterminant	0	0	0	0	0
Positive	5	0	5	0	0
Alere hCG Cassette	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	9	0	9	9
Indeterminant	0	0	0	0	0
Positive	9	0	9	0	0
Siemens Clinitek Status	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	18	0	18	18
Indeterminant	0	0	0	0	0
Positive	18	0	18	0	0
Beckman Coulter ICON II	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	2	0	2	2
Indeterminant	0	0	0	0	0
Positive	2	0	2	0	0
Beckman Coulter ICON 25	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	17	0	17	17
Indeterminant	0	0	0	0	0
Positive	17	0	17	0	0
BTNX Rapid Response	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	2	0	2	2
Indeterminant	0	0	0	0	0
Positive	2	0	2	0	0
Beckman Coulter ICON 20	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	24	0	24	24
Indeterminant	0	0	0	0	0
Positive	24	0	24	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Cardinal Health hCG Combo	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	30	0	30	30
Indeterminant	0	0	0	0	0
Positive	30	0	30	0	0
Wampole Clearview hCG Combo	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	1	0	1	1
Indeterminant	0	0	0	0	0
Positive	1	0	1	0	0
Cardinal Health hCG Cassette	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	4	0	4	4
Indeterminant	0	0	0	0	0
Positive	4	0	4	0	0
Clarity hCG Cassette	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	1	0	1	1
Indeterminant	0	0	0	0	0
Positive	1	0	1	0	0
Cen Med Elite	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	12	0	12	12
Indeterminant	1	0	0	0	0
Positive	11	0	12	0	0
Sekisui OSOM Combo (Genzyme)	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	12	0	12	12
Indeterminant	0	0	0	0	0
Positive	12	0	12	0	0
Henry Schein One Step	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	4	0	4	4
Indeterminant	0	0	0	0	0
Positive	4	0	4	0	0
Immunostics hCG Detector	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	1	0	1	1
Indeterminant	0	0	0	0	0
Positive	1	0	1	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Inverness Medical BioStar Acceava	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	1	0	1	1
Indeterminant	0	0	0	0	0
Positive	1	0	1	0	0
Jant Accutest Value+	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	2	0	2	2
Indeterminant	0	0	0	0	0
Positive	2	0	2	0	0
McKesson hCG Cassette	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	1	0	1	1
Indeterminant	0	0	0	0	0
Positive	1	0	1	0	0
Meridian Diagnostics	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	1	0	1	1
Indeterminant	0	0	0	0	0
Positive	1	0	1	0	0
Polymedco Poly stat	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	5	0	5	5
Indeterminant	0	0	0	0	0
Positive	5	0	5	0	0
Quidel QuickVue	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	63	0	63	62
Indeterminant	0	0	0	0	0
Positive	63	0	63	0	1
R. Medical One Step Strip	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	3	0	3	3
Indeterminant	0	0	0	0	0
Positive	3	0	3	0	0
Stanbio QuPID	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	15	0	15	15
Indeterminant	0	0	0	0	0
Positive	15	0	15	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

Fisher Sure-Vue serum/urine	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	38	0	38	37
Indeterminant	0	0	0	0	0
Positive	38	0	38	0	1

Fisher Sure-Vue STAT serum/urine	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	24	0	24	24
Indeterminant	0	0	0	0	0
Positive	24	0	24	0	0

Signify hCG Urine/Serum	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	3	0	3	3
Indeterminant	0	0	0	0	0
Positive	3	0	3	0	0

Siemens Immulite	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	4	0	4	4
Indeterminant	0	0	0	0	0
Positive	4	0	4	0	0

Fisher Sure-Vue (urine only)	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	8	0	8	8
Indeterminant	0	0	0	0	0
Positive	8	0	8	0	0

True Medix	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	4	0	4	4
Indeterminant	0	0	0	0	0
Positive	4	0	4	0	0

All Methods & Instruments	Specimen:U36	Specimen:U37	Specimen:U38	Specimen:U39	Specimen:U40
Negative	0	314	0	314	312
Indeterminant	1	0	0	0	0
Positive	313	0	314	0	2

New York State Department of Health – Wadsworth Center
Urinalysis/ Urine Pregnancy Proficiency Testing – 26 October 2015

hCG Summary:

[Code], hCG Kit name	<u>U36</u>			<u>U37</u>			<u>U38</u>			<u>U39</u>			<u>U40</u>		
	N	I	P	N	I	P	N	I	P	N	I	P	N	I	P
[ZZZ] Other	0	0	5	5	0	0	0	0	5	5	0	0	5	0	0
[AC] Alere hCG Cassette	0	0	9	9	0	0	0	0	9	9	0	0	9	0	0
[BA] Siemens Clinitek Status	0	0	18	18	0	0	0	0	18	18	0	0	18	0	0
[BC] Beckman Coulter ICON II	0	0	2	2	0	0	0	0	2	2	0	0	2	0	0
[BN] Beckman Coulter ICON 25	0	0	17	17	0	0	0	0	17	17	0	0	17	0	0
[BR] BTNX Rapid Response	0	0	2	2	0	0	0	0	2	2	0	0	2	0	0
[BT] Beckman Coulter ICON 20	0	0	24	24	0	0	0	0	24	24	0	0	24	0	0
[CA] Cardinal Health hCG Combo	0	0	30	30	0	0	0	0	30	30	0	0	30	0	0
[CC] Wampole Clearview hCG Combo	0	0	1	1	0	0	0	0	1	1	0	0	1	0	0
[CH] Cardinal Health hCG Cassette	0	0	4	4	0	0	0	0	4	4	0	0	4	0	0
[CL] Clarity hCG Cassette	0	0	1	1	0	0	0	0	1	1	0	0	1	0	0
[CM] Cen Med Elite	0	1	11	12	0	0	0	0	12	12	0	0	12	0	0
[GO] Sekisui OSOM Combo (Genzyme)	0	0	12	12	0	0	0	0	12	12	0	0	12	0	0
[HS] Henry Schein One Step	0	0	4	4	0	0	0	0	4	4	0	0	4	0	0
[ID] Immunostics hCG Detector	0	0	1	1	0	0	0	0	1	1	0	0	1	0	0
[IM] Inverness Medical BioStar Acceava	0	0	1	1	0	0	0	0	1	1	0	0	1	0	0
[JA] Jant Accutest Value±	0	0	2	2	0	0	0	0	2	2	0	0	2	0	0
[MC] McKesson hCG Cassette	0	0	1	1	0	0	0	0	1	1	0	0	1	0	0
[MD] Meridian Diagnostics	0	0	1	1	0	0	0	0	1	1	0	0	1	0	0
[PM] Polymedco Poly stat	0	0	5	5	0	0	0	0	5	5	0	0	5	0	0
[QQ] Quidel QuickVue	0	0	63	63	0	0	0	0	63	63	0	0	62	0	1
[RM] R. Medical One Step Strip	0	0	3	3	0	0	0	0	3	3	0	0	3	0	0
[SB] Stanbio QuPID	0	0	15	15	0	0	0	0	15	15	0	0	15	0	0
[SE] Fisher Sure-Vue serum/urine	0	0	38	38	0	0	0	0	38	38	0	0	37	0	1
[SF] Fisher Sure-Vue STAT serum/urine	0	0	24	24	0	0	0	0	24	24	0	0	24	0	0
[SG] Signify hCG Urine/Serum	0	0	3	3	0	0	0	0	3	3	0	0	3	0	0
[SI] Siemens Immulite	0	0	4	4	0	0	0	0	4	4	0	0	4	0	0
[SU] Fisher Sure-Vue (urine only)	0	0	8	8	0	0	0	0	8	8	0	0	8	0	0
[TX] True Medix	0	0	4	4	0	0	0	0	4	4	0	0	4	0	0
[99] All hCG Methods	0	1	313	314	0	0	0	0	314	314	0	0	312	0	2