

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Amphetamines

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	1000	5	5	15	0	0	0
American BioMedica Corp	1000	1	1	3	0	0	0
American Diag Prod Multi-5	1000	1	1	3	0	0	0
Bayer ADVIA	1000	4	4	12	0	0	0
Beckman EMIT	1000	1	1	3	0	0	0
Beckman	300	1	1	3	0	0	0
Beckman	1000	24	24	72	0	0	0
Roche /Hitachi enzymatic	1000	1	1	3	0	0	0
Branan QuickTox Dipcard	1000	1	0	3	0	1	0
Carolina Liquid Chemistries	1000	4	4	12	0	0	0
Microgenics CEDIA	1000	6	6	18	0	0	0
Dade Dimension	500	4	4	12	0	0	0
Dade Dimension	1000	32	32	96	0	0	0
Diagnostic Reagents (DRI)	1000	3	3	9	0	0	0
Behring SYVA EMIT	100	1	1	3	0	0	0
Behring SYVA EMIT	300	4	4	12	0	0	0
Behring SYVA EMIT	500	3	3	9	0	0	0
Behring SYVA EMIT	1000	22	22	65	1	0	0
E-Screen	1000	1	1	3	0	0	0
Abbott FPIA	500	1	1	3	0	0	0
Abbott FPIA	1000	3	3	9	0	0	0
I Screen One Step Drug Scr	1000	3	1	9	0	2	0
Roche KIMS	500	3	3	9	0	0	0
Roche KIMS	1000	16	16	48	0	0	0
Microgenics DRI	800	1	1	3	0	0	0
Microgenics DRI	1000	4	4	12	0	0	0
MedTox Sure-Screen	300	1	1	3	0	0	0
Med Tox Profile III	1000	2	1	6	0	1	0
Med Tox Profile V	500	3	3	9	0	0	0
Med Tox Profile V	1000	2	2	6	0	0	0
Med Tox Profile II	1000	4	4	12	0	0	0
RocheElectrochemiluminesc	1000	1	1	3	0	0	0
Siemens Dimension vista-EMIT	1000	1	1	3	0	0	0
Biosite Triage/ Plus TCA	1000	26	22	78	0	4	0
BiositeTriage Meter Plus/Pro	999	1	0	3	0	1	0
BiositeTriage Meter Plus/Pro	1000	29	15	87	0	14	0
BioRad Tox/See	1000	2	1	6	0	1	0
Verdict II/MedTox	1000	1	1	3	0	0	0
Vitros/Johnson & Johnson	1000	12	12	36	0	0	0
Varian OnTrak Test Cup	1000	1	0	3	0	1	0
Participant Totals:		236	211	707	1	25	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Barbiturates

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	200	5	0	20	0	0	0
American BioMedica Corp	300	1	0	4	0	0	0
Bayer ADVIA	200	5	0	20	0	0	0
Beckman EMIT	200	1	0	4	0	0	0
Beckman	200	25	0	100	0	0	0
Beckman	300	2	0	8	0	0	0
Roche /Hitachi enzymatic	200	1	0	4	0	0	0
Branan QuickTox Dipcard	300	1	0	4	0	0	0
Carolina Liquid Chemistries	200	4	0	16	0	0	0
Microgenics CEDIA	200	1	0	4	0	0	0
Microgenics CEDIA	300	2	0	8	0	0	0
Dade Dimension	200	35	0	140	0	0	0
Dade Dimension	1000	1	0	4	0	0	0
Diagnostic Reagents (DRI)	200	2	0	8	0	0	0
Diagnostic Reagents (DRI)	300	1	0	4	0	0	0
Behring SYVA EMIT	200	19	0	76	0	0	0
Behring SYVA EMIT	300	11	0	44	0	0	0
Abbott FPIA	200	2	0	8	0	0	0
Abbott FPIA	300	2	0	8	0	0	0
I Screen One Step Drug Scr	300	3	0	12	0	0	0
Roche KIMS	200	21	0	84	0	0	0
Roche KIMS	300	1	0	4	0	0	0
Microgenics DRI	200	3	0	12	0	0	0
Microgenics DRI	300	1	0	4	0	0	0
Med Tox Profile III	200	2	0	8	0	0	0
Med Tox Profile V	200	5	0	20	0	0	0
Med Tox Profile II	200	4	0	16	0	0	0
RocheElectrochemiluminesc	200	1	0	4	0	0	0
Siemens Dimension vista-EMIT	200	1	0	4	0	0	0
Biosite Triage/ Plus TCA	300	26	0	104	0	0	0
BiositeTriage Meter Plus/Pro	299	1	0	4	0	0	0
BiositeTriage Meter Plus/Pro	300	29	0	116	0	0	0
BioRad Tox/See	300	2	0	8	0	0	0
Vitros/Johnson & Johnson	200	8	0	32	0	0	0
Vitros/Johnson & Johnson	300	3	0	12	0	0	0
Varian OnTrak Test Cup	200	1	0	4	0	0	0
Participant Totals:		233	0	932	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Benzodiazepines

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	200	5	5	15	0	0	0
American BioMedica Corp	300	1	1	3	0	0	0
Bayer ADVIA	200	5	5	15	0	0	0
Beckman EMIT	200	1	1	3	0	0	0
Beckman	200	24	24	72	0	0	0
Beckman	300	3	3	9	0	0	0
Roche /Hitachi enzymatic	300	1	1	3	0	0	0
Branan QuickTox Dipcard	300	1	1	3	0	0	0
Carolina Liquid Chemistries	200	3	3	9	0	0	0
Carolina Liquid Chemistries	300	1	1	3	0	0	0
Microgenics CEDIA	200	1	1	3	0	0	0
Microgenics CEDIA	300	2	2	6	0	0	0
Dade Dimension	200	35	35	105	0	0	0
Dade Dimension	1000	1	1	3	0	0	0
Diagnostic Reagents (DRI)	200	2	2	6	0	0	0
Diagnostic Reagents (DRI)	300	1	1	3	0	0	0
Behring SYVA EMIT	100	2	2	6	0	0	0
Behring SYVA EMIT	200	16	16	48	0	0	0
Behring SYVA EMIT	300	12	12	36	0	0	0
Abbott FPIA	200	2	2	6	0	0	0
Abbott FPIA	300	2	2	6	0	0	0
I Screen One Step Drug Scr	300	3	3	9	0	0	0
Roche KIMS	100	4	4	12	0	0	0
Roche KIMS	200	9	9	27	0	0	0
Roche KIMS	300	9	9	27	0	0	0
Microgenics DRI	40	1	1	3	0	0	0
Microgenics DRI	200	2	2	6	0	0	0
Microgenics DRI	300	1	1	3	0	0	0
MedTox Sure-Screen	200	1	1	3	0	0	0
Med Tox Profile III	300	2	0	6	0	2	0
Med Tox Profile V	150	3	3	9	0	0	0
Med Tox Profile V	300	2	2	6	0	0	0
Med Tox Profile II	300	4	4	12	0	0	0
RocheElectrochemiluminesc	100	1	1	3	0	0	0
Siemens Dimension vista-EMIT	200	1	1	3	0	0	0
Biosite Triage/ Plus TCA	300	26	26	78	0	0	0
BiositeTriage Meter Plus/Pro	299	1	1	3	0	0	0
BiositeTriage Meter Plus/Pro	300	29	29	87	0	0	0
BioRad Tox/See	300	2	2	6	0	0	0
Verdict II/MedTox	300	1	1	3	0	0	0
Vitros/Johnson & Johnson	200	8	8	24	0	0	0
Vitros/Johnson & Johnson	300	4	4	12	0	0	0
Varian OnTrak Test Cup	100	1	1	3	0	0	0
Participant Totals:		236	234	708	0	2	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Cannabinoids

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	50	4	0	16	0	0	0
American BioMedica Corp	50	1	1	3	0	0	0
American Diag Prod Multi-5	50	1	1	3	0	0	0
Bayer ADVIA	50	5	1	19	0	0	0
Beckman EMIT	50	1	0	4	0	0	0
Beckman	20	2	2	6	0	0	0
Beckman	25	1	1	3	0	0	0
Beckman	50	22	9	79	0	0	0
Beckman	100	2	0	8	0	0	0
Roche /Hitachi enzymatic	50	1	1	3	0	0	0
Branan QuickTox Dipcard	50	1	1	3	0	0	0
Carolina Liquid Chemistries	50	3	2	10	0	0	0
Carolina Liquid Chemistries	439	1	0	4	0	0	0
Microgenics CEDIA	25	1	1	3	0	0	0
Microgenics CEDIA	50	3	2	10	0	0	0
Dade Dimension	25	3	3	9	0	0	0
Dade Dimension	50	31	5	119	0	0	0
Dade Dimension	1000	1	0	4	0	0	0
Diagnostic Reagents (DRI)	50	3	1	11	0	0	0
Behring SYVA EMIT	20	2	2	6	0	0	0
Behring SYVA EMIT	25	1	1	3	0	0	0
Behring SYVA EMIT	50	25	9	91	0	0	0
E-Screen	50	1	0	4	0	0	0
Abbott FPIA	25	1	1	3	0	0	0
Abbott FPIA	50	3	3	9	0	0	0
I Screen One Step Drug Scr	50	3	2	10	0	0	0
Roche KIMS	20	3	3	9	0	0	0
Roche KIMS	50	18	16	56	0	0	0
Roche KIMS	100	2	0	8	0	0	0
Microgenics DRI	20	1	1	3	0	0	0
Microgenics DRI	30	1	1	3	0	0	0
Microgenics DRI	50	3	2	10	0	0	0
MedTox Sure-Screen	100	1	0	3	0	0	0
Med Tox Profile III	50	2	1	7	0	0	0
Med Tox Profile V	50	5	5	15	0	0	0
Med Tox Profile II	50	4	2	14	0	0	0
RocheElectrochemiluminesc	20	1	1	3	0	0	0
Siemens Dimension vista-EMIT	50	1	1	3	0	0	0
Biosite Triage/ Plus TCA	50	26	26	77	1	0	0
BiositeTriage Meter Plus/Pro	50	30	24	96	0	0	0
BioRad Tox/See	50	2	2	6	0	0	0
Verdict II/MedTox	50	1	1	3	0	0	0
Vitros/Johnson & Johnson	50	12	0	48	0	0	0
Varian OnTrak Test Cup	50	1	0	4	0	0	0
Participant Totals:		237	135	811	1	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Cocaine/metab

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	300	5	5	15	0	0	0
American BioMedica Corp	300	1	1	3	0	0	0
American Diag Prod Multi-5	300	1	1	3	0	0	0
Bayer ADVIA	300	5	5	15	0	0	0
Beckman EMIT	300	1	1	3	0	0	0
Beckman	150	2	2	6	0	0	0
Beckman	300	25	25	75	0	0	0
Roche /Hitachi enzymatic	300	1	1	3	0	0	0
Branan QuickTox Dipcard	300	1	1	3	0	0	0
Carolina Liquid Chemistries	300	2	2	6	0	0	0
Carolina Liquid Chemistries	408	1	1	3	0	0	0
Microgenics CEDIA	100	1	1	3	0	0	0
Microgenics CEDIA	300	4	4	12	0	0	0
Dade Dimension	150	4	4	12	0	0	0
Dade Dimension	300	31	31	93	0	0	0
Dade Dimension	1000	1	0	3	0	0	0
Diagnostic Reagents (DRI)	300	3	3	9	0	0	0
Behring SYVAEMITRapid Scr	300	1	1	3	0	0	0
Behring SYVA EMIT	150	2	2	6	0	0	0
Behring SYVA EMIT	200	1	1	3	0	0	0
Behring SYVA EMIT	300	23	23	69	0	0	0
E-Screen	300	1	1	3	0	0	0
Abbott FPIA	300	4	4	12	0	0	0
I Screen One Step Drug Scr	300	3	3	9	0	0	0
Roche KIMS	300	22	22	66	0	0	0
Microgenics DRI	200	1	1	3	0	0	0
Microgenics DRI	300	5	5	15	0	0	0
MedTox Sure-Screen	100	1	1	3	0	0	0
LC/MS/MS	200	1	1	3	0	0	0
Med Tox Profile III	300	2	2	6	0	0	0
Med Tox Profile V	150	3	3	9	0	0	0
Med Tox Profile V	300	2	2	6	0	0	0
Med Tox Profile II	300	4	4	12	0	0	0
RocheElectrochemiluminesc	300	1	1	3	0	0	0
Siemens Dimension vista-EMIT	300	1	1	3	0	0	0
Biosite Triage/ Plus TCA	300	25	25	75	0	0	0
BiositeTriage Meter Plus/Pro	299	1	1	3	0	0	0
BiositeTriage Meter Plus/Pro	300	30	30	90	0	0	0
BioRad Tox/See	300	2	2	6	0	0	0
Verdict II/MedTox	300	1	1	3	0	0	0
Vitros/Johnson & Johnson	200	1	1	3	0	0	0
Vitros/Johnson & Johnson	300	12	12	36	0	0	0
Varian OnTrak Test Cup	300	1	1	3	0	0	0
Participant Totals:		240	239	720	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Ethanol

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	40	1	2	2	0	0	0
Alcohol Dehydrogenase	10	5	10	10	0	0	0
Alcohol Dehydrogenase	19	1	2	2	0	0	0
Alcohol Dehydrogenase	20	1	2	2	0	0	0
Alcohol Dehydrogenase	40	1	2	2	0	0	0
Abbott REA(Radiative Energy Att)	13	1	2	2	0	0	0
Bayer ADVIA	10	1	2	2	0	0	0
Bayer ADVIA	20	3	6	6	0	0	0
Beckman EMIT	10	1	2	2	0	0	0
Beckman	5	2	4	4	0	0	0
Beckman	10	5	10	10	0	0	0
Beckman	11	1	2	2	0	0	0
Beckman	20	1	2	2	0	0	0
Roche /Hitachi enzymatic	9	1	2	2	0	0	0
Roche /Hitachi enzymatic	10	3	6	6	0	0	0
Roche /Hitachi enzymatic	20	1	2	2	0	0	0
Carolina Liquid Chemistries	10	1	2	2	0	0	0
Carolina Liquid Chemistries	20	1	2	2	0	0	0
Carolina Liquid Chemistries	50	1	1	3	0	0	0
Microgenics CEDIA	20	2	4	4	0	0	0
Dade Dimension	3	1	2	2	0	0	0
Dade Dimension	50	6	7	17	0	0	0
Diagnostic Chemicals Ltd	10	1	2	2	0	0	0
Diagnostic Reagents (DRI)	20	1	2	2	0	0	0
Diagnostic Reagents (DRI)	50	1	1	3	0	0	0
Behring SYVA EMIT	10	4	8	8	0	0	0
Behring SYVA EMIT	20	6	12	12	0	0	0
Behring SYVA EMIT	50	1	1	3	0	0	0
GC Flame Ionization	10	1	2	2	0	0	0
Roche ONLINE Ethanol	50	1	1	3	0	0	0
Roche KIMS	13	1	2	2	0	0	0
Roche KIMS	20	1	2	2	0	0	0
Microgenics DRI	20	3	6	6	0	0	0
Microgenics DRI	30	1	2	2	0	0	0
Microgenics DRI	50	2	3	5	0	0	0
Roche Electrochemiluminescence	25	1	2	2	0	0	0
Siemens Dimension Alcohol	5	1	2	2	0	0	0
Vitros/Johnson & Johnson	10	1	2	2	0	0	0
Vitros/Johnson & Johnson	25	1	2	2	0	0	0
Participant Totals:		69	128	148	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Methadone

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	250	1	1	3	0	0	0
Abbott Aeroset/Architect	300	4	4	12	0	0	0
American BioMedica Corp	300	1	1	3	0	0	0
Bayer ADVIA	300	3	3	9	0	0	0
Beckman EMIT	300	1	1	3	0	0	0
Beckman	250	1	1	3	0	0	0
Beckman	300	17	17	51	0	0	0
Roche /Hitachi enzymatic	300	1	1	3	0	0	0
Branan QuickTox Dipcard	300	1	1	3	0	0	0
Carolina Liquid Chemistries	300	4	4	12	0	0	0
Microgenics CEDIA	250	1	1	3	0	0	0
Microgenics CEDIA	300	2	2	6	0	0	0
Dade Dimension	300	21	21	63	0	0	0
Dade Dimension	1000	1	1	3	0	0	0
Diagnostic Reagents (DRI)	300	2	2	6	0	0	0
Behring SYVA EMIT	50	1	1	3	0	0	0
Behring SYVA EMIT	150	2	2	6	0	0	0
Behring SYVA EMIT	200	1	1	3	0	0	0
Behring SYVA EMIT	250	1	1	3	0	0	0
Behring SYVA EMIT	300	23	23	69	0	0	0
Behring SYVA EMIT	500	1	1	3	0	0	0
Abbott FPIA	250	1	1	3	0	0	0
Abbott FPIA	300	2	2	6	0	0	0
GC Mass Spectrometry	300	1	1	3	0	0	0
I Screen One Step Drug Scr	300	3	3	9	0	0	0
Roche KIMS	300	19	19	57	0	0	0
Microgenics DRI	130	1	1	3	0	0	0
Microgenics DRI	300	3	3	9	0	0	0
MedTox Sure-Screen	200	1	1	3	0	0	0
Med Tox Profile III	300	2	2	6	0	0	0
Med Tox Profile V	200	3	3	9	0	0	0
Med Tox Profile V	300	2	2	6	0	0	0
Med Tox Profile II	300	4	4	12	0	0	0
RocheElectrochemiluminesc	300	1	1	3	0	0	0
Siemens Dimension vista-EMIT	300	1	1	3	0	0	0
Biosite Triage/ Plus TCA	300	5	5	15	0	0	0
BiositeTriage Meter Plus/Pro	299	1	1	3	0	0	0
BiositeTriage Meter Plus/Pro	300	13	13	39	0	0	0
BioRad Tox/See	300	2	2	6	0	0	0
Vitros/Johnson & Johnson	300	7	7	21	0	0	0
Participant Totals:		162	162	486	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Opiates

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	300	4	4	12	0	0	0
Abbott Aeroset/Architect	2000	1	1	3	0	0	0
American BioMedica Corp	300	1	1	3	0	0	0
American Diag Prod Multi-5	2000	1	1	3	0	0	0
Bayer ADVIA	300	5	5	15	0	0	0
Beckman EMIT	300	1	1	3	0	0	0
Beckman	200	1	1	3	0	0	0
Beckman	300	24	24	72	0	0	0
Beckman	2000	2	2	6	0	0	0
Roche /Hitachi enzymatic	300	1	1	3	0	0	0
Branan QuickTox Dipcard	300	1	1	3	0	0	0
Carolina Liquid Chemistries	300	3	3	9	0	0	0
Carolina Liquid Chemistries	399	1	1	3	0	0	0
Microgenics CEDIA	150	1	1	3	0	0	0
Microgenics CEDIA	300	2	2	6	0	0	0
Dade Dimension	300	23	23	69	0	0	0
Dade Dimension	1000	1	1	3	0	0	0
Dade Dimension	2000	12	12	36	0	0	0
Diagnostic Reagents (DRI)	300	3	3	9	0	0	0
Behring SYVA EMIT	200	2	2	6	0	0	0
Behring SYVA EMIT	300	26	26	78	0	0	0
Behring SYVA EMIT	2000	3	3	9	0	0	0
E-Screen	2000	1	1	3	0	0	0
Abbott FPIA	300	3	3	9	0	0	0
Abbott FPIA	400	1	1	3	0	0	0
I Screen One Step Drug Scr	300	3	3	9	0	0	0
Roche KIMS	300	21	21	63	0	0	0
Roche KIMS	2000	1	1	3	0	0	0
Microgenics DRI	50	1	1	3	0	0	0
Microgenics DRI	300	3	3	9	0	0	0
MedTox Sure-Screen	100	1	1	3	0	0	0
Med Tox Profile III	300	2	2	6	0	0	0
Med Tox Profile V	100	3	3	9	0	0	0
Med Tox Profile V	300	2	2	6	0	0	0
Med Tox Profile II	300	4	4	12	0	0	0
RocheElectrochemiluminesc	300	1	1	3	0	0	0
Siemens Dimension vista-EMIT	300	1	1	3	0	0	0
Biosite Triage/ Plus TCA	300	26	26	78	0	0	0
BiositeTriage Meter Plus/Pro	299	1	1	3	0	0	0
BiositeTriage Meter Plus/Pro	300	29	29	87	0	0	0
BioRad Tox/See	300	2	2	6	0	0	0
Verdict II/MedTox	300	1	1	3	0	0	0
Vitros/Johnson & Johnson	300	12	12	36	0	0	0
Vitros/Johnson & Johnson	2000	1	1	3	0	0	0
Varian OnTrak Test Cup	300	1	1	3	0	0	0
Participant Totals:		240	240	720	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

PCP

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
Abbott Aeroset/Architect	20	1	0	4	0	0	0
Abbott Aeroset/Architect	25	4	0	16	0	0	0
American BioMedica Corp	25	1	0	4	0	0	0
American Diag Prod Multi-5	25	1	0	4	0	0	0
Bayer ADVIA	25	4	0	16	0	0	0
Beckman EMIT	25	1	0	4	0	0	0
Beckman	25	21	0	84	0	0	0
Roche /Hitachi enzymatic	25	1	0	4	0	0	0
Branan QuickTox Dipcard	25	1	0	4	0	0	0
Carolina Liquid Chemistries	25	3	0	12	0	0	0
Carolina Liquid Chemistries	450	1	0	4	0	0	0
Microgenics CEDIA	25	4	0	16	0	0	0
Dade Dimension	25	32	0	128	0	0	0
Dade Dimension	1000	1	0	4	0	0	0
Diagnostic Reagents (DRI)	25	3	0	12	0	0	0
Behring SYVA EMIT	25	24	0	96	0	0	0
Behring SYVA EMIT	50	1	0	4	0	0	0
E-Screen	25	1	0	4	0	0	0
Abbott FPIA	25	3	0	12	0	0	0
I Screen One Step Drug Scr	25	3	0	12	0	0	0
Roche KIMS	25	23	0	92	0	0	0
Microgenics DRI	25	5	0	20	0	0	0
MedTox Sure-Screen	25	1	0	4	0	0	0
Med Tox Profile III	25	2	0	8	0	0	0
Med Tox Profile V	25	5	0	20	0	0	0
Med Tox Profile II	25	4	0	16	0	0	0
RocheElectrochemiluminesc	25	1	0	4	0	0	0
Siemens Dimension vista-EMIT	25	1	0	4	0	0	0
Biosite Triage/ Plus TCA	25	26	0	104	0	0	0
BiositeTriage Meter Plus/Pro	25	30	0	120	0	0	0
BioRad Tox/See	25	1	0	4	0	0	0
Vitros/Johnson & Johnson	25	11	0	44	0	0	0
Varian OnTrak Test Cup	25	1	0	4	0	0	0
Participant Totals:		222	0	888	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Propoxyphene

<i>Method Code</i>	<i>Cutoff</i>	<i>No. Labs</i>	<i>True Positive</i>	<i>True Negative</i>	<i>False Positive</i>	<i>False Negative</i>	<i>Cross-Reactant</i>
Abbott Aeroset/Architect	300	3	0	12	0	0	0
American BioMedica Corp	300	1	0	4	0	0	0
Beckman EMIT	300	1	0	4	0	0	0
Beckman	300	7	0	28	0	0	0
Roche /Hitachi enzymatic	300	1	0	4	0	0	0
Carolina Liquid Chemistries	224	1	0	4	0	0	0
Carolina Liquid Chemistries	300	2	0	8	0	0	0
Microgenics CEDIA	150	1	0	4	0	0	0
Microgenics CEDIA	300	1	0	4	0	0	0
Dade Dimension	300	2	0	8	0	0	0
Diagnostic Reagents (DRI)	300	2	0	8	0	0	0
Behring SYVA EMIT	200	1	0	4	0	0	0
Behring SYVA EMIT	300	16	0	64	0	0	0
Roche KIMS	300	9	0	36	0	0	0
Microgenics DRI	180	1	0	4	0	0	0
Microgenics DRI	300	4	0	16	0	0	0
Med Tox Profile III	300	2	0	8	0	0	0
Med Tox Profile V	300	3	0	12	0	0	0
Med Tox Profile II	300	3	0	12	0	0	0
RocheElectrochemiluminesc	300	1	0	4	0	0	0
ToxiLab	500	1	0	4	0	0	0
Biosite Triage/ Plus TCA	300	1	0	4	0	0	0
Method not listed	300	1	0	4	0	0	0
<i>Participant Totals:</i>		65	0	260	0	0	0

Please refer to the Performance Evaluation Report for the test specimen drug content.

Analytical Toxicology Proficiency Testing Program

Participant Performance Summary - Initial (Screen) Test

Clinical Toxicology - March 2010 Test Event

Tricyclic Antidepressants

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative	Cross-Reactant
American BioMedica Corp	1000	1	1	3	0	0	0
Dade Dimension	1000	1	1	3	0	0	0
Diagnostic Reagents (DRI)	150	1	1	3	0	0	0
Diagnostic Reagents (DRI)	300	1	1	3	0	0	0
Behring SYVAEMITRapid Scr	110	1	1	3	0	0	0
Behring SYVAEMITRapid Scr	1000	4	4	12	0	0	0
Behring SYVA EMIT	100	1	1	3	0	0	0
Behring SYVA EMIT	200	1	1	3	0	0	0
Behring SYVA EMIT	300	1	1	3	0	0	0
Abbott FPIA	500	1	1	3	0	0	0
I Screen One Step Drug Scr	1000	1	1	3	0	0	0
Microgenics DRI	300	3	3	9	0	0	0
Microgenics DRI	500	3	3	9	0	0	0
Med Tox Profile III	300	2	2	6	0	0	0
Med Tox Profile V	300	4	4	12	0	0	0
Med Tox Profile II	300	5	5	15	0	0	0
ToxiLab	500	2	2	6	0	0	0
Biosite Triage/ Plus TCA	1000	24	23	72	0	1	0
BiositeTriage Meter Plus/Pro	999	1	1	3	0	0	0
BiositeTriage Meter Plus/Pro	1000	26	23	78	0	3	0
BioRad Tox/See	1000	1	1	3	0	0	0
Vitros/Johnson & Johnson	300	1	1	3	0	0	0
Varian OnTrak Test Cup	1000	1	1	3	0	0	0
Participant Totals:		87	83	261	0	4	0