

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

6-Acetylmorphine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	10	1	1	0	0	0
GCMS	4	1	1	0	0	0
GCMS	5	2	2	0	0	0
GCMS	10	21	21	0	0	0
GCMS	15	1	1	0	0	0
GCMS	20	1	0	1	0	0
GCMS	300	1	0	1	0	0
LCMM	2	1	1	0	0	0
LCMM	4	1	1	0	0	0
LCMM	5	2	2	0	0	0
LCMM	10	11	11	0	0	0
LCMM	20	2	0	2	0	0
LCMM	25	1	0	1	0	0
LCMM	50	1	0	1	0	0
TOXI	200	2	0	2	0	0
TOXI	1000	1	0	1	0	0

Alprazolam

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	0	1	0	0
GCMS	100	3	0	3	0	0
GCMS	150	1	0	1	0	0
GCMS	200	2	0	2	0	0
LCDA	200	1	0	1	0	0
LCDA	400	1	0	1	0	0
LCMM	5	1	0	1	0	0
LCMM	15	2	0	2	0	0
LCMM	20	1	0	1	0	0
LCMM	25	2	0	2	0	0
LCMM	50	6	0	6	0	0
LCMM	75	2	0	2	0	0
LCMM	100	3	0	3	0	0
LCMM	200	2	0	2	0	0
LCMM	300	4	0	4	0	0

Amitriptyline

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	1	0	0	0
GCMS	100	1	1	0	0	0
GCMS	250	1	1	0	0	0
GCMS	500	1	1	0	0	0
GCNP	50	1	1	0	0	0
LCFW	25	1	1	0	0	0
LCMM	10	1	1	0	0	0
LCMM	25	1	1	0	0	0
LCMM	50	4	4	0	0	0
LCMM	100	5	5	0	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

LCMM	150	1	1	0	0	0
LCMM	200	2	2	0	0	0
TOXI	500	2	2	0	0	0

Amobarbital

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	1	0	0	0
GCMS		0	0	0	0	0
GCMS	20	2	2	0	0	0
GCMS	50	2	2	0	0	0
GCMS	100	6	6	0	0	0
GCMS	150	1	1	0	0	0
GCMS	200	17	17	0	0	0
GCMS	300	5	5	0	0	0
LCDA	200	1	1	0	0	0
LCMM	40	1	1	0	0	0
LCMM	100	2	2	0	0	0
TLCC	300	1	1	0	0	0
TLCC	1000	1	1	0	0	0

Amphetamine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
LCMM	50	1	1	0	0	0
LCMM	100	3	3	0	0	0

Barbiturates

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
TOXI	175	1	1	0	0	0
TOXI	500	1	1	0	0	0
TOXI	1000	1	1	0	0	0

Benzodiazepines

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
TLCC	1000	1	1	0	0	0
TOXI	100	2	2	0	0	0
TOXI	1000	1	1	0	0	0

Benzoyllecgonine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
LCMM	25	3	3	0	0	0
LCMM	50	1	1	0	0	0

Butalbital

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	1	0	0	0
GCMS		0	0	0	0	0
GCMS	20	1	1	0	0	0
GCMS	40	1	1	0	0	0
GCMS	50	4	4	0	0	0
GCMS	100	6	6	0	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

GCMS	150	1	1	0	0	0
GCMS	200	17	16	0	0	0
GCMS	300	6	6	0	0	0
LCDA	200	1	1	0	0	0
LCMM	20	1	1	0	0	0
LCMM	50	1	1	0	0	0
LCMM	100	6	6	0	0	0
LCMM	200	1	1	0	0	0
LCMM	300	1	1	0	0	0
TLCC	300	1	1	0	0	0
TLCC	1000	1	1	0	0	0

Codeine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	300	1	0	1	0	0
GCMS	50	1	0	1	0	0
GCMS	100	4	0	4	0	0
GCMS	150	1	0	1	0	0
GCMS	200	1	0	1	0	0
GCMS	300	13	0	13	0	0
GCMS	2000	4	0	4	0	0
LCMM	20	1	0	1	0	0
LCMM	25	2	0	2	0	0
LCMM	30	1	0	1	0	0
LCMM	50	9	0	9	0	0
LCMM	100	8	0	8	0	0
LCMM	300	3	0	3	0	0
THSC	2000	1	0	1	0	0
TLCC	300	1	0	1	0	0
TLCC	500	1	0	1	0	0
TOXI	200	2	0	2	0	0
TOXI	500	1	0	1	0	0

Doxepin

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	0	1	0	0
GCMS	100	1	0	1	0	0
GCMS	250	1	0	1	0	0
GCNP	50	1	0	1	0	0
LCFW	25	1	0	1	0	0
LCMM	50	3	0	3	0	0
LCMM	100	2	0	2	0	0
LCMM	150	1	0	1	0	0
LCMM	200	2	0	2	0	0
TOXI	500	1	0	1	0	0

EDDP

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	100	2	2	0	0	0
GCMS	150	2	2	0	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

GCMS	200	3	3	0	0	0
GCMS	250	1	1	0	0	0
GCMS	300	3	3	0	0	0
LCMM	5	1	1	0	0	0
LCMM	10	1	1	0	0	0
LCMM	25	1	1	0	0	0
LCMM	50	2	2	0	0	0
LCMM	75	1	1	0	0	0
LCMM	100	9	9	0	0	0
LCMM	200	2	2	0	0	0
LCMM	300	3	3	0	0	0
LDTM	30	1	1	0	0	0
THSD	150	1	1	0	0	0
TLCC	500	1	1	0	0	0
TOXI	1000	1	1	0	0	0

Ethanol

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
BECO	10	1	1	0	0	0
BECO	20	2	2	0	0	0
GCED	10	1	1	0	0	0
GCFL	10	14	13	1	0	0
GCFL	20	17	17	0	0	0
GCFL	40	3	3	0	0	0
GCFL	50	2	2	0	0	0
ROON	20	1	1	0	0	0
THSD	50	1	1	0	0	0

Methadone

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	1	0	0	0
GCMS	50	1	1	0	0	0
GCMS	100	5	5	0	0	0
GCMS	150	4	4	0	0	0
GCMS	200	4	4	0	0	0
GCMS	250	1	1	0	0	0
GCMS	300	13	13	0	0	0
LCMM	3	1	1	0	0	0
LCMM	10	1	1	0	0	0
LCMM	25	1	1	0	0	0
LCMM	50	3	3	0	0	0
LCMM	75	1	1	0	0	0
LCMM	100	9	9	0	0	0
LCMM	200	2	2	0	0	0
LCMM	300	3	3	0	0	0
LDTM	30	1	1	0	0	0
TLCC	300	1	1	0	0	0
TLCC	500	1	1	0	0	0
TOXI	250	1	1	0	0	0
TOXI	500	2	2	0	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

Methamphetamine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
LCMM	50	1	1	0	0	0
LCMM	100	3	3	0	0	0

Morphine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	300	1	1	0	0	0
GCMS	50	2	2	0	0	0
GCMS	100	4	4	0	0	0
GCMS	150	1	1	0	0	0
GCMS	200	1	1	0	0	0
GCMS	300	15	15	0	0	0
GCMS	2000	6	6	0	0	0
LCMM	25	2	2	0	0	0
LCMM	30	1	1	0	0	0
LCMM	50	9	9	0	0	0
LCMM	100	8	8	0	0	0
LCMM	300	3	3	0	0	0
TLCC	300	1	1	0	0	0
TLCC	500	1	1	0	0	0
TOXI	200	2	2	0	0	0
TOXI	1000	1	1	0	0	0

Morphine Free

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
LCMM	20	1	1	0	0	0

Nordiazepam

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	0	1	0	0
GCMS		0	0	1	0	0
GCMS	50	1	0	1	0	0
GCMS	100	5	0	5	0	0
GCMS	150	1	0	1	0	0
GCMS	200	5	0	5	0	0
GCMS	300	5	0	5	0	0
LCDA	200	1	0	1	0	0
LCDA	400	1	0	1	0	0
LCMM	10	1	0	1	0	0
LCMM	15	1	0	1	0	0
LCMM	20	2	0	2	0	0
LCMM	25	1	0	1	0	0
LCMM	40	2	0	2	0	0
LCMM	50	8	0	8	0	0
LCMM	75	2	0	2	0	0
LCMM	100	5	0	5	0	0
LCMM	200	3	0	3	0	0
LCMM	300	4	0	4	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

Nordoxepin

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	0	1	0	0
GCMS	100	1	0	1	0	0
GCMS	250	1	0	1	0	0
GCNP	50	1	0	1	0	0
LCFW	25	1	0	1	0	0
LCMM	50	1	0	1	0	0
LCMM	100	1	0	1	0	0
LCMM	150	1	0	1	0	0
TOXI	500	1	0	1	0	0

Norpropoxyphene

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	1	0	0	0
GCMS	100	3	3	0	0	0
GCMS	150	2	2	0	0	0
GCMS	200	9	9	0	0	0
GCMS	300	10	10	0	0	0
LCMM	10	1	1	0	0	0
LCMM	50	2	2	0	0	0
LCMM	100	6	6	0	0	0
LCMM	200	1	1	0	0	0
LCMM	300	1	1	0	0	0
TOXI	2000	2	0	2	0	0

Nortriptyline

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	1	0	0	0
GCMS	100	1	1	0	0	0
GCMS	250	1	1	0	0	0
GCMS	500	1	1	0	0	0
GCNP	50	1	1	0	0	0
LCFW	25	1	1	0	0	0
LCMM	10	1	1	0	0	0
LCMM	25	1	1	0	0	0
LCMM	50	4	4	0	0	0
LCMM	100	5	5	0	0	0
LCMM	150	1	1	0	0	0
LCMM	200	2	2	0	0	0
TOXI	500	2	2	0	0	0

Oxazepam

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	1	0	0	0
GCMS		0	0	0	0	0
GCMS	50	1	1	0	0	0
GCMS	100	6	6	0	0	0
GCMS	150	1	1	0	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

GCMS	200	5	5	0	0	0
GCMS	300	6	6	0	0	0
LCDA	200	1	1	0	0	0
LCDA	400	1	1	0	0	0
LCMM	5	1	1	0	0	0
LCMM	15	1	1	0	0	0
LCMM	20	2	2	0	0	0
LCMM	25	1	1	0	0	0
LCMM	40	2	2	0	0	0
LCMM	50	8	8	0	0	0
LCMM	75	2	2	0	0	0
LCMM	100	5	5	0	0	0
LCMM	200	3	3	0	0	0
LCMM	300	4	4	0	0	0
TLCC	300	1	1	0	0	0

Pentobarbital

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	0	1	0	0
GCMS		0	0	1	0	0
GCMS	20	1	0	1	0	0
GCMS	50	4	0	4	0	0
GCMS	60	1	0	1	0	0
GCMS	100	6	0	6	0	0
GCMS	150	1	0	1	0	0
GCMS	200	17	0	16	1	0
GCMS	300	5	0	5	0	0
LCDA	200	1	0	1	0	0
LCMM	100	3	0	3	0	0
LCMM	300	1	0	0	1	0
TLCC	300	1	0	1	0	0

Phencyclidine

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	25	1	1	0	0	0
GCMS		0	0	0	0	0
GCMS	3	1	1	0	0	0
GCMS	5	1	1	0	0	0
GCMS	10	4	4	0	0	0
GCMS	13	1	1	0	0	0
GCMS	15	1	1	0	0	0
GCMS	17	1	1	0	0	0
GCMS	20	1	1	0	0	0
GCMS	25	29	29	0	0	0
LCMM	5	2	2	0	0	0
LCMM	10	3	3	0	0	0
LCMM	25	5	5	0	0	0
LCMM	100	1	0	1	0	0
LCMS	50	1	1	0	0	0
TLCC	100	1	0	1	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

TLCC	500	1	0	1	0	0
TOXI	500	2	0	2	0	0

Propoxyphene

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMS	50	1	1	0	0	0
GCMS	100	2	2	0	0	0
GCMS	150	3	3	0	0	0
GCMS	200	3	3	0	0	0
GCMS	300	5	5	0	0	0
LCMM	3	1	1	0	0	0
LCMM	10	1	1	0	0	0
LCMM	25	1	1	0	0	0
LCMM	50	1	1	0	0	0
LCMM	75	1	1	0	0	0
LCMM	100	5	5	0	0	0
LCMM	200	1	1	0	0	0
LCMM	300	1	1	0	0	0
TOXI	500	3	3	0	0	0

Secobarbital

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	1	0	0	0
GCMS		0	0	0	0	0
GCMS	20	2	2	0	0	0
GCMS	50	4	4	0	0	0
GCMS	100	6	6	0	0	0
GCMS	150	1	1	0	0	0
GCMS	200	17	17	0	0	0
GCMS	300	6	6	0	0	0
LCDA	200	1	1	0	0	0
LCMM	50	1	1	0	0	0
LCMM	60	1	1	0	0	0
LCMM	100	5	5	0	0	0
LCMM	200	1	1	0	0	0
LCMM	300	1	1	0	0	0
TLCC	300	1	1	0	0	0
TLCC	1000	1	1	0	0	0

Temazepam

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	200	1	1	0	0	0
GCMS		0	0	0	0	0
GCMS	50	1	1	0	0	0
GCMS	100	5	5	0	0	0
GCMS	150	1	1	0	0	0
GCMS	200	2	2	0	0	0
GCMS	300	4	4	0	0	0
LCDA	200	1	1	0	0	0
LCDA	400	1	1	0	0	0

Participant Performance Evaluation - Confirmation Test

March 2013 Test Event

LCMM	5	1	1	0	0	0
LCMM	15	1	1	0	0	0
LCMM	20	2	2	0	0	0
LCMM	25	1	1	0	0	0
LCMM	40	1	1	0	0	0
LCMM	50	9	9	0	0	0
LCMM	75	2	2	0	0	0
LCMM	100	5	5	0	0	0
LCMM	200	3	3	0	0	0
LCMM	300	4	4	0	0	0

THCA

Method Code	Cutoff	No. Labs	True Positive	True Negative	False Positive	False Negative
GCMM	15	1	1	0	0	0
GCMS	3	3	3	0	0	0
GCMS	5	5	5	0	0	0
GCMS	8	1	1	0	0	0
GCMS	10	4	4	0	0	0
GCMS	15	26	26	0	0	0
LCMM	5	4	4	0	0	0
LCMM	10	4	4	0	0	0
LCMM	15	4	4	0	0	0
LCMM	20	1	1	0	0	0
LCMS	50	1	1	0	0	0
TLCC	25	1	1	0	0	0
TLCC	50	1	1	0	0	0
TLCC	100	1	1	0	0	0
TOXI	15	1	1	0	0	0
TOXI	25	1	1	0	0	0