

Urinalysis/Urine Pregnancy Proficiency Testing
14 February 2000

Specific Gravity

Method Code	# Reported	Mean $\pm$ 1 SD	Method Code	# Reported	Mean $\pm$ 1 SD
Sample U61					
1. Bayer Test Strip Manual (Multistik)	71	1.023 $\pm$ 0.003	2. Bayer Automated (Clinitek Systems)	206	1.024 $\pm$ 0.022
3. BMC Manual (Chemstrip)	24	1.019 $\pm$ 0.003	4. BMC Automated (UA Readers)	46	1.020 $\pm$ 0.001
7. Urinometer	1	1.014 $\pm$ 0.000	8. Refractometer	69	1.010 $\pm$ 0.001
10. Atlas System	53	1.010 $\pm$ 0.002	11. IRIS Mass Gravity Meter	9	1.012 $\pm$ 0.001
28. Chemstrip Criterion	8	1.020 $\pm$ 0.000	99. All methods	498	1.019 $\pm$ 0.016
Sample U62					
1. Bayer Test Strip Manual (Multistik)	71	1.011 $\pm$ 0.003	2. Bayer Automated (Clinitek Systems)	206	1.007 $\pm$ 0.071
3. BMC Manual (Chemstrip)	24	1.011 $\pm$ 0.002	4. BMC Automated (UA Readers)	46	1.011 $\pm$ 0.002
7. Urinometer	1	1.010 $\pm$ 0.000	8. Refractometer	69	1.006 $\pm$ 0.001
10. Atlas System	53	1.006 $\pm$ 0.001	11. IRIS Mass Gravity Meter	9	1.008 $\pm$ 0.001
28. Chemstrip Criterion	8	1.011 $\pm$ 0.002	99. All methods	498	1.008 $\pm$ 0.045
Sample U63					
1. Bayer Test Strip Manual (Multistik)	71	1.015 $\pm$ 0.005	2. Bayer Automated (Clinitek Systems)	206	1.005 $\pm$ 0.071
3. BMC Manual (Chemstrip)	24	1.006 $\pm$ 0.003	4. BMC Automated (UA Readers)	46	1.007 $\pm$ 0.004
7. Urinometer	1	1.013 $\pm$ 0.001	8. Refractometer	69	1.013 $\pm$ 0.001
10. Atlas System	53	1.014 $\pm$ 0.002	11. IRIS Mass Gravity Meter	9	1.014 $\pm$ 0.001
28. Chemstrip Criterion	8	1.006 $\pm$ 0.003	99. All methods	498	1.009 $\pm$ 0.046
Sample U64					
1. Bayer Test Strip Manual (Multistik)	71	1.020 $\pm$ 0.003	2. Bayer Automated (Clinitek Systems)	206	1.017 $\pm$ 0.075
3. BMC Manual (Chemstrip)	24	1.015 $\pm$ 0.002	4. BMC Automated (UA Readers)	46	1.017 $\pm$ 0.002
7. Urinometer	1	1.011 $\pm$ 0.000	8. Refractometer	69	1.010 $\pm$ 0.002
10. Atlas System	53	1.010 $\pm$ 0.002	11. IRIS Mass Gravity Meter	9	1.011 $\pm$ 0.000
28. Chemstrip Criterion	8	1.016 $\pm$ 0.002	99. All methods	498	1.016 $\pm$ 0.048
Sample U65					
1. Bayer Test Strip Manual (Multistik)	71	1.023 $\pm$ 0.004	2. Bayer Automated (Clinitek Systems)	206	1.021 $\pm$ 0.075
3. BMC Manual (Chemstrip)	24	1.017 $\pm$ 0.003	4. BMC Automated (UA Readers)	46	1.018 $\pm$ 0.003
7. Urinometer	1	1.016 $\pm$ 0.000	8. Refractometer	69	1.015 $\pm$ 0.002
10. Atlas System	53	1.015 $\pm$ 0.002	11. IRIS Mass Gravity Meter	9	1.017 $\pm$ 0.001
28. Chemstrip Criterion	8	1.018 $\pm$ 0.003	99. All methods	498	1.019 $\pm$ 0.048

Urinalysis/Urine Pregnancy Proficiency Testing
14 February 2000

pH

Method Code	# Reported	Mean $\pm$ 1 SD	Method Code	# Reported	Mean $\pm$ 1 SD
Sample U61					
1. Bayer Test Strip Manual (Multistik)	104	5.7 $\pm$ 0.51	2. Bayer Automated (Clinitek Systems)	218	5.2 $\pm$ 0.30
3. BMC Manual (Chemstrip)	42	5.1 $\pm$ 0.24	4. BMC Automated (UA Readers)	51	5.5 $\pm$ 0.52
9. pH Meter	2	5.5 $\pm$ 0.00	10. Atlas Systems	53	5.9 $\pm$ 0.23
12. IRIS Clinitek Reader	8	5.1 $\pm$ 0.23	13. Iris BMC UA Reader	5	5.4 $\pm$ 0.42
17. Hitachi 747	2	5.3 $\pm$ 0.35	28. Chemstrip Criterion	10	5.9 $\pm$ 0.32
99. All methods	509	5.4 $\pm$ 0.45			
Sample U62					
1. Bayer Test Strip Manual (Multistik)	104	6.8 $\pm$ 0.38	2. Bayer Automated (Clinitek Systems)	218	6.8 $\pm$ 0.59
3. BMC Manual (Chemstrip)	42	6.9 $\pm$ 0.21	4. BMC Automated (UA Readers)	51	7.0 $\pm$ 0.10
9. pH Meter	2	6.9 $\pm$ 0.07	10. Atlas Systems	53	7.0 $\pm$ 0.10
12. IRIS Clinitek Reader	8	6.9 $\pm$ 0.18	13. Iris BMC UA Reader	5	6.9 $\pm$ 0.22
17. Hitachi 747	2	6.9 $\pm$ 0.14	28. Chemstrip Criterion	10	7.0 $\pm$ 0.00
99. All methods	509	6.9 $\pm$ 0.44			
Sample U63					
1. Bayer Test Strip Manual (Multistik)	104	5.4 $\pm$ 0.49	2. Bayer Automated (Clinitek Systems)	218	5.1 $\pm$ 0.43
3. BMC Manual (Chemstrip)	42	5.4 $\pm$ 0.57	4. BMC Automated (UA Readers)	51	6.3 $\pm$ 0.78
9. pH Meter	2	6.7 $\pm$ 0.07	10. Atlas Systems	53	5.5 $\pm$ 0.31
12. IRIS Clinitek Reader	8	5.9 $\pm$ 0.50	13. Iris BMC UA Reader	5	5.6 $\pm$ 0.65
17. Hitachi 747	2	5.3 $\pm$ 0.49	28. Chemstrip Criterion	10	5.9 $\pm$ 0.69
99. All methods	509	5.4 $\pm$ 0.61			
Sample U64					
1. Bayer Test Strip Manual (Multistik)	104	6.3 $\pm$ 0.29	2. Bayer Automated (Clinitek Systems)	218	6.2 $\pm$ 0.50
3. BMC Manual (Chemstrip)	42	6.0 $\pm$ 0.15	4. BMC Automated (UA Readers)	51	6.3 $\pm$ 0.25
9. pH Meter	2	6.3 $\pm$ 0.07	10. Atlas Systems	53	6.5 $\pm$ 0.10
12. IRIS Clinitek Reader	8	6.2 $\pm$ 0.26	13. Iris BMC UA Reader	5	6.2 $\pm$ 0.45
17. Hitachi 747	2	6.0 $\pm$ 0.07	28. Chemstrip Criterion	10	6.3 $\pm$ 0.24
99. All methods	509	6.2 $\pm$ 0.39			
Sample U65					
1. Bayer Test Strip Manual (Multistik)	104	6.3 $\pm$ 0.67	2. Bayer Automated (Clinitek Systems)	218	6.2 $\pm$ 0.49
3. BMC Manual (Chemstrip)	42	5.9 $\pm$ 0.86	4. BMC Automated (UA Readers)	51	6.2 $\pm$ 0.24
9. pH Meter	2	6.3 $\pm$ 0.00	10. Atlas Systems	53	6.5 $\pm$ 0.13
12. IRIS Clinitek Reader	8	6.0 $\pm$ 0.46	13. Iris BMC UA Reader	5	5.9 $\pm$ 0.65
17. Hitachi 747	2	6.1 $\pm$ 0.07	28. Chemstrip Criterion	10	6.1 $\pm$ 0.16
99. All methods	509	6.2 $\pm$ 0.54			

Urinalysis/Urine Pregnancy Proficiency Testing
14 February 2000

Bilirubin

Method Code	1	2	3	4	5	6	7	Method Code	1	2	3	4	5	6	7
Sample U61															
1. Bayer Test Strip Manual (Multistik)	3	1	9	31	37	0	18	2. Bayer Automated (Clinitek Systems)	0	0	126	75	16	0	0
3. BMC Manual (Chemstrip)	3	1	12	15	5	0	0	4. BMC Automated (UA Readers)	9	0	40	0	0	0	1
10. Atlas System	0	0	2	50	0	0	0	12. Iris Clinitek Reader	3	0	5	0	0	0	0
13. Iris BMC UA Reader	2	0	3	0	0	0	0	15. Ictotest	1	0	11	9	4	146	0
28. Chemstrip Criterion	9	0	1	0	0	0	0	99. All methods	31	3	213	181	63	149	19
Sample U62															
1. Bayer Test Strip Manual (Multistik)	97	0	0	0	0	0	0	2. Bayer Automated (Clinitek Systems)	216	0	0	0	0	0	0
3. BMC Manual (Chemstrip)	35	0	0	0	0	0	0	4. BMC Automated (UA Readers)	50	0	0	0	0	0	0
10. Atlas System	51	0	0	0	0	0	1	12. Iris Clinitek Reader	8	0	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	0	15. Ictotest	12	0	0	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	0	0	99. All methods	492	0	0	0	0	0	1
Sample U63															
1. Bayer Test Strip Manual (Multistik)	20	3	20	31	9	0	16	2. Bayer Automated (Clinitek Systems)	19	2	158	37	1	0	0
3. BMC Manual (Chemstrip)	12	5	12	6	1	0	0	4. BMC Automated (UA Readers)	37	0	12	0	0	0	1
10. Atlas System	8	2	17	25	0	0	0	12. Iris Clinitek Reader	6	0	2	0	0	0	0
13. Iris BMC UA Reader	4	1	0	0	0	0	0	15. Ictotest	1	1	13	7	0	118	0
28. Chemstrip Criterion	10	0	0	0	0	0	0	99. All methods	120	14	238	107	11	121	17
Sample U64															
1. Bayer Test Strip Manual (Multistik)	97	0	0	0	0	0	0	2. Bayer Automated (Clinitek Systems)	216	0	0	0	0	0	0
3. BMC Manual (Chemstrip)	35	0	0	0	0	0	0	4. BMC Automated (UA Readers)	50	0	0	0	0	0	0
10. Atlas System	51	0	0	0	0	0	1	12. Iris Clinitek Reader	8	0	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	0	15. Ictotest	12	0	0	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	0	0	99. All methods	492	0	0	0	0	0	1
Sample U65															
1. Bayer Test Strip Manual (Multistik)	96	0	1	0	0	0	0	2. Bayer Automated (Clinitek Systems)	216	0	0	0	0	0	0
3. BMC Manual (Chemstrip)	35	0	0	0	0	0	0	4. BMC Automated (UA Readers)	50	0	0	0	0	0	0
10. Atlas System	51	0	0	0	0	0	1	12. Iris Clinitek Reader	8	0	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	0	15. Ictotest	12	0	0	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	0	0	99. All methods	491	0	1	0	0	0	1

Urinalysis/Urine Pregnancy Proficiency Testing
14 February 2000

Hemoglobin

Method Code	1	2	3	4	5	Method Code	1	2	3	4	5
Sample U61											
1. Bayer Test Strip Manual (Multistik)	0	0	1	16	89	2. Bayer Automated (Clinitek Systems)	0	0	0	87	132
3. BMC Manual (Chemstrip)	0	0	1	7	36	4. BMC Automated (UA Readers)	0	0	1	13	37
10. Atlas Systems	0	0	0	14	39	12. Iris Clinitek Reader	0	0	0	1	7
13. Iris BMC UA Reader	0	0	0	2	3	17. Hitachi 747	0	0	0	1	1
28. Chemstrip Criterion	0	0	0	0	10	29. Hemastix	0	0	0	4	9
99. All methods	0	0	4	148	372						
Sample U62											
1. Bayer Test Strip Manual (Multistik)	105	0	0	0	0	2. Bayer Automated (Clinitek Systems)	218	0	0	1	0
3. BMC Manual (Chemstrip)	35	5	2	0	0	4. BMC Automated (UA Readers)	50	1	0	0	0
10. Atlas Systems	52	0	0	0	0	12. Iris Clinitek Reader	8	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	17. Hitachi 747	2	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	29. Hemastix	2	0	0	0	0
99. All methods	498	6	2	1	0						
Sample U63											
1. Bayer Test Strip Manual (Multistik)	102	3	0	0	0	2. Bayer Automated (Clinitek Systems)	217	1	0	0	1
3. BMC Manual (Chemstrip)	40	2	0	0	0	4. BMC Automated (UA Readers)	51	0	0	0	0
10. Atlas Systems	51	0	0	1	0	12. Iris Clinitek Reader	8	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	17. Hitachi 747	2	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	29. Hemastix	2	0	0	0	0
99. All methods	499	6	0	1	1						
Sample U64											
1. Bayer Test Strip Manual (Multistik)	1	0	0	9	96	2. Bayer Automated (Clinitek Systems)	1	0	2	43	173
3. BMC Manual (Chemstrip)	0	0	0	6	38	4. BMC Automated (UA Readers)	1	0	0	1	49
10. Atlas Systems	0	0	0	0	53	12. Iris Clinitek Reader	0	0	0	0	8
13. Iris BMC UA Reader	0	0	0	0	5	17. Hitachi 747	0	0	0	0	2
28. Chemstrip Criterion	0	0	0	0	10	29. Hemastix	0	0	0	2	11
99. All methods	3	0	3	62	456						
Sample U65											
1. Bayer Test Strip Manual (Multistik)	86	17	2	0	1	2. Bayer Automated (Clinitek Systems)	181	35	1	0	2
3. BMC Manual (Chemstrip)	33	6	2	1	0	4. BMC Automated (UA Readers)	39	4	8	0	0
10. Atlas Systems	37	15	0	0	0	12. Iris Clinitek Reader	7	0	1	0	0
13. Iris BMC UA Reader	5	0	0	0	0	17. Hitachi 747	1	1	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	29. Hemastix	2	0	0	0	0
99. All methods	411	79	14	1	3						

Urinalysis/Urine Pregnancy Proficiency Testing
14 February 2000

Ketones

Method Code	1	2	3	4	5	6	Method Code	1	2	3	4	5	6
Sample U61													
1. Bayer Test Strip Manual (Multistik)	105	1	0	0	0	0	2. Bayer Automated (Clinitek Systems)	219	0	0	0	0	0
3. BMC Manual (Chemstrip)	43	0	0	0	0	0	4. BMC Automated (UA Readers)	51	0	0	0	0	0
10. Atlas Systems	53	0	0	0	0	0	12. Iris Clinitek Reader	8	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	16. Acetest	6	0	0	0	0	0
23. Chemstrip UGK	0	0	0	0	0	0	28. Chemstrip Criterion	10	0	0	0	0	0
32. Ketostik	0	0	0	0	0	0	33. Keto-Diastik	0	0	0	0	0	0
99. All methods	510	1	0	0	0	0							
Sample U62													
1. Bayer Test Strip Manual (Multistik)	0	2	10	85	10	0	2. Bayer Automated (Clinitek Systems)	0	0	128	91	0	0
3. BMC Manual (Chemstrip)	0	0	13	28	3	1	4. BMC Automated (UA Readers)	0	0	1	49	1	0
10. Atlas Systems	0	0	23	30	0	0	12. Iris Clinitek Reader	0	0	2	6	0	0
13. Iris BMC UA Reader	0	0	4	1	0	0	16. Acetest	0	0	20	40	4	5
23. Chemstrip UGK	0	0	1	3	0	0	28. Chemstrip Criterion	0	0	2	8	0	0
32. Ketostik	0	0	5	4	0	0	33. Keto-Diastik	0	0	1	4	0	0
99. All methods	0	2	212	358	18	6							
Sample U63													
1. Bayer Test Strip Manual (Multistik)	1	2	20	77	8	0	2. Bayer Automated (Clinitek Systems)	0	1	170	48	0	0
3. BMC Manual (Chemstrip)	0	2	21	21	0	1	4. BMC Automated (UA Readers)	0	0	12	38	1	0
10. Atlas Systems	0	0	6	46	1	0	12. Iris Clinitek Reader	0	0	5	3	0	0
13. Iris BMC UA Reader	0	0	5	0	0	0	16. Acetest	0	0	35	27	2	5
23. Chemstrip UGK	0	0	2	2	0	0	28. Chemstrip Criterion	0	0	8	2	0	0
32. Ketostik	0	0	5	4	0	0	33. Keto-Diastik	0	0	2	3	0	0
99. All methods	1	4	298	275	12	6							
Sample U64													
1. Bayer Test Strip Manual (Multistik)	106	0	0	0	0	0	2. Bayer Automated (Clinitek Systems)	219	0	0	0	0	0
3. BMC Manual (Chemstrip)	43	0	0	0	0	0	4. BMC Automated (UA Readers)	51	0	0	0	0	0
10. Atlas Systems	52	0	0	0	0	0	12. Iris Clinitek Reader	8	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	16. Acetest	6	0	0	0	0	0
23. Chemstrip UGK	0	0	0	0	0	0	28. Chemstrip Criterion	10	0	0	0	0	0
32. Ketostik	0	0	0	0	0	0	33. Keto-Diastik	0	0	0	0	0	0
99. All methods	510	0	0	0	0	0							
Sample U65													
1. Bayer Test Strip Manual (Multistik)	105	0	0	1	0	0	2. Bayer Automated (Clinitek Systems)	219	0	0	0	0	0
3. BMC Manual (Chemstrip)	43	0	0	0	0	0	4. BMC Automated (UA Readers)	51	0	0	0	0	0
10. Atlas Systems	52	0	0	0	0	0	12. Iris Clinitek Reader	8	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	16. Acetest	7	0	0	0	0	0
23. Chemstrip UGK	0	0	0	0	0	0	28. Chemstrip Criterion	10	0	0	0	0	0
32. Ketostik	0	0	0	0	0	0	33. Keto-Diastik	0	0	0	0	0	0
99. All methods	510	0	0	1	0	0							

Protein

Method Code	1	2	3	4	5	6	7	8	Method Code	1	2	3	4	5	6	7	8
Sample U61																	
1. Bayer Test Strip Manual (Multistik)	0	0	0	5	76	11	2	14	2. Bayer Automated (Clinitek Systems)	0	1	0	16	196	2	2	0
3. BMC Manual (Chemstrip)	0	0	0	3	18	24	0	0	4. BMC Automated (UA Readers)	0	0	0	3	2	46	0	0
10. Atlas Systems	0	0	0	0	50	2	0	1	12. Iris Clinitek Reader	0	0	0	0	5	3	0	0
13. Iris BMC UA Reader	0	0	0	0	5	0	0	0	14. Sulfosalicylic Acid	0	0	4	9	98	22	0	4
17. Hitachi 747	0	0	0	1	2	0	0	0	25. Bayer Albustix	0	0	0	1	2	0	0	0
28. Chemstrip Criterion	0	0	0	0	0	10	0	0	30. Uristik	0	0	0	0	3	0	0	0
99. All methods	0	1	4	38	464	129	4	19									
Sample U62																	
1. Bayer Test Strip Manual (Multistik)	104	2	0	0	0	0	0	0	2. Bayer Automated (Clinitek Systems)	217	0	0	0	0	0	0	0
3. BMC Manual (Chemstrip)	42	1	0	0	0	0	0	0	4. BMC Automated (UA Readers)	50	1	0	0	0	0	0	0
10. Atlas Systems	52	0	0	0	0	0	1	0	12. Iris Clinitek Reader	8	0	0	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	0	0	14. Sulfosalicylic Acid	19	1	0	0	0	0	0	0
17. Hitachi 747	2	0	0	0	0	0	0	0	25. Bayer Albustix	0	0	0	0	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	0	0	0	30. Uristik	2	1	0	0	0	0	0	0
99. All methods	523	6	0	0	0	0	1	0									
Sample U63																	
1. Bayer Test Strip Manual (Multistik)	103	2	0	0	1	0	0	0	2. Bayer Automated (Clinitek Systems)	217	0	0	0	0	0	0	0
3. BMC Manual (Chemstrip)	42	1	0	0	0	0	0	0	4. BMC Automated (UA Readers)	51	0	0	0	0	0	0	0
10. Atlas Systems	52	0	0	0	0	0	1	0	12. Iris Clinitek Reader	8	0	0	0	0	0	0	0
13. Iris BMC UA Reader	5	0	0	0	0	0	0	0	14. Sulfosalicylic Acid	18	1	0	0	0	0	0	0
17. Hitachi 747	2	0	0	0	0	0	0	0	25. Bayer Albustix	0	0	0	0	0	0	0	0
28. Chemstrip Criterion	10	0	0	0	0	0	0	0	30. Uristik	3	0	0	0	0	0	0	0
99. All methods	523	4	0	0	1	0	1	0									
Sample U64																	
1. Bayer Test Strip Manual (Multistik)	0	0	0	0	44	24	6	34	2. Bayer Automated (Clinitek Systems)	0	0	0	1	211	2	2	0
3. BMC Manual (Chemstrip)	0	0	0	0	16	29	0	0	4. BMC Automated (UA Readers)	0	0	0	0	2	49	0	0
10. Atlas Systems	0	0	0	0	50	2	0	1	12. Iris Clinitek Reader	0	0	0	0	5	3	0	0
13. Iris BMC UA Reader	0	0	0	0	5	0	0	0	14. Sulfosalicylic Acid	0	0	1	6	81	40	4	6
17. Hitachi 747	0	0	0	0	1	2	0	0	25. Bayer Albustix	0	0	0	0	3	0	0	0
28. Chemstrip Criterion	0	0	0	0	0	10	0	0	30. Uristik	0	0	0	0	1	2	0	0
99. All methods	0	0	1	7	425	173	12	41									
Sample U65																	
1. Bayer Test Strip Manual (Multistik)	1	0	0	4	67	18	1	17	2. Bayer Automated (Clinitek Systems)	1	0	2	12	198	2	2	0
3. BMC Manual (Chemstrip)	0	0	0	3	16	26	0	0	4. BMC Automated (UA Readers)	0	0	0	8	2	41	0	0
10. Atlas Systems	1	0	0	0	49	2	0	1	12. Iris Clinitek Reader	0	0	0	1	5	2	0	0
13. Iris BMC UA Reader	0	0	0	0	5	0	0	0	14. Sulfosalicylic Acid	0	0	2	14	95	23	0	4
17. Hitachi 747	0	0	0	2	1	0	0	0	25. Bayer Albustix	0	0	0	0	3	0	0	0
28. Chemstrip Criterion	0	0	0	0	0	10	0	0	30. Uristik	0	0	0	0	3	0	0	0
99. All methods	3	0	4	46	450	132	3	22									

Urinalysis/Urine Pregnancy Proficiency Testing
14 February 2000

In Sample #	hCG Content in Sample	<u>Reported Results</u>		
		Positive	Indeterminate	Negative
U61	100 mIU/mL	436	0	3
U62	0 mIU/mL	1	0	438
U63	0 mIU/mL	2	0	437
U64	150 mIU/mL	436	1	2
U65	75 mIU/mL	432	2	4