

Participant Performance Summary - Initial (Screen) Test
Clinical Toxicology - March 2016 Test Event

Amphetamines Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
LC MS/MS	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Abbott Architect	1000	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Alere (Biosite) Triage Meter/PRO	1000	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	1000	Negative	11	11	100	Positive	11	11	100	Negative	11	11	100	Negative	11	11	100
Alere (Biosite) Triage TOX Drug Screen	1000	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
American Screening, Multi-panel Drug Test Cup	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	1000	Negative	18	18	100	Positive	18	18	100	Negative	18	18	100	Negative	18	18	100
BMC QuickTox	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
BioRad Tox/See	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Carolina Liquid Chemistries	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Diatron Pictus	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
E-Screen	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Ortho Vitros	1000	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
Roche Cobas	1000	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100	Negative	6	6	100
Roche KIMS	1000	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100
Siemens ADVIA	1000	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100
Siemens Dimension	1000	Negative	24	24	100	Positive	24	24	100	Negative	24	24	100	Negative	24	24	100
Siemens Syva Emit / II Plus	1000	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Thermo Scientific / Microgenics DRI	1000	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
Other	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	378	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Abbott Architect	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
AMS Diagnostics Liasys 450 device	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	500	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
MedTox Profile	500	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Roche Cobas	500	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100	Negative	2	2	100
Roche KIMS	500	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
Siemens Dimension	500	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100	Negative	6	6	100
Siemens Syva Emit / II Plus	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	500	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
All Labs	All Labs	Negative	141	141	100	Positive	141	141	100	Negative	141	141	100	Negative	141	141	100

Barbiturates Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Dimension	1000	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Abbott Architect	200	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100
AMS Diagnostics Liasys 450 device	200	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	200	Negative	21	21	100	Negative	21	21	100	Negative	21	21	100	Negative	21	21	100
Carolina Liquid Chemistries	200	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Diatron Pictus	200	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
MedTox Profile	200	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100
Ortho Vitros	200	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100
Roche Cobas	200	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100
Roche KIMS	200	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100
Siemens ADVIA	200	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Siemens Dimension	200	Negative	30	30	100	Negative	30	30	100	Negative	30	30	100	Negative	30	30	100
Siemens Syva Emit / II Plus	200	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100
Thermo Scientific / Microgenics DRI	200	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100
Other	200	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Alere (Biosite) Triage Meter/PRO	300	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	300	Negative	11	11	100	Negative	11	11	100	Negative	11	11	100	Negative	11	11	100
Alere (Biosite) Triage TOX Drug Screen	300	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100
American Screening, Multi-panel Drug Test Cup	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	300	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
BMC QuickTox	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
BioRad Tox/See	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens ADVIA	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Siemens Syva Emit / II Plus	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Thermo Scientific / Microgenics DRI	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Beckman Coulter	367	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	142	142	100	Negative	142	142	100	Negative	142	142	100	Negative	142	142	100

Benzodiazepines Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
LC MS/MS	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Roche Cobas	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Roche KIMS	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics CEDIA	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Dimension	1000	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Abbott Architect	150	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
MedTox Profile	150	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Abbott Architect	200	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
AMS Diagnostics Liasys 450 device	200	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	200	Negative	19	19	100	Negative	19	19	100	Positive	19	19	100	Negative	19	19	100
Carolina Liquid Chemistries	200	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Diatron Pictus	200	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ortho Vitros	200	Negative	4	4	100	Negative	4	4	100	Negative	4	1	25	Negative	4	4	100
										Positive	4	3	75				
Roche Cobas	200	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100
Roche KIMS	200	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Siemens ADVIA	200	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens Dimension	200	Negative	30	30	100	Negative	30	30	100	Positive	30	30	100	Negative	30	30	100
Siemens Syva Emit / II Plus	200	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
Thermo Scientific / Microgenics DRI	200	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Other	200	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Abbott Architect	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Alere (Biosite) Triage Meter/PRO	300	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	300	Negative	11	11	100	Negative	11	11	100	Positive	11	11	100	Negative	11	11	100
Alere (Biosite) Triage TOX Drug Screen	300	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
American Screening, Multi-panel Drug Test Cup	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	300	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
BMC QuickTox	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
BioRad Tox/See	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ortho Vitros	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Roche Cobas	300	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Roche KIMS	300	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens ADVIA	300	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens Syva Emit / II Plus	300	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Thermo Scientific / Microgenics DRI	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	350	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	142	142	100	Negative	142	142	100	Negative	142	1	0.7	Negative	142	142	100
										Positive	142	141	99				

Cocaine/metab		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
LC MS/MS	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Dimension	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Abbott Architect	150	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
AMS Diagnostics Liasys 450 device	150	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	150	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
MedTox Profile	150	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Roche Cobas	150	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Dimension	150	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
Thermo Scientific / Microgenics DRI	150	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100
Abbott Architect	300	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
Alere (Biosite) Triage Meter/PRO	300	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	300	Negative	10	10	100	Positive	10	10	100	Negative	10	10	100	Negative	10	10	100
Alere (Biosite) Triage TOX Drug Screen	300	Negative	8	8	100	Positive	8	8	100	Negative	8	8	100	Negative	8	8	100
American Screening, Multi-panel Drug Test Cup	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	300	Negative	20	20	100	Positive	20	20	100	Negative	20	20	100	Negative	20	20	100
BMC QuickTox	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
BioRad Tox/See	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Carolina Liquid Chemistries	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Diatron Pictus	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
E-Screen	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Ortho Vitros	300	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
Roche Cobas	300	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Roche KIMS	300	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Siemens ADVIA	300	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100
Siemens Dimension	300	Negative	25	25	100	Positive	25	25	100	Negative	25	25	100	Negative	25	25	100
Siemens Syva Emit / II Plus	300	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100	Negative	7	7	100
Thermo Scientific / Microgenics DRI	300	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
Other	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	399	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	143	143	100	Positive	143	143	100	Negative	143	143	100	Negative	143	143	100

Ethyl glucuronide Screen		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Carolina Liquid Chemistries	1000	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Immunalysis HEIA	250	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	250	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	500	Negative	4	2	50	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
		Positive	4	2	50												
Immunalysis HEIA	500	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens ADVIA	500	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	500	Negative	5	4	80	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100
		Positive	5	1	20												
All Labs	All Labs	Negative	14	7	50	Negative	14	14	100	Negative	14	14	100	Negative	14	14	100
		Positive	14	7	50												

Ethanol Screen (mg/dL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Abbott Architect	10	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Alcohol Dehydrogenases/UV or Visible Spectrome	10	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	10	Positive	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Diatron Pictus	10	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ortho Vitros	10	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Roche Cobas	10	Positive	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens ADVIA	10	Positive	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens Dimension	10	Positive	3	3	100	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100
Siemens Syva Emit / II Plus	10	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Alcohol Dehydrogenases/UV or Visible Spectrome	20	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	20	Positive	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
Siemens ADVIA	20	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Dimension	20	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	20	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	20	Positive	5	5	100	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100
Beckman Coulter	40	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens ADVIA	40	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Dimension	5	Positive	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Alcohol Dehydrogenases/UV or Visible Spectrome	50	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens Dimension	50	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100
Carolina Liquid Chemistries	70	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	44	7	16	Negative	44	44	100	Positive	44	44	100	Negative	44	44	100
		Positive	44	37	84												

Fentanyl Screen		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
IMEL - Immunalysis ELISA	0.7	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
LC MS/MS	10	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
AMS Diagnostics Liasys 450 device	2	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	2	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Immunalysis HEIA	2	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100
LC MS/MS	5	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	11	11	100	Negative	11	11	100	Negative	11	11	100	Negative	11	11	100

MDMA		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
LC MS/MS	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100	Negative	2	2	100
Alere (Biosite) Triage Meter/PRO	1000	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Alere (Biosite) Triage TOX Drug Screen	1500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	300	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
American Screening, Multi-panel Drug Test Cup	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	500	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
BioRad Tox/See	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Carolina Liquid Chemistries	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
MedTox Profile	500	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100	Negative	3	3	100
Siemens ADVIA	500	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	500	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100	Negative	2	2	100
Thermo Scientific / Microgenics DRI	500	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100	Negative	5	5	100
All Labs	All Labs	Negative	22	22	100	Positive	22	22	100	Negative	22	22	100	Negative	22	22	100

Methadone Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
LC MS/MS	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	100	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Dimension	1000	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
MedTox Profile	200	Negative	6	1	17	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100
		Positive	6	5	83												
Abbott Architect	300	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
Alere (Biosite) Triage Meter/PRO	300	Positive	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Alere (Biosite) Triage TOX Drug Screen	300	Positive	6	6	100	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100
American Screening, Multi-panel Drug Test Cup	300	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	300	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
AMS Diagnostics Liasys 450 device	300	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	300	Positive	21	21	100	Negative	21	21	100	Negative	21	21	100	Negative	21	21	100
Carolina Liquid Chemistries	300	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Diatron Pictus	300	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Ortho Vitros	300	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
Roche Cobas	300	Positive	6	6	100	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100
Roche KIMS	300	Positive	4	4	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
Siemens ADVIA	300	Positive	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Siemens Dimension	300	Positive	21	21	100	Negative	21	21	100	Negative	21	21	100	Negative	21	21	100
Siemens Syva Emit / II Plus	300	Positive	6	6	100	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100
Thermo Scientific / Microgenics DRI	300	Positive	6	6	100	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100
Other	300	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Ortho Vitros	3002	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	450	Positive	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	99	1	1	Negative	99	99	100	Negative	99	99	100	Negative	99	99	100
		Positive	99	98	99												

Opiates Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
LC MS/MS	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
MedTox Profile	100	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100
Siemens Syva Emit / II Plus	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
Siemens Dimension	1000	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
AMS Diagnostics Liasys 450 device	1999	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Alere (Biosite) Triage TOX Drug Screen	2000	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
American Screening, Multi-panel Drug Test Cup	2000	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Beckman Coulter	2000	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
E-Screen	2000	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Siemens ADVIA	2000	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Siemens Dimension	2000	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
Siemens Syva Emit / II Plus	2000	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
Beckman Coulter	270	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Abbott Architect	300	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100	Positive	8	8	100
Alere (Biosite) Triage Meter/PRO	300	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100
Alere (BioSite) Triage Plus TCA	300	Negative	11	11	100	Negative	11	11	100	Negative	11	11	100	Positive	11	11	100
Alere (Biosite) Triage TOX Drug Screen	300	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100
Ameditech ImmuTest Multi-Drug Screen	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Beckman Coulter	300	Negative	23	23	100	Negative	23	23	100	Negative	23	23	100	Positive	23	23	100
BMC QuickTox	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
BioRad Tox/See	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Carolina Liquid Chemistries	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Diatron Pictus	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Ortho Vitros	300	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100	Positive	5	5	100
Roche Cobas	300	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100	Positive	8	8	100
Roche KIMS	300	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100
Siemens ADVIA	300	Negative	3	3	100	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100
Siemens Dimension	300	Negative	28	28	100	Negative	28	28	100	Negative	28	28	100	Positive	28	28	100
Siemens Syva Emit / II Plus	300	Negative	5	5	100	Negative	5	5	100	Negative	5	5	100	Positive	5	5	100
Thermo Scientific / Microgenics DRI	300	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100
Other	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
All Labs	All Labs	Negative	143	143	100	Negative	143	143	100	Negative	143	143	100	Positive	143	143	100

Oxycodone Screen		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Abbott Architect	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
American Screening, Multi-panel Drug Test Cup	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
AMS Diagnostics Liasys 450 device	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Beckman Coulter	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100	Negative	4	4	100
BioRad Tox/See	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Immunalysis HEIA	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Innovacon Oxycodone Single Drug Test	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
LC MS/MS	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
MedTox Profile	100	Negative	9	9	100	Negative	9	9	100	Negative	9	9	100	Negative	9	9	100
Ortho Vitros	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Roche Cobas	100	Negative	3	3	100	Negative	3	3	100	Negative	3	3	100	Negative	3	3	100
Roche KIMS	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens ADVIA	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Dimension	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics CEDIA	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	100	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100	Negative	7	7	100
Other	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Alere (Biosite) Triage TOX Drug Screen	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Beckman Coulter	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
Carolina Liquid Chemistries	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Immunalysis HEIA	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Siemens Dimension	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100
All Labs	All Labs	Negative	46	46	100	Negative	46	46	100	Negative	46	46	100	Negative	46	45	98
														Positive	46	1	2.2

Phencyclidine (PCP) Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Dimension	1000	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Abbott Architect	25	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Alere (Biosite) Triage Meter/PRO	25	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	25	Negative	11	11	100	Negative	11	11	100	Positive	11	11	100	Negative	11	11	100
Alere (Biosite) Triage TOX Drug Screen	25	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
American Screening, Multi-panel Drug Test Cup	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
AMS Diagnostics Liasys 450 device	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	25	Negative	23	23	100	Negative	23	23	100	Positive	23	23	100	Negative	23	23	100
BMC QuickTox	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
BioRad Tox/See	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Diatron Pictus	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
E-Screen	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Immunalysis HEIA	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
LC MS/MS	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
MedTox Profile	25	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Ortho Vitros	25	Negative	5	5	100	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100
Roche Cobas	25	Negative	8	8	100	Negative	8	8	100	Positive	8	8	100	Negative	8	8	100
Roche KIMS	25	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Siemens ADVIA	25	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Siemens Dimension	25	Negative	27	27	100	Negative	27	27	100	Positive	27	27	100	Negative	27	27	100
Siemens Syva Emit / II Plus	25	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
Thermo Scientific / Microgenics DRI	25	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
Other	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	467	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	132	132	100	Negative	132	132	100	Positive	132	132	100	Negative	132	132	100

Propoxyphene Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Abbott Architect	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
American Screening, Multi-panel Drug Test Cup	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Beckman Coulter	300	Negative	8	8	100	Negative	8	8	100	Negative	8	8	100	Positive	8	8	100
Diatron Pictus	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
MedTox Profile	300	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100
Ortho Vitros	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Roche Cobas	300	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
Roche KIMS	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
Siemens ADVIA	300	Negative	2	2	100	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100
Siemens Syva Emit / II Plus	300	Negative	3	3	100	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100
Thermo Scientific / Microgenics DRI	300	Negative	6	6	100	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100
Beckman Coulter	361	Negative	1	1	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100
All Labs	All Labs	Negative	35	35	100	Negative	35	35	100	Negative	35	35	100	Positive	35	35	100

Cannabinoids (THCA) Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
Siemens Syva Emit / II Plus	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Dimension	1000	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	20	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Syva Emit / II Plus	20	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	20	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Siemens Dimension	25	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100
Thermo Scientific / Microgenics CEDIA	25	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	409	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Abbott Architect	50	Negative	8	8	100	Negative	8	8	100	Positive	8	8	100	Negative	8	8	100
Alere (Biosite) Triage Meter/PRO	50	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	50	Negative	11	11	100	Negative	11	11	100	Positive	11	11	100	Negative	11	11	100
Alere (Biosite) Triage TOX Drug Screen	50	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
American Screening, Multi-panel Drug Test Cup	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	50	Negative	24	24	100	Negative	24	24	100	Positive	24	24	100	Negative	24	24	100
BMC QuickTox	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
BioRad Tox/See	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Carolina Liquid Chemistries	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Diatron Pictus	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
E-Screen	50	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
MedTox Profile	50	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Ortho Vitros	50	Negative	5	5	100	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100
Roche Cobas	50	Negative	8	8	100	Negative	8	8	100	Positive	8	8	100	Negative	8	8	100
Roche KIMS	50	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Siemens ADVIA	50	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Siemens Dimension	50	Negative	27	27	100	Negative	27	27	100	Positive	27	27	100	Negative	27	27	100
Siemens Syva Emit / II Plus	50	Negative	7	7	100	Negative	7	7	100	Positive	7	7	100	Negative	7	7	100
Thermo Scientific / Microgenics DRI	50	Negative	5	5	100	Negative	5	5	100	Positive	5	5	100	Negative	5	5	100
All Labs	All Labs	Negative	141	141	100	Negative	141	141	100	Positive	141	141	100	Negative	141	141	100

Tricyclic Antidepressants Screen (ng/mL)		UTC 1				UTC 2				UTC 3				UTC 4			
Method Code	Cutoff	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%	Response	No. labs	No. Resp.	%
IMEL - Immunalysis ELISA	100	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Abbott Architect	1000	Negative	3	3	100	Negative	3	3	100	Negative	3	3	100	Negative	3	3	100
Alere (Biosite) Triage Meter/PRO	1000	Negative	4	4	100	Negative	4	4	100	Positive	4	4	100	Negative	4	4	100
Alere (BioSite) Triage Plus TCA	1000	Negative	10	10	100	Negative	10	10	100	Positive	10	10	100	Negative	10	10	100
Alere (Biosite) Triage TOX Drug Screen	1000	Negative	6	6	100	Negative	6	6	100	Positive	6	6	100	Negative	6	6	100
American Screening, Multi-panel Drug Test Cup	1000	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Ameditech ImmuTest Multi-Drug Screen	1000	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
BioRad Tox/See	1000	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Siemens Syva Rapid Test	1000	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100
DIPR - Diagnostic Products	110	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Abbott Architect	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Alere (BioSite) Triage Plus TCA	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Beckman Coulter	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Carolina Liquid Chemistries	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
MedTox Profile	300	Negative	7	7	100	Negative	7	7	100	Negative	7	1	14	Negative	7	7	100
										Positive	7	6	86				
Ortho Vitros	300	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
Thermo Scientific / Microgenics DRI	300	Negative	3	3	100	Negative	3	3	100	Positive	3	3	100	Negative	3	3	100
Beckman Coulter	500	Negative	2	2	100	Negative	2	2	100	Positive	2	2	100	Negative	2	2	100
Roche Cobas	500	Negative	1	1	100	Negative	1	1	100	Positive	1	1	100	Negative	1	1	100
All Labs	All Labs	Negative	50	50	100	Negative	50	50	100	Negative	50	4	8	Negative	50	50	100
										Positive	50	46	92				