

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Specific Gravity

Specimen: U76 -----	Specimen: U77 -----	Specimen: U78 -----	Specimen: U79 -----	Specimen: U80 -----	Number -----	Instrument or Reagent System -----
1.018 ± 0.004 1.020	1.019 ± 0.004 1.020	1.015 ± 0.003 1.015	1.013 ± 0.002 1.014	1.019 ± 0.003 1.020	n = 382 n = 382	All Methods & Instruments Median value, All Methods & Instruments
<Instruments>						
1.014 ± 0.003	1.019 ± 0.004	1.016 ± 0.004	1.010 ± 0.003	1.019 ± 0.004	n = 34	Visual examination
1.020 ± 0.002	1.020 ± 0.002	1.013 ± 0.002	1.015 ± 0.002	1.020 ± 0.002	n = 87	Siemens Atlas
1.022 ± 0.003	1.021 ± 0.002	1.017 ± 0.003	1.017 ± 0.003	1.021 ± 0.002	n = 29	Siemens Clinitek Status
1.019 ± 0.002	1.018 ± 0.003	1.016 ± 0.002	1.010 ± 0.001	1.018 ± 0.003	n = 10	Siemens Clinitek 50
1.022 ± 0.002	1.021 ± 0.002	1.018 ± 0.002	1.014 ± 0.003	1.021 ± 0.002	n = 13	Siemens Clinitek 100
1.020 ± 0.002	1.020 ± 0.002	1.018 ± 0.003	1.012 ± 0.003	1.020 ± 0.003	n = 70	Siemens Clinitek 500
1.020 ± 0.001	1.019 ± 0.002	1.017 ± 0.003	1.012 ± 0.003	1.020 ± 0.002	n = 29	Siemens Clinitek Advantus
1.011 ± 0.002	1.016 ± 0.002	1.017 ± 0.003	1.011 ± 0.002	1.016 ± 0.002	n = 5	Roche UA Reader/Super UA Reader
1.011 ± 0.003	1.017 ± 0.002	1.017 ± 0.002	1.011 ± 0.002	1.015 ± 0.001	n = 13	Roche Chemstrip Criterion
1.016 ± 0.004	1.019 ± 0.002	1.015 ± 0.004	1.014 ± 0.002	1.018 ± 0.002	n = 36	Roche Urisys
1.016 ± 0.000	1.013 ± 0.001	1.010 ± 0.000	1.011 ± 0.001	1.016 ± 0.001	n = 29	IRIS AUTION MAX
1.019 ± 0.001	1.019 ± 0.002	1.011 ± 0.001	1.014 ± 0.001	1.019 ± 0.001	n = 16	Refractometer
1.015 ± 0.003	1.017 ± 0.003	1.016 ± 0.003	1.012 ± 0.002	1.017 ± 0.002	n = 8	Other
<Reagents>						
1.020 ± 0.003	1.020 ± 0.003	1.017 ± 0.003	1.013 ± 0.004	1.020 ± 0.003	n = 172	Siemens Reagent Strips
1.020 ± 0.001	1.020 ± 0.001	1.013 ± 0.002	1.015 ± 0.001	1.020 ± 0.001	n = 75	Siemens Atlas Reagent
1.012 ± 0.003	1.017 ± 0.002	1.017 ± 0.003	1.012 ± 0.003	1.016 ± 0.002	n = 48	Roche Chemstrips
1.019 ± 0.001	1.019 ± 0.001	1.012 ± 0.001	1.014 ± 0.001	1.019 ± 0.001	n = 17	Roche Urisys Cassette
1.016 ± 0.001	1.013 ± 0.001	1.010 ± 0.000	1.011 ± 0.001	1.017 ± 0.001	n = 16	ARKRAY Aution Sticks
1.015 ± 0.000	1.013 ± 0.001	1.010 ± 0.000	1.011 ± 0.002	1.016 ± 0.002	n = 8	IRIS Reagent Strips
1.017 ± 0.005	1.020 ± 0.003	1.015 ± 0.004	1.013 ± 0.005	1.020 ± 0.003	n = 15	Bayer Clinitest
1.018 ± 0.002	1.017 ± 0.003	1.011 ± 0.002	1.013 ± 0.002	1.019 ± 0.002	n = 22	None, no reagent or strip used
1.017 ± 0.004	1.020 ± 0.006	1.015 ± 0.005	1.013 ± 0.003	1.019 ± 0.004	n = 9	Other

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

pH

Specimen: U76	Specimen: U77	Specimen: U78	Specimen: U79	Specimen: U80	Number	Instrument or Reagent System
7.15 ± 0.25	6.12 ± 0.30	5.70 ± 0.39	7.00 ± 0.14	5.89 ± 0.33	n = 390	All Methods & Instruments
7.00	6.00	6.00	7.00	6.00	n = 390	Median value, All Methods & Instruments
<Instruments>						
7.20 ± 0.31	6.09 ± 0.27	5.60 ± 0.48	7.06 ± 0.16	5.78 ± 0.41	n = 45	Visual examination
7.29 ± 0.26	6.35 ± 0.23	5.96 ± 0.14	7.00 ± 0.00	6.01 ± 0.12	n = 86	Siemens Atlas
7.02 ± 0.09	6.17 ± 0.27	5.61 ± 0.25	7.00 ± 0.00	5.92 ± 0.22	n = 32	Siemens Clinitek Status
7.05 ± 0.15	5.95 ± 0.15	5.64 ± 0.23	7.00 ± 0.00	5.91 ± 0.20	n = 11	Siemens Clinitek 50
7.50 ± 0.00	6.00 ± 0.00	5.54 ± 0.24	7.18 ± 0.25	5.82 ± 0.32	n = 14	Siemens Clinitek 100
7.15 ± 0.28	6.01 ± 0.12	5.90 ± 0.22	6.95 ± 0.17	5.97 ± 0.15	n = 74	Siemens Clinitek 500
7.15 ± 0.27	6.00 ± 0.13	5.90 ± 0.24	6.97 ± 0.13	5.90 ± 0.24	n = 30	Siemens Clinitek Advantus
7.00 ± 0.00	6.20 ± 0.27	5.70 ± 0.67	7.00 ± 0.00	6.20 ± 0.27	n = 5	Roche UA Reader/Super UA Reader
7.00 ± 0.00	6.14 ± 0.23	5.43 ± 0.51	7.00 ± 0.00	6.11 ± 0.21	n = 14	Roche Chemstrip Criterion
7.00 ± 0.00	6.04 ± 0.67	5.13 ± 0.37	7.00 ± 0.00	5.56 ± 0.62	n = 36	Roche Urisys
7.00 ± 0.00	6.00 ± 0.00	5.52 ± 0.09	6.95 ± 0.28	5.76 ± 0.34	n = 29	IRIS AUTION MAX
7.06 ± 0.17	6.06 ± 0.17	5.44 ± 0.39	7.00 ± 0.00	5.78 ± 0.26	n = 9	Other
<Reagents>						
7.17 ± 0.27	6.06 ± 0.19	5.81 ± 0.28	7.00 ± 0.16	5.92 ± 0.21	n = 186	Siemens Reagent Strips
7.29 ± 0.26	6.36 ± 0.23	5.97 ± 0.12	7.00 ± 0.00	6.01 ± 0.08	n = 76	Siemens Atlas Reagent
7.04 ± 0.13	6.02 ± 0.35	5.25 ± 0.45	7.00 ± 0.00	5.78 ± 0.51	n = 54	Roche Chemstrips
7.00 ± 0.00	6.15 ± 0.84	5.09 ± 0.36	7.00 ± 0.00	5.47 ± 0.67	n = 17	Roche Urisys Cassette
7.00 ± 0.00	6.00 ± 0.00	5.52 ± 0.11	7.00 ± 0.00	5.71 ± 0.25	n = 21	ARKRAY Aution Sticks
7.00 ± 0.00	6.00 ± 0.00	5.50 ± 0.00	6.83 ± 0.50	5.83 ± 0.50	n = 9	IRIS Reagent Strips
7.23 ± 0.26	6.23 ± 0.26	5.73 ± 0.37	7.03 ± 0.13	5.97 ± 0.40	n = 15	Bayer Clinitest
6.98 ± 0.23	6.01 ± 0.03	5.62 ± 0.38	7.03 ± 0.17	5.82 ± 0.25	n = 11	Other

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Bilirubin	Specimen: U76	Negative	Small 1+ 1 mg/dL	Moderate 2+ 3 mg/dL	Large 3+ 6 mg/dL	Positive
<Instruments>						
Visual examination	73	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	28	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0
Siemens Clinitek 500	71	1	0	0	0	0
Siemens Clinitek Advantus	31	0	0	0	0	0
Roche UA Reader/Super UA Read	4	0	0	0	0	0
Roche Mini UA Reader	2	0	0	0	0	0
Roche Chemstrip Criterion	13	0	0	0	0	0
Roche Urisys	35	0	0	0	0	0
IRIS AUTION MAX	28	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0
Other	11	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	176	1	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0
Roche Chemstrips	51	0	0	0	0	0
Roche Urisys Cassette	15	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0
Bayer Clinitest	17	0	0	0	0	0
Siemens Ictotest	36	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Other	11	0	0	0	0	0
All Instruments and Reagents	413	1	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Bilirubin	Specimen: U77	Negative	Small 1+ 1 mg/dL	Moderate 2+ 3 mg/dL	Large 3+ 6 mg/dL	Positive
<Instruments>						
Visual examination	72	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	28	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0
Siemens Clinitek Advantus	31	0	0	0	0	0
Roche UA Reader/Super UA Read	4	0	0	0	0	0
Roche Mini UA Reader	2	0	0	0	0	0
Roche Chemstrip Criterion	13	0	0	0	0	0
Roche Urisys	35	0	0	0	0	0
IRIS AUTION MAX	28	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0
Other	11	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	177	0	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0
Roche Chemstrips	51	0	0	0	0	0
Roche Urisys Cassette	15	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0
Bayer Clinitest	17	0	0	0	0	0
Siemens Ictotest	35	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Other	11	0	0	0	0	0
All Instruments and Reagents	413	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Bilirubin	Specimen: U78	Negative	Small 1+ 1 mg/dL	Moderate 2+ 3 mg/dL	Large 3+ 6 mg/dL	Positive
<Instruments>						
Visual examination	73	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	28	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0
Siemens Clinitek 500	71	1	0	0	0	0
Siemens Clinitek Advantus	31	0	0	0	0	0
Roche UA Reader/Super UA Read	4	0	0	0	0	0
Roche Mini UA Reader	2	0	0	0	0	0
Roche Chemstrip Criterion	13	0	0	0	0	0
Roche Urisys	35	0	0	0	0	0
IRIS AUTION MAX	28	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0
Other	11	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	176	1	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0
Roche Chemstrips	51	0	0	0	0	0
Roche Urisys Cassette	15	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0
Bayer Clinitest	17	0	0	0	0	0
Siemens Ictotest	36	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Other	11	0	0	0	0	0
All Instruments and Reagents	413	1	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Bilirubin	Specimen: U79	Negative	Small 1+ 1 mg/dL	Moderate 2+ 3 mg/dL	Large 3+ 6 mg/dL	Positive
<Instruments>						
Visual examination	72	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0
Siemens Clinitek Status	28	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0
Siemens Clinitek Advantus	31	0	0	0	0	0
Roche UA Reader/Super UA Read	4	0	0	0	0	0
Roche Mini UA Reader	2	0	0	0	0	0
Roche Chemstrip Criterion	13	0	0	0	0	0
Roche Urisys	35	0	0	0	0	0
IRIS AUTION MAX	27	0	1	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0
Other	11	0	0	0	0	0
<Reagents>						
Siemens Reagent Strips	177	0	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0
Roche Chemstrips	51	0	0	0	0	0
Roche Urisys Cassette	15	0	0	0	0	0
ARKRAY Aution Sticks	20	0	0	0	0	0
IRIS Reagent Strips	8	0	1	0	0	0
Bayer Clinitest	17	0	0	0	0	0
Siemens Ictotest	35	0	0	0	0	0
Siemens Hemastix	1	0	0	0	0	0
Other	11	0	0	0	0	0
All Instruments and Reagents	412	0	1	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Bilirubin	Specimen: U80	Negative	Small 1+ 1 mg/dL	Moderate 2+ 3 mg/dL	Large 3+ 6 mg/dL	Positive
<Instruments>						
Visual examination		3	10	31	22	134
Siemens Atlas		0	45	41	0	3
Siemens Clinitek Status		0	10	18	0	0
Siemens Clinitek 50		0	1	10	0	0
Siemens Clinitek 100		0	0	11	4	1
Siemens Clinitek 500		0	14	57	0	3
Siemens Clinitek Advantus		0	1	30	0	1
Roche UA Reader/Super UA Read		2	2	0	0	0
Roche Mini UA Reader		0	2	0	0	0
Roche Chemstrip Criterion		9	0	3	0	1
Roche Urisys		1	33	0	0	1
IRIS AUTION MAX		1	1	26	0	1
Other Chemistry Analyzer		0	0	2	0	0
Other		0	4	6	0	3
<Reagents>						
Siemens Reagent Strips		1	27	132	16	4
Siemens Atlas Reagent		0	41	35	0	2
Roche Chemstrips		12	25	11	2	2
Roche Urisys Cassette		1	14	0	0	0
ARKRAY Aution Sticks		0	1	19	0	0
IRIS Reagent Strips		1	0	8	0	1
Bayer Clinitest		0	4	9	4	0
Siemens Ictotest		0	9	15	2	137
Siemens Hemastix		0	0	0	0	1
Other		1	2	6	2	1
All Instruments and Reagents		16	123	235	26	148

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Glucose	Specimen: U76	Negative	Trace 50-100 mg/dL	1/4% 150-300 mg/dL	1/2% 500 mg/dL	3/4% 750 mg/dL	1% >=1000 mg/dL	2% >=2000 mg/dL
<Instruments>								
Visual examination		97	2	1	0	0	0	0
Siemens Atlas		88	0	0	0	0	0	0
Siemens Clinitek Status		32	0	0	0	0	0	0
Siemens Clinitek 50		11	0	0	0	0	0	0
Siemens Clinitek 100		14	2	0	0	0	0	0
Siemens Clinitek 500		72	0	1	0	0	0	0
Siemens Clinitek Advantus		30	0	0	0	0	0	0
Roche UA Reader/Super UA Read		5	0	0	0	0	0	0
Roche Mini UA Reader		1	0	0	0	0	0	0
Roche Chemstrip Criterion		14	0	0	0	0	0	0
Roche Urisys		37	0	0	0	0	0	0
IRIS AUTION MAX		29	0	0	0	0	0	0
Hitachi/Roche Analyzer		2	0	0	0	0	0	0
Other Chemistry Analyzer		3	0	0	0	0	0	0
Other		11	0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips		187	1	1	0	0	0	0
Siemens Atlas Reagent		76	0	0	0	0	0	0
Sulfosalicylic Acid		1	0	0	0	0	0	0
Roche Chemstrips		55	0	0	0	0	0	0
Roche Urisys Cassette		17	0	0	0	0	0	0
ARKRAY Aution Sticks		21	0	0	0	0	0	0
IRIS Reagent Strips		9	0	0	0	0	0	0
Siemens Acetest		1	0	0	0	0	0	0
Bayer Clinitest		63	3	1	0	0	0	0
Hitachi/Roche Reagents		2	0	0	0	0	0	0
Other		14	0	0	0	0	0	0
All Instruments and Reagents		446	4	2	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Glucose	Specimen: U77	Negative	Trace 50-100 mg/dL	1/4% 150-300 mg/dL	1/2% 500 mg/dL	3/4% 750 mg/dL	1% >=1000 mg/dL	2% >=2000 mg/dL
<Instruments>								
Visual examination		2	0	0	9	11	114	23
Siemens Atlas		0	0	0	3	9	76	0
Siemens Clinitek Status		0	0	1	20	1	10	0
Siemens Clinitek 50		0	0	0	0	0	11	0
Siemens Clinitek 100		0	0	0	0	4	12	0
Siemens Clinitek 500		0	0	1	8	2	63	0
Siemens Clinitek Advantus		0	0	0	2	1	27	0
Roche UA Reader/Super UA Read		0	0	0	0	0	5	0
Roche Mini UA Reader		0	0	0	0	0	1	0
Roche Chemstrip Criterion		0	0	2	2	0	10	0
Roche Urisys		0	0	2	0	2	33	0
IRIS AUTION MAX		0	0	0	0	0	29	0
Hitachi/Roche Analyzer		0	0	0	0	0	2	0
Other Chemistry Analyzer		0	0	0	0	0	3	0
Other		0	0	0	1	2	8	0
<Reagents>								
Siemens Reagent Strips		0	0	1	30	10	142	8
Siemens Atlas Reagent		0	0	0	3	6	67	0
Sulfosalicylic Acid		0	0	0	1	0	0	0
Roche Chemstrips		0	0	3	4	1	51	1
Roche Urisys Cassette		0	0	1	0	1	15	0
ARKRAY Aution Sticks		0	0	0	0	0	21	0
IRIS Reagent Strips		0	0	0	0	0	9	0
Siemens Acetest		0	0	0	0	0	1	0
Bayer Clinitest		1	0	1	5	13	87	10
Bayer Diastix		0	0	0	0	0	3	1
Hitachi/Roche Reagents		0	0	0	0	0	2	0
Other		1	0	0	2	1	6	3
All Instruments and Reagents		2	0	6	45	32	404	23

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Glucose	Specimen: U78	Negative	Trace 50-100 mg/dL	1/4% 150-300 mg/dL	1/2% 500 mg/dL	3/4% 750 mg/dL	1% >=1000 mg/dL	2% >=2000 mg/dL
<Instruments>								
Visual examination	97	0	0	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0	0
Siemens Clinitek Status	32	0	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0	0	0
Siemens Clinitek 500	72	1	0	0	0	0	0	0
Siemens Clinitek Advantus	30	0	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0	0
Roche Urisys	37	0	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0	0
Hitachi/Roche Analyzer	2	0	0	0	0	0	0	0
Other Chemistry Analyzer	3	0	0	0	0	0	0	0
Other	11	0	0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips	188	1	0	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0	0	0
Sulfosalicylic Acid	1	0	0	0	0	0	0	0
Roche Chemstrips	55	0	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0	0
Siemens Acetest	1	0	0	0	0	0	0	0
Bayer Clinitest	65	0	0	0	0	0	0	0
Hitachi/Roche Reagents	2	0	0	0	0	0	0	0
Other	13	0	0	0	0	0	0	0
All Instruments and Reagents	448	1	0	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Glucose	Specimen: U79	Negative	Trace 50-100 mg/dL	1/4% 150-300 mg/dL	1/2% 500 mg/dL	3/4% 750 mg/dL	1% >=1000 mg/dL	2% >=2000 mg/dL
<Instruments>								
Visual examination		2	9	113	29	0	1	1
Siemens Atlas		0	2	85	0	1	0	0
Siemens Clinitek Status		0	1	27	4	0	0	0
Siemens Clinitek 50		0	0	5	5	0	1	0
Siemens Clinitek 100		0	0	16	0	0	0	0
Siemens Clinitek 500		0	2	66	6	0	0	0
Siemens Clinitek Advantus		0	0	28	2	0	0	0
Roche UA Reader/Super UA Read		0	0	5	0	0	0	0
Roche Mini UA Reader		0	0	1	0	0	0	0
Roche Chemstrip Criterion		0	0	13	1	0	0	0
Roche Urisys		0	2	32	3	0	0	0
IRIS AUTION MAX		2	0	19	6	2	0	0
Hitachi/Roche Analyzer		0	0	2	0	0	0	0
Other Chemistry Analyzer		0	0	3	0	0	0	0
Other		0	0	7	4	0	0	0
<Reagents>								
Siemens Reagent Strips		0	4	165	20	0	1	0
Siemens Atlas Reagent		0	2	74	0	0	0	0
Sulfosalicylic Acid		0	0	1	0	0	0	0
Roche Chemstrips		0	3	53	4	0	0	0
Roche Urisys Cassette		0	0	15	2	0	0	0
ARKRAY Aution Sticks		1	0	14	5	1	0	0
IRIS Reagent Strips		1	0	6	1	1	0	0
Siemens Acetest		0	0	1	0	0	0	0
Bayer Clinitest		2	7	77	26	1	1	0
Bayer Diastix		0	0	4	0	0	0	0
Hitachi/Roche Reagents		0	0	2	0	0	0	0
Other		0	0	10	2	0	0	1
All Instruments and Reagents		4	16	422	60	3	2	1

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Glucose	Specimen: U80	Negative	Trace 50-100 mg/dL	1/4% 150-300 mg/dL	1/2% 500 mg/dL	3/4% 750 mg/dL	1% >=1000 mg/dL	2% >=2000 mg/dL
<Instruments>								
Visual examination		89	26	0	0	0	0	0
Siemens Atlas		88	0	0	0	0	0	0
Siemens Clinitek Status		19	13	0	0	0	0	0
Siemens Clinitek 50		1	10	0	0	0	0	0
Siemens Clinitek 100		16	0	0	0	0	0	0
Siemens Clinitek 500		41	33	0	0	0	0	0
Siemens Clinitek Advantus		22	8	0	0	0	0	0
Roche UA Reader/Super UA Read		2	3	0	0	0	0	0
Roche Mini UA Reader		0	1	0	0	0	0	0
Roche Chemstrip Criterion		3	11	0	0	0	0	0
Roche Urisys		17	20	0	0	0	0	0
IRIS AUTION MAX		27	1	1	0	0	0	0
Hitachi/Roche Analyzer		0	1	0	0	1	0	0
Other Chemistry Analyzer		1	2	0	0	0	0	0
Other		8	3	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips		123	66	0	0	0	0	0
Siemens Atlas Reagent		76	0	0	0	0	0	0
Sulfosalicylic Acid		1	0	0	0	0	0	0
Roche Chemstrips		12	47	0	0	0	0	0
Roche Urisys Cassette		14	3	0	0	0	0	0
ARKRAY Aution Sticks		20	1	0	0	0	0	0
IRIS Reagent Strips		8	0	1	0	0	0	0
Siemens Acetest		1	0	0	0	0	0	0
Bayer Clinitest		71	7	0	0	0	0	0
Bayer Diastix		0	2	0	0	0	0	0
Hitachi/Roche Reagents		0	1	0	0	1	0	0
Other		8	5	0	0	0	0	0
All Instruments and Reagents		334	132	1	0	1	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Hemoglobin	Specimen: U76	Negative	Trace 5-10 ery/uL	Small 1+ 25 ery/uL	Moderate 2+ 50 ery/uL	150 ery/uL	Large 3+ 250 ery/uL
<Instruments>							
Visual examination		0	0	0	1	1	54
Siemens Atlas		1	0	0	4	0	83
Siemens Clinitek Status		0	0	0	1	0	32
Siemens Clinitek 50		0	0	0	0	0	11
Siemens Clinitek 100		0	0	0	0	0	15
Siemens Clinitek 500		1	0	0	1	0	70
Siemens Clinitek Advantus		0	0	0	0	0	29
Roche UA Reader/Super UA Read		0	0	0	0	0	5
Roche Mini UA Reader		0	0	0	0	0	1
Roche Chemstrip Criterion		1	0	0	0	1	12
Roche Urisys		0	0	0	5	1	30
IRIS AUTION MAX		0	0	0	16	0	13
Hitachi/Roche Analyzer		0	0	0	0	0	1
Other Chemistry Analyzer		0	0	0	0	0	2
Other		0	0	0	1	0	8
<Reagents>							
Siemens Reagent Strips		1	0	0	3	1	183
Siemens Atlas Reagent		1	0	0	3	0	72
Roche Chemstrips		1	0	0	4	2	51
Roche Urisys Cassette		0	0	0	2	0	15
ARKRAY Aution Sticks		0	0	0	14	0	7
IRIS Reagent Strips		0	0	0	3	0	6
Bayer Clinitest		0	0	0	0	0	15
Siemens Hemastix		0	0	0	0	0	4
None, no reagent or strip use		0	0	0	0	0	1
Other		0	0	0	0	0	12
All Instruments and Reagents		3	0	0	29	3	366

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Hemoglobin	Specimen: U77	Negative	Trace 5-10 ery/uL	Small 1+ 25 ery/uL	Moderate 2+ 50 ery/uL	150 ery/uL	Large 3+ 250 ery/uL
<Instruments>							
Visual examination	47	2	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	32	1	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	15	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	28	1	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	35	1	0	0	0	0	0
IRIS AUTION MAX	28	1	0	0	0	0	0
Hitachi/Roche Analyzer	0	0	0	0	0	1	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	185	2	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0	0
Roche Chemstrips	52	3	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	20	1	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Hemastix	2	0	0	0	0	0	0
None, no reagent or strip use	1	0	0	0	0	0	0
Other	10	0	0	0	0	1	0
All Instruments and Reagents	387	6	0	0	0	1	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Hemoglobin	Specimen: U78	Negative	Trace 5-10 ery/uL	Small 1+ 25 ery/uL	Moderate 2+ 50 ery/uL	150 ery/uL	Large 3+ 250 ery/uL
<Instruments>							
Visual examination	49	0	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	33	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	15	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	29	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	36	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Hitachi/Roche Analyzer	0	0	0	0	0	1	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	187	0	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0	0
Roche Chemstrips	55	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Hemastix	2	0	0	0	0	0	0
None, no reagent or strip use	1	0	0	0	0	0	0
Other	10	0	0	0	0	1	0
All Instruments and Reagents	393	0	0	0	0	1	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Hemoglobin	Specimen: U79	Negative	Trace 5-10 ery/uL	Small 1+ 25 ery/uL	Moderate 2+ 50 ery/uL	150 ery/uL	Large 3+ 250 ery/uL
<Instruments>							
Visual examination	47	2	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	33	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	15	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	29	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	35	1	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Hitachi/Roche Analyzer	0	0	0	0	0	1	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	187	0	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0	0
Roche Chemstrips	52	3	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Hemastix	2	0	0	0	0	0	0
None, no reagent or strip use	1	0	0	0	0	0	0
Other	10	0	0	0	0	1	0
All Instruments and Reagents	390	3	0	0	0	1	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Hemoglobin	Specimen: U80	Negative	Trace 5-10 ery/uL	Small 1+ 25 ery/uL	Moderate 2+ 50 ery/uL	150 ery/uL	Large 3+ 250 ery/uL
<Instruments>							
Visual examination	49	0	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	30	3	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	15	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	29	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	36	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Hitachi/Roche Analyzer	0	0	0	0	0	0	1
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	184	3	0	0	0	0	0
Siemens Atlas Reagent	76	0	0	0	0	0	0
Roche Chemstrips	55	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Hemastix	2	0	0	0	0	0	0
None, no reagent or strip use	1	0	0	0	0	0	0
Other	10	0	0	0	0	0	1
All Instruments and Reagents	390	3	0	0	0	0	1

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Ketones	Specimen: U76	Negative	Trace 5 mg/dL	Small 1+ 15 mg/dL	Moderate 2+ 40-50 mg/dL	Large 3+ >=80 mg/dL	Positive
<Instruments>							
Visual examination		0	0	17	62	10	13
Siemens Atlas		0	1	9	78	0	0
Siemens Clinitek Status		0	0	2	28	2	0
Siemens Clinitek 50		0	1	0	10	0	0
Siemens Clinitek 100		0	0	6	9	0	1
Siemens Clinitek 500		1	0	45	26	0	0
Siemens Clinitek Advantus		0	0	21	9	0	0
Roche UA Reader/Super UA Read		0	0	0	5	0	0
Roche Mini UA Reader		0	0	1	0	0	0
Roche Chemstrip Criterion		0	1	0	13	0	0
Roche Urisys		1	0	0	33	2	0
IRIS AUTION MAX		0	0	0	29	0	0
Other Chemistry Analyzer		0	0	0	2	0	0
Other		0	0	0	9	0	0
<Reagents>							
Siemens Reagent Strips		0	1	71	110	3	0
Siemens Atlas Reagent		0	1	8	68	0	0
Roche Chemstrips		0	1	3	49	3	0
Roche Urisys Cassette		1	0	0	15	1	0
ARKRAY Aution Sticks		0	0	0	21	0	0
IRIS Reagent Strips		0	0	0	9	0	0
Siemens Acetest		1	0	14	17	4	14
Bayer Clinitest		0	0	3	11	1	0
Siemens Ictotest		0	0	0	1	0	0
Bayer Diastix		0	0	0	4	0	0
Other		0	0	2	8	2	0
All Instruments and Reagents		2	3	101	313	14	14

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

Ketones	Specimen: U77	Negative	Trace 5 mg/dL	Small 1+ 15 mg/dL	Moderate 2+ 40-50 mg/dL	Large 3+ >=80 mg/dL	Positive
<Instruments>							
Visual examination	61	0	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	32	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	30	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	36	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	184	0	0	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0	0
Roche Chemstrips	54	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Siemens Acetest	16	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Ictotest	1	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Other	11	0	0	0	0	0	0
All Instruments and Reagents	406	0	0	0	0	0	0

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

Ketones	Specimen: U78	Negative	Trace 5 mg/dL	Small 1+ 15 mg/dL	Moderate 2+ 40-50 mg/dL	Large 3+ >=80 mg/dL	Positive
<Instruments>							
Visual examination	61	0	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	32	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	30	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	36	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	184	0	0	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0	0
Roche Chemstrips	54	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Siemens Acetest	16	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Ictotest	1	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Other	11	0	0	0	0	0	0
All Instruments and Reagents	406	0	0	0	0	0	0

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

Ketones	Specimen: U79	Negative	Trace 5 mg/dL	Small 1+ 15 mg/dL	Moderate 2+ 40-50 mg/dL	Large 3+ >=80 mg/dL	Positive
<Instruments>							
Visual examination	61	0	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	32	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	30	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	36	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	184	0	0	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0	0
Roche Chemstrips	54	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Siemens Acetest	16	0	0	0	0	0	0
Bayer Clinitest	15	0	0	0	0	0	0
Siemens Ictotest	1	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Other	11	0	0	0	0	0	0
All Instruments and Reagents	406	0	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Ketones	Specimen: U80	Negative	Trace 5 mg/dL	Small 1+ 15 mg/dL	Moderate 2+ 40-50 mg/dL	Large 3+ >=80 mg/dL	Positive
<Instruments>							
Visual examination	60	1	0	0	0	0	0
Siemens Atlas	88	0	0	0	0	0	0
Siemens Clinitek Status	32	0	0	0	0	0	0
Siemens Clinitek 50	11	0	0	0	0	0	0
Siemens Clinitek 100	16	0	0	0	0	0	0
Siemens Clinitek 500	72	0	0	0	0	0	0
Siemens Clinitek Advantus	30	0	0	0	0	0	0
Roche UA Reader/Super UA Read	5	0	0	0	0	0	0
Roche Mini UA Reader	1	0	0	0	0	0	0
Roche Chemstrip Criterion	14	0	0	0	0	0	0
Roche Urisys	36	0	0	0	0	0	0
IRIS AUTION MAX	29	0	0	0	0	0	0
Other Chemistry Analyzer	2	0	0	0	0	0	0
Other	9	0	0	0	0	0	0
<Reagents>							
Siemens Reagent Strips	184	0	0	0	0	0	0
Siemens Atlas Reagent	77	0	0	0	0	0	0
Roche Chemstrips	54	0	0	0	0	0	0
Roche Urisys Cassette	17	0	0	0	0	0	0
ARKRAY Aution Sticks	21	0	0	0	0	0	0
IRIS Reagent Strips	9	0	0	0	0	0	0
Siemens Acetest	16	0	0	0	0	0	0
Bayer Clinitest	14	1	0	0	0	0	0
Siemens Ictotest	1	0	0	0	0	0	0
Bayer Diastix	1	0	0	0	0	0	0
Other	11	0	0	0	0	0	0
All Instruments and Reagents	405	1	0	0	0	0	0

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Protein	Specimen: U76	Negative	Trace 10-15 mg/dL	1+ 25-30 mg/dL	2+ 75-150 mg/dL	3+ 300-600 mg/dL	4+ >=1000 mg/dL	Positive
<Instruments>								
Visual examination		0	0	1	11	66	22	36
Siemens Atlas		0	1	0	1	82	5	0
Siemens Clinitek Status		0	0	0	0	31	1	0
Siemens Clinitek 50		0	0	0	0	11	0	0
Siemens Clinitek 100		0	0	0	0	15	0	1
Siemens Clinitek 500		1	0	0	1	69	1	0
Siemens Clinitek Advantus		0	0	0	0	30	0	0
Roche UA Reader/Super UA Read		0	0	0	1	4	0	0
Roche Mini UA Reader		0	0	0	0	1	0	0
Roche Chemstrip Criterion		0	0	0	1	13	0	0
Roche Urisys		0	0	0	6	30	0	0
IRIS AUTION MAX		0	0	0	12	17	0	0
Hitachi/Roche Analyzer		0	0	0	2	0	0	0
Other Chemistry Analyzer		0	0	0	0	2	0	0
Other		0	0	0	1	9	0	0
<Reagents>								
Siemens Reagent Strips		1	0	0	3	176	8	0
Siemens Atlas Reagent		0	1	0	1	72	5	0
Sulfosalicylic Acid		0	0	0	7	32	9	37
Roche Chemstrips		0	0	0	8	46	1	0
Roche Urisys Cassette		0	0	0	2	15	0	0
ARKRAY Aution Sticks		0	0	0	8	13	0	0
IRIS Reagent Strips		0	0	0	4	5	0	0
Bayer Clinitest		0	0	1	0	13	0	0
Other		0	0	0	3	8	6	0
All Instruments and Reagents		1	1	1	36	380	29	37

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Protein	Specimen: U77	Negative	Trace 10-15 mg/dL	1+ 25-30 mg/dL	2+ 75-150 mg/dL	3+ 300-600 mg/dL	4+ >=1000 mg/dL	Positive
<Instruments>								
Visual examination		27	42	4	0	0	0	8
Siemens Atlas		85	3	1	0	0	0	0
Siemens Clinitek Status		9	15	8	0	0	0	0
Siemens Clinitek 50		3	6	2	0	0	0	0
Siemens Clinitek 100		10	5	0	0	0	0	1
Siemens Clinitek 500		70	2	0	0	0	0	0
Siemens Clinitek Advantus		29	1	0	0	0	0	0
Roche UA Reader/Super UA Read		4	0	1	0	0	0	0
Roche Mini UA Reader		0	1	0	0	0	0	0
Roche Chemstrip Criterion		9	5	0	0	0	0	0
Roche Urisys		15	7	13	1	0	0	0
IRIS AUTION MAX		18	11	0	0	0	0	0
Hitachi/Roche Analyzer		0	2	0	0	0	0	0
Other Chemistry Analyzer		2	0	0	0	0	0	0
Other		4	5	1	0	0	0	0
<Reagents>								
Siemens Reagent Strips		134	43	11	0	0	0	0
Siemens Atlas Reagent		76	2	1	0	0	0	0
Sulfosalicylic Acid		11	11	2	0	0	0	9
Roche Chemstrips		18	21	13	1	0	0	0
Roche Urisys Cassette		12	3	2	0	0	0	0
ARKRAY Aution Sticks		14	7	0	0	0	0	0
IRIS Reagent Strips		5	4	0	0	0	0	0
Bayer Clinitest		9	5	0	0	0	0	0
Other		6	9	1	0	0	0	0
All Instruments and Reagents		285	105	30	1	0	0	9

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Protein	Specimen: U78	Negative	Trace 10-15 mg/dL	1+ 25-30 mg/dL	2+ 75-150 mg/dL	3+ 300-600 mg/dL	4+ >=1000 mg/dL	Positive
<Instruments>								
Visual examination		46	21	2	0	0	0	2
Siemens Atlas		88	1	0	0	0	0	0
Siemens Clinitek Status		28	3	1	0	0	0	0
Siemens Clinitek 50		11	0	0	0	0	0	0
Siemens Clinitek 100		15	0	0	0	0	0	1
Siemens Clinitek 500		72	0	0	0	0	0	0
Siemens Clinitek Advantus		30	0	0	0	0	0	0
Roche UA Reader/Super UA Read		4	0	1	0	0	0	0
Roche Mini UA Reader		1	0	0	0	0	0	0
Roche Chemstrip Criterion		14	0	0	0	0	0	0
Roche Urisys		30	1	5	0	0	0	0
IRIS AUTION MAX		29	0	0	0	0	0	0
Hitachi/Roche Analyzer		2	0	0	0	0	0	0
Other Chemistry Analyzer		2	0	0	0	0	0	0
Other		9	1	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips		179	8	1	0	0	0	0
Siemens Atlas Reagent		78	1	0	0	0	0	0
Sulfosalicylic Acid		15	2	2	0	0	0	3
Roche Chemstrips		37	11	6	0	0	0	0
Roche Urisys Cassette		17	0	0	0	0	0	0
ARKRAY Aution Sticks		21	0	0	0	0	0	0
IRIS Reagent Strips		9	0	0	0	0	0	0
Bayer Clinitest		13	1	0	0	0	0	0
Other		12	4	0	0	0	0	0
All Instruments and Reagents		381	27	9	0	0	0	3

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Protein	Specimen: U79	Negative	Trace 10-15 mg/dL	1+ 25-30 mg/dL	2+ 75-150 mg/dL	3+ 300-600 mg/dL	4+ >=1000 mg/dL	Positive
<Instruments>								
Visual examination		3	37	60	2	0	0	30
Siemens Atlas		1	8	80	0	0	0	0
Siemens Clinitek Status		0	0	19	13	0	0	0
Siemens Clinitek 50		0	0	8	3	0	0	0
Siemens Clinitek 100		0	1	14	0	0	0	1
Siemens Clinitek 500		2	39	31	0	0	0	0
Siemens Clinitek Advantus		0	10	20	0	0	0	0
Roche UA Reader/Super UA Read		0	1	4	0	0	0	0
Roche Mini UA Reader		0	0	1	0	0	0	0
Roche Chemstrip Criterion		0	10	4	0	0	0	0
Roche Urisys		3	6	23	4	0	0	0
IRIS AUTION MAX		2	21	6	0	0	0	0
Hitachi/Roche Analyzer		0	2	0	0	0	0	0
Other Chemistry Analyzer		0	2	0	0	0	0	0
Other		0	4	5	1	0	0	0
<Reagents>								
Siemens Reagent Strips		2	55	113	18	0	0	0
Siemens Atlas Reagent		1	6	72	0	0	0	0
Sulfosalicylic Acid		3	27	20	0	0	0	31
Roche Chemstrips		1	19	31	4	0	0	0
Roche Urisys Cassette		2	3	12	0	0	0	0
ARKRAY Aution Sticks		1	15	5	0	0	0	0
IRIS Reagent Strips		1	7	1	0	0	0	0
Bayer Clinitest		0	4	9	1	0	0	0
Other		0	5	12	0	0	0	0
All Instruments and Reagents		11	141	275	23	0	0	31

New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010

Protein	Specimen: U80	Negative	Trace 10-15 mg/dL	1+ 25-30 mg/dL	2+ 75-150 mg/dL	3+ 300-600 mg/dL	4+ >=1000 mg/dL	Positive
<Instruments>								
Visual examination		64	4	0	0	0	0	0
Siemens Atlas		89	0	0	0	0	0	0
Siemens Clinitek Status		32	0	0	0	0	0	0
Siemens Clinitek 50		11	0	0	0	0	0	0
Siemens Clinitek 100		16	0	0	0	0	0	0
Siemens Clinitek 500		72	0	0	0	0	0	0
Siemens Clinitek Advantus		30	0	0	0	0	0	0
Roche UA Reader/Super UA Read		5	0	0	0	0	0	0
Roche Mini UA Reader		1	0	0	0	0	0	0
Roche Chemstrip Criterion		14	0	0	0	0	0	0
Roche Urisys		32	1	3	0	0	0	0
IRIS AUTION MAX		28	1	0	0	0	0	0
Hitachi/Roche Analyzer		2	0	0	0	0	0	0
Other Chemistry Analyzer		2	0	0	0	0	0	0
Other		10	0	0	0	0	0	0
<Reagents>								
Siemens Reagent Strips		188	0	0	0	0	0	0
Siemens Atlas Reagent		79	0	0	0	0	0	0
Sulfosalicylic Acid		20	0	0	0	0	0	0
Roche Chemstrips		48	2	3	0	0	0	0
Roche Urisys Cassette		17	0	0	0	0	0	0
ARKRAY Aution Sticks		21	0	0	0	0	0	0
IRIS Reagent Strips		8	1	0	0	0	0	0
Bayer Clinitest		12	2	0	0	0	0	0
Other		15	1	0	0	0	0	0
All Instruments and Reagents		408	6	3	0	0	0	0

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

Qualitative hCG report:

Bayer Clinitek Status	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	4	4	4	4	4
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Beckman Coulter ICON II	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	2	2	2	2	2
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Beckman Coulter ICON 25	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	29	29	29	29	29
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Beckman Coulter ICON 20	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	21	21	21	21	21
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Cardinal Health hCG Combo	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	10	10	10	10	10
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Chembio	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	1	1	1	1	1
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Wampole Clearview hCG Combo	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	4	4	4	4	4
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

Cen Med Elite	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	20	20	20	20	20
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Genzyme Contrast hCG	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	1	1	1	1	1
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Wampole Clearview hCG II	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	8	8	8	8	8
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Genzyme OSOM Combo	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	19	19	19	19	19
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Inverness Medical BioStar Acceava	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	2	2	2	2	2
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Lifesign Status (serum/urine)	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	1	1	1	1	1
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Polymedco Poly stat	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	7	7	7	7	7
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

Quidel QuickVue	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	72	72	72	71	69
Indeterminant	0	0	0	0	2
Positive	0	0	0	1	1
SA Scientific	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	2	2	2	2	2
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Stanbio QuPID	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	28	28	28	27	28
Indeterminant	0	0	0	0	0
Positive	0	0	0	1	0
Fisher Sure-Vue serum/urine	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	39	39	39	39	39
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Fisher Sure-Vue STAT serum/urine	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	29	29	29	29	28
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	1
Signify hCG Urine/Serum	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	4	4	4	4	4
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Sacks Medical RefuAH	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	5	5	5	5	5
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0

**New York State Department of Health – Wadsworth Center
Urinalysis/Urine Pregnancy Proficiency Testing – 19 April 2010**

SureStep (Applied Biotech)	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	1	1	1	1	1
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Fisher Sure-Vue (urine only)	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	7	7	7	7	7
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
True Medix	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	6	6	6	6	6
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
Method not listed	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	27	27	27	27	27
Indeterminant	0	0	0	0	0
Positive	0	0	0	0	0
All Methods & Instruments	Specimen:U76	Specimen:U77	Specimen:U78	Specimen:U79	Specimen:U80
Negative	349	349	349	347	345
Indeterminant	0	0	0	0	2
Positive	0	0	0	2	2