

<b>Transplant Monitoring</b>	
<b>Standard</b>	<b>Guidance</b>
The following specialty sustaining standards of practices shall be incorporated into the laboratory's quality management system, where applicable to the scope of services provided.	
<b>Engraftment Monitoring Standard 1 (EM S1)</b>  The laboratory shall include a sensitivity control in each patient run.	For this control, it is suggested that a small amount of a positive sample be mixed with an excess of a negative sample, e.g., 1:20 for a 5% sensitivity.
<b>Engraftment Monitoring Standard 2 (EM S2)</b>  The final report shall include, at minimum, a summary of the method that was used, the DNA loci tested, and the objective findings in a readily interpretable format.	

